EtOH AQ 0,3 g/L – In vitro diagnosticum

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AQ03-030 (10 x 3,0 ml)

Version: 4 – 201510

EtOH AQ 0.3 g/l – For in vitro diagnostic use

Aqueous Ethanol Standards

Application

This aqueous ethanol standard is intended for the calibration of analytical techniques used to determine the concentration of ethanol.

User guide

This ACQ Science EtOH AQ standard requires no additional preparation and is ready for use. This material should be used as defined by the individual laboratory's operating procedures and applied for instrument calibration or as a control material.

Assigned value

The assigned ethanol concentration was determined by 3 independent laboratories, each accredited to DIN EN 17025. Repeat determinations were carried out daily on 5 days using two independent analytical methods (Gas Chromatography and Enzymatic determination (ADH)). The material is traceable to primary ethanol standards (ERM¹ ethanol standard in aqueous solution).

Storage and stability

Storage: 2 ° to 8 ° C

Stability:

- Sealed container, stored in the dark: see expiration date on the package.
- Stored tightly capped in the dark: see expiration date on package.

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Confidence range (of the target value)

The confidence interval indicates the range in which the target value is located with a significance level of 95%.

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ACQ Science GmbH