RIDASCREEN® Histamin

REF R1601

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RIDASCREEN® Histamine

Brief information

RIDASCREEN® Histamine (Art. No. R1601) is a competitive enzyme immunoassay for the quantitative analysis of histamine in foodstuffs validated for the method (see 1. Intended use).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: different depending on the sample (see 9.)

Time requirement: sample preparation (for 10 samples) approx. 15 min
test implementation (incubation time) 90 min
test implementation (incubation with shaker) 70 min

Limit of detection: sparkling wine, white and red wine 250 ppb
milk 100 ppb
cheese, fresh and canned fish 2.5 ppm
fish meal 100 ppm

Recovery rate: sparkling wine, white and red wine 95 %
milk 109 %
cheese, fresh and canned fish 100 %
fish meal 96 %

Specificity: Histamine 100 %
3-Methyl-histamine approx. 0.01 %
Tyramine n.d.
L-Phenylalanine n.d.
L-Histidine n.d.
L-Tyrosine n.d.
Tryptamine n.d.
5-Hydroxy-Indole-Acetic Acid n.d.
Serotonin n.d.

Related products for histamine determination

RIDASCREEN® Histamine (enzymatic) (Art. No. 1601)

1. Intended use

RIDASCREEN® Histamin (Art. No. R1601) is a competitive enzyme immunoassay for the quantitative analysis of histamine in white and red wine, sparkling wine, milk, cheese, fresh fish, canned fish and fish meal. It can be assumed that the test is also suitable for the analysis of other foodstuffs; this must be checked by the user.

2. General

Histamine is a product of the decomposition of histidine caused by the growth of certain bacteria in protein rich food like fish, cheese, meat as well as sparkling wine, wine and beer. The amount of histamine formed depends on the bacterial species, the temperature and time of exposure, and may exceed 1000 ppm (mg/kg). Fish containing high levels of histamine have been associated with poisoning commonly referred to as "scombroid poisoning". Quality control measures have to be taken to minimize the occurrence of histamine contaminated fish. Good quality fish contains less than 10 ppm histamine. The limit value for fish and products thereof is 100 mg/kg, in some countries even 50 mg/kg. The Swiss limit value for histamine in wine is 10 mg/l.

3. Test principle

After sample preparation histamine is derivatized quantitatively by an acylation reagent into N-acylhistamine. The competitive ELISA is based on an antigen-antibody reaction. The wells of the microtiter strips are coated with histamine. Standards or samples and anti-histamine antibodies are added. Free acylated histamine and bound histamine compete for the antibody binding sites (competitive enzyme immunoassay). After washing, secondary antibodies labeled with peroxidase (enzyme conjugate) are added. These antibodies bind to the antibody-histamine complexes. Any unbound enzyme conjugated antibody is then removed in a washing step. Enzyme substrate and chromogen are added to the wells and incubated. Bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the histamine concentration in the sample.

The results is expressed in ng/ml or µg/kg (ppb) histamine.

4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format	Concentration	Volume
Reagents for sample acylation				
Acylation plate	-	ready to use		96 wells
Acylation tubes	-	ready to use		96 tubes (plastic)
Standard 1	white	ready to use	0 ng/ml	4 ml
Standard 2	light yellow	ready to use	0.5 ng/ml	4 ml
Standard 3	light red	ready to use	1.5 ng/ml	4 ml
Standard 4	black	ready to use	5.0 ng/ml	4 ml
Standard 5	grey	ready to use	15.0 ng/ml	4 ml
Standard 6	black	ready to use	50.0 ng/ml	4 ml
Control 1	yellow green	ready to use	see bottle label	4 ml
Control 2	dark red	ready to use	see bottle label	4 ml
Acylation reagent	green	ready to use		3 ml
Acylation buffer	brown	ready to use, red-colored		22 ml
Reagents for the enzyme immunoassay				
Microtiter plate	-	ready to use		96 wells
Wash buffer	rose	concentrate	50x	20 ml
Conjugate	red	ready to use		12 ml
Anti-histamine Antibody	blue	ready to use		12 ml
Substrate/Chromogen	black	ready to use		12 ml
Stop solution	grey	ready to use		12 ml

5. Reagents required but not provided

5.1. Equipment

- Microtiter plate spectrophotometer (450 nm)
- Gloves
- Scale
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge + centrifugal vials (plastic)
- Shaker/rotator
- Graduated pipettes
- Graduated cylinder
- Variable 20 - 200 µl and 200 - 1000 µl micropipettes
- To shorten incubation time: microtiter plate (MTP) shaker (shaking amplitude 3 mm, approx. 600 rpm)
- If necessary: 8-channel pipette for 100 µl
- Optional: RIDASOFT® Win.NET (Z9996FF)

5.2. Reagents

- Distilled water (dist. water) or deionized water
- For milk samples: 10 mM PBS, 0.05 % Tween 20 (e.g. Sigma P3563)

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the plates. Use separate pipet tips for each standard and each sample extract to avoid cross contamination.

Ensure the proper and responsible disposal of all reagents and materials after their use. For disposal, please adhere to national regulations.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Any bluish coloration of the substrate/chromogen prior to test implementation.
- A value of less than 0.9 absorbance units ($A_{450\text{ nm}} < 0.9$) for standard 1.

9. Preparation and extraction of samples

Before starting and during the assay wear gloves. Impure laboratory equipment may lead to contamination of the assay. Hence, it is recommended to take the following precautionary measures:

- Clean surfaces, glass vials, mincers and other equipment before and after each sample preparation,
- Carry out the sample preparation in a room isolated from the assay procedure,

Acylation tubes or reaction wells of acylation plate are for single use only to avoid contaminations.

For sample preparation use only plastic vials!!!

9.1. Sparkling wine, white wine, red wine

- dilute 20 µl wine with 10 ml of dist. water
- apply 100 µl per acylation tube or per well in the acylation plate

9.2. Milk

- centrifuge 4 ml milk for 10 min / 3,000 x g at room temperature (20 - 25 °C / 68 - 77 °F)
- remove the lipid layer
- mix 20 µl defatted milk with 4 ml PBS buffer (see 5.2.)
- apply 100 µl per acylation tube or per well in the acylation plate

9.3. Cheese, fish (fresh and canned fish)

- homogenize 10 g of sample
- add 9 ml of distilled water to 1g of homogenate and mix well
- centrifuge: 5 min / 2,500 x g at room temperature (20 - 25 °C / 68 - 77 °F)
- remove lipid layer or cut through lipid layer
- mix 1 ml of supernatant with 9 ml of water and mix well
- dilute 200 µl of this solution with 9.8 ml of dist. water
- apply 100 µl per acylation tube or per well in the acylation plate

9.4. Fish meal

- mix 1 g of fish meal with 200 ml of dist. water and stir for 15 min
- centrifuge aliquot: 5 min / 2,500 x g at room temperature (20 - 25 °C / 68 - 77 °F)
- dilute 200 µl of the supernatant with 200 ml of dist. water
- apply 100 µl per acylation tube or per well in the acylation plate

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **washing buffer** is provided as a 50fold concentrate. The buffer concentrate (20 - 25 °C / 68 - 77 °F) is diluted 1:50 (1+49) with dist. water before use (i.e. 10 ml buffer concentrate + 490 ml distilled water). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. 4 weeks.

The **acylation reagent** is provided ready to use, it has a freezing point of 18.5 °C (65 °F). When being used, it must be ensured that the acylation reagent has reached room temperature and forms a homogeneous, crystal-free solution. The acylation reagent can be stored at room temperature (20 - 25 °C / 68 - 77 °F) or brought to room temperature before use.

10.2. Test procedure for the acylation

Pipette carefully and avoid splashes.

1. All standards, controls and samples should be run in duplicate. Record standard and sample positions on the acylation plate.
Alternatively, insert required number of acylation tubes into the acylation plate and use as reaction vials.
2. Add 100 µl of each standard solution, control or prepared sample to acylation tubes or separate wells of the acylation plate.
3. Add 25 µl of the acylation reagent to each acylation tube or well.
4. Add 200 µl of the acylation buffer to each acylation tube or well, mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Use 25 µl for the ELISA.

10.3. Test procedure for ELISA

Pipette carefully and avoid splashes. Careful washing is very important. Avoid drying of the wells between the steps.

It is recommended to pipette the antibody, conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards, controls and samples to be run in duplicate. Record standard and sample positions.
2. Add 25 µl of acylated standard solution, control or prepared sample to separate wells.
3. Add 100 µl of the anti-histamine antibody solution to each well, mix gently by shaking the plate manually and incubate for 40 min at room temperature (20 - 25 °C / 68 - 77 °F). Alternatively: 30 min with a MTP-shaker (approx. 600 rpm)
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl of washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
5. Add 100 µl of the conjugate solution to of each well. Mix gently by shaking the plate manually and incubate for 20 min at room temperature (20 - 25 °C / 68 - 77 °F). Alternatively: 10 min with a MTP-shaker (approx. 600 rpm)

6. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl of washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
7. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark. Alternatively: 15 min with a MTP-shaker (approx. 600 rpm).
8. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm against an air blank. Read within 10 min after addition of stop solution.

11. Results

A special software, the **RIDASOFT® Win.NET (Art. Nr. Z9996FF)**, is available for evaluation of the RIDASCREEN® enzyme immunoassays.

For the evaluation it must be clarified, that all quality criteria are fulfilled for the current test run. The criteria as well as the course of the standard curve can be taken from the attached certificate of analysis.

The histamine concentration in µg/kg (ppb) corresponding to the absorbance of each sample is read from the calibration curve and then further multiplied by the corresponding dilution factor. The dilution factor depends on the sample preparation and is listed in the table below.

Sample preparation	Sample	Dilution factor
9.1.	red wine, white wine, sparkling wine	500
9.2.	milk	200
9.3.	cheese, fresh and canned fish	5,000
9.4.	fish meal	200,000

12. Limits of the method

An incorrect weight of the sample to be analysed has an influence on the measurement results. A sufficient accuracy is given with a fluctuation of max. +/- 1 %.

A further dilution and new detection of samples is recommended for absorbance values ($A_{450\text{nm}}$) > 0.4. In case of a further dilution, the additional dilution factor must be taken into account when calculating the histamine concentration.

For the present ELISA, only individual, exemplary foods from different product groups could be validated due to the large number of foods. When analysing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spike experiments.

Detection and determination limits depend on the respective sample matrix, the degree of processing and the extraction method.

13. Recommendation

In order to ensure a high analytical performance we recommend:

- Follow the general quality assurance requirements for laboratories (e.g. performing duplicate determinations).
- Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors).
- Rinse pipette tips with standard or sample extract before pipetting.
- Test controls should be carried along for quality control. The controls contained in the kit should be used for this purpose.
- The controls should be measured with a maximum deviation of +/- 30 % (otherwise repeat the assay).
- The optimal temperature for using the assay is at room temperature (20 - 25 °C / 68 - 77 °F).

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2010-06-14	Version 1
2020-08-28	Version 2 Applied changes: 4. Reagents provided General linguistic revision

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG.
R-Biopharm AG