# RIDASCREEN® Chinolone/Quinolones

**REF R3113**

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# RIDASCREEN® Chinolone/Quinolones

## Brief information

RIDASCREEN® Chinolone/Quinolones (R3113) is a competitive enzyme immunoassay for the quantitative analysis of quinolones in egg, meat, fish and shrimp.

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for max. 48 duplicate determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: Egg, shrimps: homogenization, extraction, centrifugation  
Meat, fish: homogenization, extraction, centrifugation, dilution

Time requirement: Sample preparation (for 10 samples).... approx. 30 min  
Test implementation (incubation time) ..... 1 h 15 min

Limit of detection: Shrimps ..... 6 µg/kg (ppb)  
(corresponding to the standard substance) Fish..... 8 µg/kg  
Egg ..... 9 µg/kg  
Meat..... 10 µg/kg

Recovery rate: 80 - 110 %  
(corresponding to the Standard substance) Further information is contained in the validation report.

Specificity: Ciprofloxacin (standard substance).....100 %  
Norfloxacin, Enrofloxacin,  
Marbofloxacin, Danofloxacin,  
Difloxacin, Flumequine, Ofloxacin .....> 100 %  
Sarafloxacin .....43 %  
Oxolinic acid .....24 %

The specificity of the RIDASCREEN® Chinolone/Quinolones test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice manual. It lists minimum standards concerning the framework conditions when using test kits of R-Biopharm AG and performing ELISA analyses with them. The manual can be retrieved, printed and downloaded from [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

## **Related products**

RIDA® Ciprofloxacin Spiking Solution (R3198)

### **1. Intended use**

RIDASCREEN® Chinolone/Quinolones is a competitive enzyme immunoassay for the quantitative analysis of quinolones in egg, meat (beef, pork, sheep, chicken, turkey), fish and shrimps.

### **2. General information**

Gyrase inhibitors are divided into four subgroups. The majority of quinolones belong to the subgroup of fluoroquinolones, which have a fluoro-group attached at the central ring system, typically at the 6<sup>th</sup> position. The fluoroquinolones belong to the so called second generation quinolones. Some of these are admitted as antibiotics in veterinary medicine for food producing animals. Fluoroquinolones are broad spectrum antibiotics against a lot of bacterial species. They are used frequently in veterinary medicine especially for cattle, pigs and chicken. The usage of fluoroquinolone has increased in the past years because large amounts of fluoroquinolones were applied to prevent infectious diseases, especially in chicken, swine and fish/shrimp farming.

The wide use in food producing animals has generated microbial resistances. Therefore, in the European Union, maximum residue limits (MRL) apply to some substances of this group, as described in COMMISSION REGULATION (EU) No. 37/2010.

### 3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-quinolone antibodies. Standards or sample solution, ciprofloxacin enzyme conjugate and anti-quinolone antibodies are added. Free quinolones and ciprofloxacin conjugate compete for the quinolone antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-quinolone antibodies are also bound by the immobilized capture antibodies. Any unbound conjugate is then removed in a washing step. Substrate/Chromogen is added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue end product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the quinolone concentration in the sample.

### 4. Reagents provided

Each kit contains sufficient materials for max. 48 duplicate measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format	Volume
<b>Microtiter plate</b>	-	Ready to use	96 wells
<b>Standard 1</b>	White	Ready to use	0 µg/l 1.3 ml
<b>Standard 2</b>	White	Ready to use	0.5 µg/l 1.3 ml
<b>Standard 3</b>	White	Ready to use	1.5 µg/l 1.3 ml
<b>Standard 4</b>	White	Ready to use	3 µg/l 1.3 ml
<b>Standard 5</b>	White	Ready to use	6 µg/l 1.3 ml
<b>Standard 6</b>	White	Ready to use	18 µg/l 1.3 ml
<b>Wash buffer salt Tween</b>		Salt for Dissolving	
<b>Conjugate</b>	Red	Ready to use	6 ml
<b>Antibody</b>	Black	Ready to use	6 ml
<b>Substrate/Chromogen</b> Red Chromogen Pro	Brown	Ready to use	10 ml
<b>Stop solution</b>	Yellow	Ready to use	14 ml

## 5. Reagents required but not provided

### 5.1 Equipment

- Microtiter plate spectrophotometer (450 nm)
- Centrifuge
- Up-side-down shaker or vortex
- Graduated pipettes
- Variable 20 - 200 µl and 200 - 1000 µl micropipettes

### 5.2 Reagents

- Methanol
- Shrimp extraction buffer: dilute methanol/water (70:30, v/v) 1:2 (1+1) with wash buffer (see chapter 10.1.)

## 6. Warnings and precautions for the users

This test should be carried out only by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. Please refer to the component safety information in the material safety data sheets (SDS) for this product, available online at [www.r-biopharm.de](http://www.r-biopharm.de).

Ensure the proper and responsible cleanup and / or disposal of all reagents and materials after their use. For disposal, please adhere to national regulations.

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to addition in the wells
- Extinction less than 0.6 ( $E_{450\text{ nm}} < 0.6$ ) for standard 1

## 9. Sample preparation

Working devices such as a mill, glass vials or spatulas must be cleaned before.

The samples should be stored in a cool place, protected against light.

### 9.1. Egg

- Homogenize the sample
- Use 1 g of the homogenized sample and add 9 ml methanol/water (35/65, v/v)
- Mix vigorously for 10 min (with up-side-down shaker or vortex)
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- Use 50 µl of the supernatant per well in the test

### 9.2. Meat (beef, pork, sheep, chicken, turkey) and fish

- Homogenize the sample
- Use 1 g of the homogenized sample and add 4 ml methanol/water (70/30, v/v)
- Mix vigorously for 10 min (with up-side-down shaker or vortex)
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- Dilute the supernatant 1:2 (1+1) with wash buffer (see chapter 10.1.)
- Use 50 µl per well in the test

### 9.3. Shrimp

- Homogenize the sample
- Use 1 g of the homogenized sample and add 4 ml shrimp extraction buffer (see chapter 5.2.)
- Mix vigorously for 10 min (with up-side-down shaker or vortex)
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- Use 50 µl of the supernatant per well in the test

## 10. Test procedure

### 10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

A PBS-Tween buffer is needed as **wash buffer**, please use the wash buffer salt (pouch) contained in the kit (see chapter 4.). Dissolve the total content of the pouch in one liter of distilled water. The ready to use wash buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternatively: Dissolve the contents of the pouch in 100 ml of distilled water to obtain a 10fold concentrated wash buffer. This solution expires after approx. 8 - 12 weeks, stored at room temperature (20 - 25 °C / 68 - 77 °F). Use one part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use wash buffer (e.g. 100 ml concentrate + 900 ml distilled water).

Unused reagents should be immediately stored at 2 - 8 °C.

### 10.2 Test procedure

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow wells to dry between work steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µl of each standard or prepared sample to separate duplicate wells. Use a new pipette tip for each standard or sample.
3. Add 50 µl of the conjugate (red cap) to the bottom of each well.
4. Add 50 µl of the antibody (black cap) to each well, mix gently by shaking the plate manually and incubate for 1 h in a refrigerator (2 - 8 °C / 35 - 46 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see chapter 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl of substrate/chromogen (brown cap) to each well. Mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.

7. Add 100 µl of the stop solution (yellow cap) to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 minutes after addition of stop solution.

## 11. Evaluation

Special software, **RIDASOFT® Win.NET** (Z9996FF), is optionally available for evaluating RIDASCREEN® enzyme immunoassays. The course of the standard curve is shown in the Quality Assurance Certificate (certificate of analysis) enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = B/B_0 (\%)$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a semi-logarithmic plot against the quinolone concentration [µg/kg].

## 12. Result interpretation

In order to obtain the quinolone concentration in µg/kg (ppb) actually contained in a sample, the concentration read from the standard curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factors are as follows:

Egg .....	10
Meat, fish .....	10
Shrimp .....	5

## 13. Recommendation

In order to ensure a high analytical performance we recommend:

- To analyze each sample material in duplicates
- To use also quinolone-free and quinolone containing (spiked) samples as test controls
- To perform spike experiments to ensure an accurate and correct test procedure

**Further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

## Version overview

Version number	Chapter and title
2008-05-20	Release version
2015-08-21	General revision
2021-02-15	General revision Correction in “Related products” NEW: Chapter 13. Recommendation, Version overview, Explanation of symbols

## Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

## Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.