RIDA® QUICK CIS

REF R4303

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RIDA®QUICK CIS

Brief information

RIDA®QUICK CIS (R4303) is an immunochromatographic test for the determination of cow's milk in milk or cheese of other species (sheep and goat).

The test kit contains 25 test strips, each of which can be used for one analysis. All reagents required for the assay are contained in the test kit. Results are read and evaluated visually.

Sample preparation:	Milk: use directly in the test Cheese: homogenization and filtration
Time requirement:	Sample preparation (for 10 samples) Milk:none Cheese: approx. 10 min Test implementation (incubation time) 5 min
Limit of detection:	approx. 0.5 % cow's milk in sheep's and goat's milk approx. 0.5 % cow's milk in sheep's and goat's cheese
Specificity:	The use of monoclonal antibodies in the elaboration of the RIDA®QUICK CIS test assures its high degree of specificity for bovine IgG.

1. Intended use

RIDA®QUICK CIS test is an immunochromatographic test for the detection of cow's milk in milk or cheese of other species (sheep and goat).

2. General information

Sheep's milk and goat's milk are much more expensive than cow's milk which is often used to extend sheep's and goat's milk. This means an adulteration of the product according to EU-law.

Adulterations of milk can be detected chromatographically or by electrophoresis. These methods are laborious and require sophisticated equipment. Therefore, they are of limited value in routine screening of milk or cheese.

3. Test principle

The basis of the assay is the antigen-antibody reaction. The presence of cow's milk in a sample is determined by the immunological detection of bovine IgG (an antibody class) which is a natural constituent of cow's milk. In ultra high heated products the antibodies will be denatured. Therefore an adulteration with cow's milk cannot be detected with the RIDA®QUICK CIS test in these kind of products.

Each RIDA®QUICK CIS reaction strip contains a membrane, which is coated with specific antibodies directed against bovine IgG. The reagents are migrating through the membrane of the reaction strip. In case of falsification of the sample with cow's milk, bovine IgG (antigen) will bind to red colored beads, coated with antibodies against bovine IgG. The bead-antigen complex will be carried to the part of the membrane, which is surfaced with antibodies against bovine IgG. The red colored bead-antigen complex will be bound by the antibodies against bovine IgG, and a red band becomes visible in the field of reaction of the test strip. Unbound colored beads pass the field of reaction and will be absorbed in the absorption pad at the end of the membrane. A second blue upper band becomes visible as a control band for a correct working test strip.

4. Reagents provided

Each kit contains sufficient materials for max. 25 determinations. Each test kit contains:

Component	Volume
25 x Test strips	One for each determination
25 x Disposable pipettes	
25 x Test tubes	
6 x Tube holder	
Buffer	10 ml, ready to use

5. Reagents required but not provided

5.1 Equipment

- Stomacher or mortar (for cheese samples only)
- Graduated pipettes
- Variable 20 - 200 µl micropipettes

5.2 Reagents

- Distilled or deionized water
- PBS-buffer, pH 7.2: 0.55 g $\text{NaH}_2\text{PO}_4 \times \text{H}_2\text{O}$ + 2.85 g $\text{Na}_2\text{HPO}_4 \times 2 \text{H}_2\text{O}$ + 9 g NaCl, ad 1000 distilled water

6. Warnings and precautions for the users

This test should be carried out only by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. Please refer to the component safety information in the material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Ensure the proper and responsible disposal of all reagents and materials after their use. For disposal, please adhere to national regulations.

7. Storage instructions

Store the reagents dry at 2 - 30 °C (36 - 86 °F). Do not freeze the test strips or any test kit components.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. For this reason, it is important to protect the strips against moisture and to take them out of their container just for **direct use**. Return any unused reaction strips to their original container and close the container always tightly. Store the container at room temperature (20 - 25 °C / 68 - 77 °F).

Do not use RIDA®QUICK CIS after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Sample preparation

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

The samples should be stored in a cool place.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

8.1. Milk samples

Milk samples can be stored in the refrigerator for 1 or 2 days. For longer storage they must be kept frozen at -20 °C (-4 °F). Frozen samples should be totally thawed and brought to room temperature (20 - 25 °C / 68 - 77 °F) before testing.

- Use 25 µl of the milk sample per test strip

8.2. Cheese samples

- Mix and homogenize 1 g of cheese with 5 ml PBS-buffer (see chapter 5.2.)
- Filter the homogenate and remove gross pieces
- Use 25 µl of the filtrate per test strip

9. Test procedure

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature. Return all test strips into the container and store them at room temperature (20 - 25 °C / 68 - 77 °F), (see chapter 7.)

1. Introduce as many test tubes as samples to be analyzed in the test tube holder.
2. Using the provided disposable pipette, place 5 drops (alternatively 200 µl) of the buffer solution in the test tube.
3. Place 25 µl of the sample in the test tube and shake slightly.
4. Place the test strip inside vertically with the arrow end into the tube.
5. Read the result after 5 minutes without further action.

10. Evaluation

The upper band (blue line) in the reaction field is a control band (control line) and must be present after each test procedure. If this band is missing, the test was not performed correctly or the reagents were not functional. In this case, the test should be repeated with a new strip. In case of repeated absence of the upper band, please inform your local distributor.

The lower band is the test band (test line). The sample contains bovine IgG and therefore cow's milk if both the control band (control line) and the test band (test line) are present.

To obtain the maximum possible sensitivity of the test, the result should be read after exactly **5 minutes**. The result should not be interpreted after this time to avoid incorrect conclusions.

Any line or staining that should appear after 15 minutes is of no significance.

Remark:

For documentation the upper part of the strip marked with „bovine IgG“ can be cut off and the strip can be stored. By cutting the upper end of the strip the reaction is also stopped.

10.1 Positive samples

Two colored bands (the upper blue band and the lower reddish band) appearing in the field of reaction on the strip mean that the sample contains bovine IgG, the sample shows a falsification with cow's milk. The color intensity of the spot depends on the concentration of cow's milk adulteration in the sample.

With the RIDA®QUICK CIS it is possible to detect cow's milk adulteration up to 0.5 %. A positive result means that the sample has a cow's milk content of 0.5% or more.

10.2 Negative samples

If only **one colored band** (the upper blue colored control band) appears, the sample used does not contain cow's milk (< 0.5 %) and has therefore not been adulterated.

11. Limits of the method

The presence of cow's milk in a sample is determined by the immunological detection of bovine IgG (an antibody class) which is a natural constituent of cow's milk. In ultra high heated products the antibodies will be denatured. Therefore, an adulteration with cow's milk cannot be detected with the RIDA®QUICK CIS test in these kind of products.

Furthermore, a negative result does not exclude the possibility of cow's milk product adulteration below the detection limit of this test.

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2001-05-02	Release version
2010-03-04	General revision
2021-02-17	General revision In chapter 3.: product reference removed Removed: chapters “11. Sensitivity” and “12. Specificity” New chapter “11. Limits of the method”

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.