# RIDASCREEN®FAST Casein

**REF** R4612

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# RIDASCREEN®FAST Casein

## Brief information

RIDASCREEN®FAST Casein (Art. No. R4612) is a sandwich enzyme immunoassay for the quantitative analysis of casein in food validated for the method (see chapter 1. Intended Use).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization and extraction

Time requirement: sample preparation:  
for 10 samples..... approx. 20 min  
for 10 heated samples..... approx. 30 min  
test implementation (incubation time)..... 30 min

Standard material: Casein from SIGMA Aldrich (Art. No. C5890)

Limit of detection: Allergen extraction buffer 0.12 mg/kg (ppm) Casein\*  
(depending on matrix) ..... 0.19 mg/kg (ppm) Casein  
RIDA® Extractor 2..... 0.71 mg/kg (ppm) Casein\*  
..... 0.41 - 0.95 mg/kg (ppm) Casein  
\*mean value

Limit of quantification: Allergen extraction buffer 0.5 mg/kg (ppm) Casein  
RIDA® Extractor 2..... 2.5 mg/kg (ppm) Casein

Specificity: The antibodies used specifically detect caseins of cow's milk. There is a cross reactivity to sheep's, goat's and buffalo's milk. Any further cross reactivity is not known. Especially, a cross reactivity to  $\beta$ -lactoglobulin does not exist. Further information regarding cross reactivity is contained in the validation report.

Cross reactivities of the antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spiking experiments (see chapter 13. Limits of the method).

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

## **Related products for milk determination**

RIDA® Extractor 2 (Art. No. R4613)

RIDASCREEN®FAST Milk (Art. No. R4652)

RIDASCREEN®FAST  $\beta$ -Lactoglobulin (Art. No. R4912)

RIDASCREEN®  $\beta$ -Lactoglobulin (Art. No. R4901)

Bioavid Lateral Flow Milch / Milk (Art. No. BL613-10, BL613-25)

### **1. Intended use**

RIDASCREEN®FAST Casein (Art. No. R4612) is a sandwich enzyme immunoassay test developed for the quantitative analysis of casein in food. Due to the large number of different food products, the following samples were examined as representative for different food product categories within the scope of the test development: ice cream, wine, chocolate, sausage and cookies. It can be assumed that the assay is also suitable for the analysis of other foods, however, this must be verified by the user before applying the test kit to that food category.

For detailed results and further information on validation data with other food matrices, please refer to the validation report. Further applications are regularly validated in our laboratories, which we make available in our application notes (see chapter 15. Further application notes).

## 2. General information

Casein refers to a family of proteins commonly found in mammalian milk. The term “casein” is derived from the Latin word “*caseus*”, which means “cheese”. Cow’s milk contains 3.2 % proteins, which consist of approx. 10 %  $\beta$ -lactoglobulin (leading protein of whey), approx. 80 % caseins and approx. 10 % other proteins. The most important allergen for children is  $\beta$ -lactoglobulin while the caseins become the dominant allergen later in adults. Caseins are also heat-stable allergens.

Cow's milk is one of the most important allergenic components of food, especially for children. That is why the labeling of casein or milk is mandatory in many countries worldwide. Although there are no legal limits for casein, food producers are strongly advised to test for very low casein concentrations to protect people with allergies and to avoid allergen-related recalls.

Furthermore, Casein is a protein that is often used as a fining agent in wine production. In the fining process, casein is added to bind the suspended particles that settle out of the wine.

Casein is a rough flaked curdling protein, which forms micelles in the milk and precipitates under acidic conditions. The group of caseins consists of the  $\alpha$ -caseins,  $\beta$ -casein,  $\kappa$ -casein and  $\gamma$ -casein (proteolytic protein fragment of  $\beta$ -casein by the milk protease plasmin).  $\kappa$ -casein can be cleaved into a hydrophobic (para  $\kappa$ -casein) and into a water soluble polar component (macropeptid) by proteolysis, e.g. by using lab ferment. Caseins can be present as an ingredient or as a contaminant in raw and processed food products. Skim milk powder or caseins / caseinates are often added to food products (e.g. in sausages) to increase the protein content or as thickening agent. According to the **regulation (EU) No. 1169/2011**, milk and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

## 3. Test principle

The principle of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies against casein. By adding the standard or sample solution to the wells, casein present in the sample will bind to the specific capture antibodies resulting in the formation of an antibody-antigen-complex. Components not bound by the antibodies are then removed in a washing step. Following the washing step, a solution containing antibodies conjugated to peroxidase is added. This conjugate binds to the Ab-Ag-complex and an antibody-antigen-antibody (sandwich) complex is formed. Any unbound

conjugate is then removed in another washing step. A substrate/chromogen solution is added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue product. A stop solution is added which results in a color change from blue to yellow. The absorbance of the solution which is proportional to the casein concentration in the sample is measured photometrically at 450 nm and expressed as mg/kg casein.

#### 4. Reagents provided

Each kit contains sufficient materials for a maximum of 48 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
<b>Microtiter plate</b>	-	Ready to use		48 wells
<b>Extractor 2</b>	Blue	<b>Concentrate</b>	<b>2x</b>	30 ml
<b>Allergen extraction buffer</b>	Green	<b>Concentrate</b>	<b>10x</b>	100 ml
<b>Additive 1</b>	Blue			2 g
<b>Standard 1*</b>	Transparent	Ready to use	0 mg/kg	1.3 ml
<b>Standard 2*</b>	Transparent	Ready to use	0.5 mg/kg	1.3 ml
<b>Standard 3*</b>	Transparent	Ready to use	1.5 mg/kg	1.3 ml
<b>Standard 4*</b>	Transparent	Ready to use	4.5 mg/kg	1.3 ml
<b>Standard 5*</b>	Transparent	Ready to use	13.5 mg/kg	1.3 ml
<b>Wash buffer</b>	Brown	<b>Concentrate</b>	<b>10x</b>	100 ml
<b>Conjugate</b>	Red	<b>Concentrate</b>	<b>11x</b>	0.7 ml
<b>Conjugate buffer</b>	Black	Ready to use		7 ml
<b>Substrate/Chromogen</b> Red Chromogen Pro	Brown	Ready to use		13 ml
<b>Stop solution</b>	Yellow	Ready to use		14 ml

\* The concentration values of the standards already take into account the dilution factor 20, which results from the sample preparation according to chapter 9.2. Thus, the casein concentrations of the samples can directly be read from the standard curve.

If sample preparation with Extractor 2 is used (see chapter 9.1), a dilution factor of 100 is obtained. The values read from the standard curve must be multiplied by factor 5 in addition.

## 5. Reagents required but not provided

### 5.1. Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of  $\pm 0.01$  g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials (e.g. 50 ml centrifuge tubes from Greiner Art. No. 227261)
- Shaker
- Water bath (60 °C / 140 °F and 100 °C / 212 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Fluted filter (pore size 8 - 12  $\mu\text{m}$ )
- Graduated pipettes
- Graduated cylinders
- Variable 20 - 200  $\mu\text{l}$  and 200 - 1000  $\mu\text{l}$  micropipettes
- If necessary: a further microtiter plate (e.g. universal binding, breakable MTP from Thermo Fisher Scientific Art. No. 95029390 or low binding Greiner bio-one Art. No. 655901)
- If necessary: 8-channel pipette for 100  $\mu\text{l}$
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET (Art. No. Z9996FF)

### 5.2. Reagents

- Distilled water (dist. water) or deionized water
- 1 M sodium hydroxide (NaOH)
- 1 M hydrochloric acid (HCl)
- If necessary: bovine serum albumin (BSA, e.g. Serva, fraction V, protease free, Art. No. 11926)

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

Do not reuse wells of the microtiter strips (coated microtiter plate and pre-plate, if necessary (see chapter 10.2. Test procedure). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

The Extractor 2 is harmful to health. It contains mercaptoethanol. It should be worked **under a chemical hood** and skin contact should be avoided (use gloves!).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation
- Value of less than 1.2 absorbance units ( $A_{450\text{ nm}} < 1.2$ ) for standard 5

## 9. Preparation of Samples

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment before and after each sample preparation.
- Carry out the sample preparation in a room isolated from the ELISA procedure.

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The Allergen extraction buffer (AEB) is provided as a 10-fold concentrate and must be diluted prior use. Before dilution of the buffer concentrate, dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that, dilute the heated buffer concentrate 1:10 (1+9) with dist. water (e.g. 100 ml buffer concentrate + 900 ml dist. water). The **diluted Allergen extraction buffer (AEB)** is stable at 20 - 25 °C (68 - 77 °F) for approx. 4 weeks or at 2 - 8 °C (35 - 46 °F) for 12 weeks.

The AEB is required for the dilution of the extracts for use in ELISA. For the dilution of matrices containing pine, sunflower or pumpkin seeds, the diluted AEB must be additionally mixed with BSA (final concentration: 2.5 %; e.g. 10 ml diluted AEB + 0.25 g BSA (= BSA-AEB)).

For the preparation of the final **Allergen extraction buffer containing Additive 1 (A-AEB)**, which is required for the extraction according to chapter 9.1., weigh 1.35 g of Additive 1 in a glass beaker and add 15 ml 1 M NaOH. Stir until the Additive 1 is solved. Then, fill 700 ml diluted AEB (see above) in a measuring cylinder. Add the 15 ml Additive 1 solution by stirring constantly. Transfer any residues of the Additive 1 solution into the measuring cylinder by rinsing with diluted AEB. Adjust the Additive 1 containing Allergen extraction buffer (A-AEB) to pH 9 with 1 M HCl and fill up to 750 ml with diluted AEB.

750 ml A-AEB are sufficient for 45 samples. The buffer can be used for approx. 3 weeks at room temperature 20 - 25 °C (68 - 77 °F) (do **not** store in the refrigerator). Discard the buffer if crystals are present. Use clean bottles when preparing the buffer. Particles of dust can initiate crystallization.

Make sure to heat the A-AEB in time in the 60 °C (140 °F) water bath (the 100 °C / 212 °F water bath is used for the extraction of the samples).

The **Extractor 2** is provided as 2fold concentrate and has to be diluted 1:2 (1+1) with dist. water (e.g. 30 ml Extractor 2 + 30 ml dist. water). The complete diluted Extractor 2 is sufficient for 15 samples and can be used for approx. 3 month at room temperature (20 - 25 °C / 68 - 77 °F). Additional Extractor 2 can be ordered at R-Biopharm with the article number R4613.

In the following section, the following abbreviations are used:

- AEB: final diluted Allergen extraction buffer
- A-AEB: AEB with addition of Additive 1

## 9.1. Extraction with Extractor 2 and A-AEB as universal sample extraction

### **Heat the A-AEB to 60 °C (140 °F) before sample extraction.**

- Homogenize (grind thoroughly to powder and mix well or mix well a solution respectively) well a sufficient amount (e.g. 50 g or 50 ml) to ensure taking a representative test portion of sample.
- Transfer 1 g (or in case of liquid samples 1 ml) of homogenized sample to a new vial, add 4 ml prepared Extractor 2 (see chapter 9.), close the vial and mix well.
- Cook it for 10 min at 100 °C (212 °F) in a water bath (if a lump forms, it must be dissolved in the still warm liquid by shaking the vial).
- Let the sample cool down shortly on the lab bench (1 - 2 min).
- Add 16 ml (or 15 ml in case of liquid samples) pre-heated A-AEB (see chapter 9.) to the cooked sample.
- Mix vigorously (e.g. vortexer).
- Extract for 10 min at 60 °C (140 °F) in a water bath.
- Let the sample cool down shortly in ice water (3 - 5 min).
- Filter sample or centrifuge for 10 min at > 2,500 x g; if possible at 4 °C (39 °F).  
(Alternatively: transfer 2 ml of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 10 min in a microcentrifuge).
- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The extract (supernatant of centrifugation step or filtrate) can be stored undiluted in a well-sealed container at 2 - 8 °C until used in the test (shelf life approx. 3 days). Unused extracts can also be stored undiluted at -20 °C (- 4 °F) for several months.
- For use in the test, extracts prepared with Extractor 2 (supernatant of centrifugation step or filtrate) must be diluted with AEB (for dilution see chapter 10.2.). The diluted sample extracts have a limited shelf life and must be used in the test within 30 minutes.
- If further dilutions are required (samples with absorbance values ( $A_{450\text{ nm}}$ ) > Standard 5), these extract dilutions must be made BEFORE dilution with AEB for test use (see chapter 10.2.).
- For these dilutions, the following buffer should be used to keep the composition of the extract constant:

A-AEB	16 ml
Extractor 2	2 ml
Dest. water	2 ml

This is followed by the normal dilution with AEB (see chapter 10.2.).

## Remark

**The extracts prepared with Extractor 2 and A-AEB according to chapter 9.1. can be also used with the RIDASCREEN®FAST  $\beta$ -Lactoglobulin (Art. No. R4912) and RIDASCREEN®FAST Milk (Art. No. R4652) assays.**

9.2. Extraction with AEB for food like non-hydrolyzed infant formula, ice cream, wine, chocolate, beverages, sausage, rice crispies

*The following extraction procedure is only suited for the mentioned food since they can also be prepared without Extractor 2, Additive 1 and cooking step.*

### Heat the AEB to 60 °C (140 °F) before sample extraction

- Homogenize (grind thoroughly to powder and mix well or mix well a solution respectively) well a sufficient amount (e.g. 50 g or 50 ml) to ensure taking a representative test portion of sample.
- Transfer 1 g (or in case of liquid samples 1 ml) homogenized sample to a new vial and add 20 ml (or 19 ml in case of liquid samples) pre-heated AEB (see chapter 9.).

*1 ml wine can also be extracted with 9 ml diluted and heated AEB. Doing this increases the method's sensitivity (LOD: 0.12 mg/l; LOQ: 0.25 mg/l).*

- Mix vigorously (e.g. vortexer).
- Extract for 10 min at 60 °C (140 °F) in a water bath by shaking occasionally.
- Let the sample cool down shortly in ice water (3 - 5 min).
- Filter sample or centrifuge for 10 min at > 2,500 g; if possible at 4 °C (39 °F).

(Alternatively: transfer 2 ml of the extract into a reaction vial and centrifuge at high speed (> 10,000 g) for 10 min in a microcentrifuge).

- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- Use the sample extracts (supernatant of centrifugation step or filtrate) immediately (within 30 minutes) in the assay. A longer time period may influence the recovery. Alternatively, the extracts can be stored in a well-sealed vial at 2 - 8 °C until it is used (shelf life approx. 1 day). Unused extracts can be stored undiluted at -20 °C (-4 °F) for several months.

## 10. Test implementation

### 10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **conjugate** (bottle with red cap) is provided as a 11fold concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is needed should be diluted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in conjugate buffer (e.g. 2 ml conjugate buffer + 200 µl conjugate, sufficient for 2 microtiter strips).

The **wash buffer** is provided as a 10-fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with dist. water (i.e. 900 ml dist. water + 100 ml buffer concentrate). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for 4 weeks.

Components should be stored immediately at 2 - 8 °C (35 - 46 °F) when no longer required.

### 10.2. Test procedure

**The extracts prepared according to chapter 9.1. (with Extractor 2 and A-AEB) must be diluted 1:5 (1+4) with AEB (see chapter 9. Sample preparation) before use (e.g. 400 µl AEB + 100 µl sample).**

**Individual samples extracted with Extractor 2 and A-AEB (e.g. corn or corn products, pine, sunflower or pumpkin seeds) may produce non-specific matrix effects in the ELISA. To avoid these, such sample extracts can be diluted 1:5 with BSA-AEB (see chapter 9.) before use in the test (e.g. 400 µl BSA-AEB + 100 µl sample). For this purpose, the diluted AEB is mixed with BSA (final concentration 2.5 %): e.g. 10 ml diluted AEB + 0.25 g BSA.**

**Use the diluted sample extracts immediately (within 30 minutes) in the assay. A longer time period may influence the recovery.**

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at a time. If more than three strips are needed, a second uncoated plate (see chapter 5.1.) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then exactly 100 µl are quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard or sample (prepared according to chapter 9. Sample preparation and - if necessary - diluted according to chapter 10.2. Test procedure) in duplicate to the wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see chapter 10.1.) and pour out the liquid as before. Repeat three more times (a total of four wash cycles).
4. Add 100 µl of the diluted conjugate (see chapter 10.1.) to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times in a row) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see chapter 10.1.) and pour out the liquid as before. Repeat three more times (a total of four wash cycles).
6. Add 100 µl of substrate/chromogen to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 min after addition of stop solution.

## 11. Evaluation

A special software, **RIDASOFT® Win.NET (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done using the 4-parameter or cubic spline function. Once an evaluation method has been selected, it should be retained and not be switched between the two functions. Due to its mathematical background, the cubic spline function cannot calculate concentrations outside the measuring range (< standard 2 and > standard 5). The 4-parameter function enables the calculation of values between standard 1 and standard 2 (refer to chapter 13. Limits of the method).

For the evaluation it should be clarified, that the quality criteria are fulfilled for the current test run. The course of the standard curve can be taken from the attached analysis certificate.

The assay can be also evaluated when running in single determinations. This has no influence on the function of the test kit. A special assay evaluation must be written in the RIDASOFT® Win.NET software for this purpose. It is not present by default. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. However, it is not consistent with standards like EN 15633-1 and EN 15842. It should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

### **Extraction with Extractor 2 and A-AEB (see chapter 9.1.):**

When working according with this extraction, the sample dilution factor is 100. A dilution factor of 20 is already taken into account with the standard concentrations (see chapter 4\*). Hence, values read from the standard curve need to be multiplied with factor 5 only.

### **Extraction with AEB (see chapter 9.2.):**

When working in accordance with this extraction, the sample dilution factor is 20. Hence, a dilution factor of 20 is already taken into account with the standard concentrations (see chapter 4.\*), the casein concentration can directly be read from the standard curve.

If wine is extracted with 9 ml AEB only, a lower dilution of 1:10 needs to be respected when calculating the casein concentration. Values read from the standard curve need to be multiplied with factor 0.5 in this case.

## 12. Result interpretation

The test result is calculated in mg casein per kg food and thus indicates the protein concentration.

Results between LOD and LOQ indicate a low allergen concentration in the sample. Calculated results show a high uncertainty in this area due to the method's high variation below LOQ. Therefore, such results should not be reported with a quantitative value, but qualitative as "< LOQ".

A result below the LOD does not exclude an allergen contamination below the detection limit of the assay, or that other allergen components, such as lipids, may be present in a sample. The result should be reported accordingly.

A further dilution and new detection of samples is recommended for absorbance values ( $A_{450\text{ nm}}$ ) > standard 5. In case of a further dilution, the additional dilution factor must be taken into account when calculating the allergen concentration.

Compared to the certificate, higher absorbance values ( $A_{450\text{ nm}}$ ) for the standard curve, especially for the zero standard, may be a result of insufficient washing or allergen contamination.

## 13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

Technical limits of the test method are approached outside the designated measurement range resulting in higher variation. This may cause a switch of results between the different areas of the calibration curve especially at the test characteristic boundaries (LOD, LOQ, upper limit of measurement range).

Raw fish strongly binds casein. Therefore, the recovery may only be 10 %. This is also true for cooked fish if casein is added before cooking.

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max  $\pm 1$  %.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and inter-laboratory comparisons.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

The addition of foreign protein (depending on the test e.g. BSA, gelatine, skim milk powder) during extraction or test procedure may suppress matrix effects.

The recommendation for the use of BSA mentioned under chapter 10.2. Test procedure should only be applied if a matrix effect is actually present. It is also not limited to the samples mentioned, but is to be understood as a basic recommendation for samples with matrix effect.

In processed foods (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented and this may have an impact on the recovery and assay results.

Cross reactivities are side reactions of the antibody used for preparing the test kit with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. Therefore, cross reactivities may also appear after food processing in single case or are lost.

For evaluation of the cross reactivity only one representative sample was analyzed, other samples may show a different result. All analyzed cross reactivities are described in the validation report.

## 14. Recommendation

In order to ensure a high analytical performance we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in standards like EN 15633-1 and EN 15842 (e.g. performing duplicate determinations).
- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- Carry along test controls for quality control and to ensure an accurate and correct test procedure. Casein-free and casein-containing (naturally contaminated or spiked) samples should be used. An example of a spike experiment is given in the validation report.
- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction may be necessary.
- To contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de) if automates (e.g. ThunderBolt® / GEMINI) are used.

## 15. Further application notes

- Sample extraction for heavily liquid absorbing matrices using RIDA® Extractor 2 (Art. No. R4613).
- More sensitive allergen analysis for wine (dilution 1:10).
- RIDASCREEN®FAST Allergen - Swabbing method for the qualitative analysis of allergens in the production line or for laboratory equipment.
- Sample extraction of pine nuts, sunflower and pumpkin seeds.
- Preparation of food samples with the RIDA® Extraction Solution (colorless) (Art. No. R7098).

*This sample preparation shows similar results as the method with A-AEB and Extractor 2. However, the so prepared extracts cannot be used with the RIDASCREEN®FAST  $\beta$ -Lactoglobulin (Art. No. R4912) and in the RIDASCREEN®FAST Milk (Art. No. R4652).*

**Further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

## Version overview

Version number	Chapter and title
2019-05-03	Release version
2021-06-16	Current version Changes made: <ul style="list-style-type: none"><li>– General linguistic revision</li><li>– Adaptations in chapter 1. Intended use</li><li>– Modification of the disposal clause in chapter 6.</li><li>– Changes in chapter 9: Addition on further dilution for results of <math>A_{450\text{ nm}} &gt; \text{standard 5}</math> for extracted samples (chapter 9.1.)</li><li>– Simplified sample processing: extraction in chapter 9.2. removed and included in chapter 10.2. via special dilution</li><li>– Revision of chapter 11., 12., 13. and 14.</li><li>– Addition in chapter 11. Evaluation (moved from chapter 14.)</li></ul>

## Explanation of symbols

General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM)
	Manufacturer + address

## Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.