RIDASCREEN[®] Aflatoxin Total

Art. No. R4701

Enzymimmunoassay zur quantitativen Bestimmung von
Aflatoxin

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RIDASCREEN[®] Aflatoxin Total

Brief information

RIDASCREEN[®] Aflatoxin Total (Art. No. R4701) is a competitive enzyme immunoassay for the quantitative determination of aflatoxins in cereals and feed. All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: grinding, extraction, filtration/centrifugation and dilution

Time requirement: sample preparation (for 10 samples)approx. 30 min
test implementation (incubation time)45 min

Detection limit: cereals and feed: approx. 1.75 ppb
(corresponding to the standard substance)

Recovery rate: approx. 85 %
(corresponding to the standard substance) coefficient of variation: 15 %

Cross reactivity: Aflatoxin B₁ 100 %
Aflatoxin B₂ approx. 48 %
Aflatoxin G₁ approx. 75 %
Aflatoxin G₂ approx. 18 %

The specificity of the RIDASCREEN[®] Aflatoxin Total test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the limit of detection and the recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our good ELISA practice (GEP) – manual in the respective version. Such lists minimum standards concerning the framework conditions when using test kits of R-Biopharm AG and performing ELISA analysis. The manual can be retrieved, printed and downloaded from the website www.r-biopharm.com/products/food-feed-analysis.

1. Intended use

RIDASCREEN® Aflatoxin Total is a competitive enzyme immunoassay for the quantitative determination of aflatoxins in cereals and feed.

2. General

Aflatoxins are secondary metabolites of the fungi species *Aspergillus flavus*, *parasiticus* und *nomius*. These fungi occur in humid tropical areas and the contamination of vegetable food takes place in the cultivable countries. Aflatoxins belong to the strongest natural occurring cancerogenic substances.

Aflatoxin B₁ which is mostly found together with the aflatoxins B₂, G₁ and G₂ is the one with the highest toxic importance. It is found above all in corn, peanuts, brazil nuts, cotton seed and pistachios.

Due to the toxicity of these mycotoxins maximum levels for aflatoxin B₁ and total aflatoxins for food and feed apply in EU countries.

3. Test principle

The basis of the test is the antigen-antibody reaction. The wells in the microtiter strips are coated with capture antibodies directed against anti-aflatoxin antibodies. Standards or the sample solutions, aflatoxin-enzyme conjugate and anti-aflatoxin antibodies are added. Free and enzyme conjugated aflatoxin compete for the aflatoxin antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-aflatoxin antibodies are also bound by the immobilized capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen solution is added to the wells and incubated. Bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is performed photometrically at 450 nm; the absorption is inversely proportional to the aflatoxin concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Content
Microtiter plate M	-	Ready to use		96 wells
Standard 1*	White	Ready to use	0 µg/l	1.3 ml
Standard 2*	White	Ready to use	0,05 µg/l	1.3 ml
Standard 3*	White	Ready to use	0,15 µg/l	1.3 ml
Standard 4*	White	Ready to use	0,45 µg/l	1.3 ml
Standard 5*	White	Ready to use	1,35 µg/l	1.3 ml
Standard 6*	White	Ready to use	4,05 µg/l	1.3 ml
Wash buffer salt Tween		Dissolve the salt		
Conjugate	Red	Ready to use		6 ml
Antibody	Black	Ready to use		6 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 ml
Stop solution	Yellow	Ready to use		14 ml

5. Materials required but not provided

5.1. Equipment

- grinder (mill)
- graduated cylinder (plastic or glass) 100 ml
- graduated pipettes
- filter funnel and 50 ml flask
- filter paper: Whatman No. 1 or equivalent
- variable 20 µl - 200 µl- and 200 - 1000 µl micropipettes
- microtiter plate spectrophotometer (450 nm)
- optional: shaker, centrifuge

5.2. Reagents

- 70 % methanol solution: mix 70 ml methanol (100 %) with 30 ml distilled water
- distilled water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The standard solutions contain aflatoxin B₁, particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

Decontamination of the glassware and aflatoxin solutions is best carried out using a sodium hypochlorite (bleach) solution (10 %; v/v) overnight (adjust solution with HCl to pH 7).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

Aflatoxins are light-sensitive, therefore, avoid exposure to direct light.

The red stained substrate/chromogen solution is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen solution prior to test implementation
- a value of less than 0.6 absorbance units ($A_{450\text{ nm}} < 0.6$) for the zero standard

9. Preparation of Samples

The samples should be stored in a cool place, protected from light.

Bring all reagents and samples to room temperature (20 - 25 °C / 68 - 77 °F) before use and perform the sample preparation at room temperature.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- weigh 2 g of the ground and homogenized sample into a suitable container and add 10 ml of 70 % methanol*)
- mix for 10 min at room temperature using a shaker
- filter the extract through a Whatman No. 1 filter (or equivalent) or centrifuge (10 min / 3500 g / room temperature)
- dilute 100 µl of the filtrate with 600 µl distilled water
- use 50 µl of the diluted filtrate per well in the test

*) sample size may be increased if required, but the volume of methanol/water must be adapted accordingly, e.g.: 10 g in 50 ml of 70 % methanol

Remark:

If the aflatoxin concentration is expected to exceed 120 ppb further dilutions are necessary. For this use distilled water containing 10 % methanol (i.e. 9 ml distilled water + 1 ml methanol (100 %)). Please note that any kind of sample being used in the assay has to be provided in distilled water with 10 % methanol.

For the analysis of nut, herb, spice or tea leave samples please perform a sample preparation using immunoaffinity columns. For this please request the RIDA[®] Aflatoxin column (R5001 / R5002) product information.

An application note for green coffee in combination with RIDA[®] Aflatoxin column (R5001 / R5002) is available on request. Please contact your local distributor.

Please note, that for the sample preparations with RIDA[®] Aflatoxin column (R5001 / R5002) the respective eluates have to be diluted 1:10 with distilled water (see section 9.1. of the instructions for use of RIDA[®] Aflatoxin column (R5001 / R5002)) . For all further dilutions the use of distilled water with 10 % methanol is absolutely essential.

10. Test implementation

10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use and perform the test at room temperature.

Return all reagents to 2 - 8 °C (35 - 46 °F) immediately after use.

10.2. Wash buffer

As **wash buffer** a PBS tween buffer is needed. Please use the wash buffer salt contained in the kit (see 4.). Dissolve the entire buffer salt in one liter of distilled water. The ready to use wash buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternative: Dissolve the contents of the envelope in 100 ml of distilled water to obtain a 10fold concentrated wash buffer. This solution expires after approx. 8 - 12 weeks, store at room temperature (20 - 25 °C / 68 - 77 °F).

Use 1 part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use wash buffer.

10.3. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry up and avoid prolonged intervals between the working steps.

Avoid direct sunlight during all incubation steps and incubate plates in the dark.

1. Insert a sufficient number of microtiter wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Add 50 µl of the standard or prepared sample to separate duplicate wells.
3. Add 50 µl of the conjugate to each well.
4. Add 50 µl of the antibody to each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously against absorbent paper (three times in a row) to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.2.). Empty the wells again. Repeat two more times.
6. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C; 68 - 77 °F).
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 minutes after addition of stop solution.

11. Results

A specific software, the RIDA[®]SOFT Win.net (Art. No. Z9996), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays.

For sample analysis we recommend to perform double or multiple determinations. Please use the RIDA[®]SOFT Win.net cubic spline function for evaluation of double or multiple determinations.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the aflatoxin concentration [ng/kg].

In order to obtain the aflatoxin concentration in ng/kg actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factor is as follows:

cereals and feed 35

For further information and applications please contact your local distributor of R-Biopharm AG (sales@r-biopharm.de).

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG