r-biopharm®

RIDASCREEN® FAST Ochratoxin A

Art. No. R5402

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RIDASCREEN®FAST Ochratoxin A

Brief information

RIDASCREEN®FAST Ochratoxin A (Art. No. R5402) is a competitive enzyme immunoassay for the quantitative analysis of ochratoxin A in corn, wheat, barley, oats (grains), and feed.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification. Please note that the following method must be used for the evaluation: Ridawin.NET \ FOOD \ Mycotoxins FAST \ R5402 FAST Ochratoxin A ECO.met.

Sample preparation: extraction, centrifugation, dilution

Time requirement: sample preparation (for 10 samples)
Corn, wheat, barley, oatsapprox. 20 min
Feed.....approx. 20 min
Test implementation (incubation time).....8 min

Limit of detection: Corn 1.3 µg/kg (ppb)
Wheat..... 1.5 µg/kg (ppb)
Barley..... 1.5 µg/kg (ppb)
Oats 2.0 µg/kg (ppb)
Feed..... 2.8 µg/kg (ppb)

Recovery rate: Corn 100 +/- 30 %
(naturally contaminated samples - TRILOGY®) Wheat..... 100 +/- 30 %
Feed..... 80 -130 %

Recovery rate: Corn (FAPAS) 60 - 130 %
(spiked samples) Wheat (FAPAS)..... 60 - 100 %
Barley (FAPAS)..... 60 - 100 %
Oats (grains) 60 - 130 %
Feed (FAPAS)..... 100 +/- 30 %

Note: The assay was checked using naturally contaminated samples. Deviations in the recovery of spiked samples are possible.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

RIDASCREEN® Ochratoxin A (R1312)

RIDA® Ochratoxin A column (R1303)

TRILOGY® Certified Liquid Standard Ochratoxin A (CTSL-520-5)

TRILOGY® Liquid Standard Ochratoxin A (TSL-504-5)

TRILOGY® Dried Standard Ochratoxin A (TS-503-5)

1. Intended use

RIDASCREEN®FAST Ochratoxin A is a competitive enzyme immunoassay for the quantitative analysis of ochratoxin A in corn, wheat, barley, oats (grains), and feed.

2. General

The mycotoxin ochratoxin A is formed by fungi of the species *Aspergillus* and *Penicillium*. Apart from a marked nephrotoxicity, ochratoxin A displays hepatotoxic, teratogenic, carcinogenic and immunosuppressive properties.

There is a risk to human health not only through the intake of contaminated foods of vegetable origin, but also through foods of animal origin. Ochratoxin A has been detected in pig blood and kidneys, as well as in human blood and mother's milk.

3. Test principle

The basis of the enzyme immunoassay is the antigen-antibody reaction. The wells in the microtiter strips are coated with specific antibodies against ochratoxin A. Ochratoxin A standards or the sample solutions and enzyme conjugate are added. Free and enzyme-conjugated ochratoxin A compete for the ochratoxin A antibody binding sites (competitive enzyme immunoassay). Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells and incubated. Bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the ochratoxin A concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for as many as 42 analyses (plus 6 standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
ECO extractor	transparent	Concentrate 10fold		2 x 120 ml
Standard 1*	white	Ready to use	0 µg/L	1.3 ml
Standard 2*	white	Ready to use	1 µg/L	1.3 ml
Standard 3*	white	Ready to use	3 µg/L	1.3 ml
Standard 4*	white	Ready to use	10 µg/L	1.3 ml
Standard 5*	white	Ready to use	30 µg/L	1.3 ml
Standard 6*	white	Ready to use	100 µg/L	1.3 ml
Wash buffer salt Tween**		Dissolve the salt		
Conjugate	red	Ready to use		3 ml
Substrate/Chromogen Red Chromogen Pro	brown	Ready to use		10 ml
Stop solution	yellow	Ready to use		14 ml

- *) The dilution factor 10 for the sample preparation has already been considered. Therefore, the ochratoxin A concentrations (µg/kg or µg/L) of samples can be read directly from the standard curve.
- ***) The ready-to-use wash buffer (see preparation 10.1) is used as sample and wash buffer.

5. Materials required but not provided

5.1. Equipment

- laboratory or grain mill
- graduated cylinder (plastic or glass) 100 ml
- centrifuge ($\geq 3,500$ g)
- microtiter plate spectrophotometer (450 nm)
- (horizontal) shaker
- Vortex mixer
- variable 50 μ l - 100 μ l and 1000 μ l micropipettes

5.2. Reagents

- distilled or deionized water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The standards contain ochratoxin A. Particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

Decontamination of the glassware and ochratoxin A solutions is best carried out using a sodium hypochlorite (bleach) solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen solution is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

The test kit can be regularly used at least up to the expiry date (indicated on the kit package), if stored correctly.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Any bluish coloration of the reddish substrate/chromogen solution prior to test implementation
- A value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for the zero standard

9. Preparation of Samples

The samples should be stored in a cool place, protected against light.

9.1. Extraction buffer

For extraction, the diluted ECO extractor is needed. To obtain the ready-to-use ECO extractor, dilute the ECO Extractor (10x Concentrate) 1:10 with distilled or deionised water. The diluted ECO Extractor is stable for one week at 2 - 8 ° C (36 – 46 °F). If turbidity occurs in the diluted ECO extractor (possibly caused by contamination), discard or do not use it.

9.2. Extraction of corn, wheat, barley, oats (grains) and feed

Bring all reagents and samples to room temperature (20 – 25 °C / 68 – 77 °F) before use and perform the sample preparation at room temperature.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure (recommended particle size: 500 µm).

We recommend to produce an additional replica per unknown sample.

- Weigh 10 g of ground and homogenized sample into a suitable container (e.g. 125 ml bottle) and add 50 ml of diluted ECO extractor
- Vortex the sample briefly (10 seconds)
- Shake the sample vigorously for 5 minutes (manually or with shaker at 420 rpm)
- Centrifuge for 5 min at 3,500 g and room temperature (20 – 25 °C; 36 – 46 °F)
- Dilute 1 ml of the supernatant with 1 ml of ready-to-use wash buffer (sample buffer, see 10.1.)
- Add 50 µl of the diluted supernatant per well in the assay

Note: Samples measured outside the measurement range of > 100 µg/kg (ppb) are recommended to be further diluted 1:10 or maximum 1:20 with ready-to-use wash buffer (see 10.1.). The resulting dilution factor of 10 (1:10 dilution) or 20 (1:20 dilution) should be taken into account in the calculation (i.e. the results must still be multiplied by 10 and 20, respectively).

10. Test implementation

10.1. Preliminary comments

1. Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use and perform the test at room temperature.
2. Return all reagents to 2 - 8 °C (35 - 46 °F) immediately after use.

As **wash buffer** and sample buffer a PBS tween buffer is needed. Please use the buffer salt contained in the kit (see 4.). Dissolve the entire buffer salt in one liter of distilled water. The ready-to-use wash buffer expires after approx. 4 - 6 weeks when stored at 2 - 8 °C (36 - 46 °F).

Alternative: Dissolve the contents of the envelope in 100 ml of distilled water to obtain a 10fold concentrated washing buffer. The 10fold concentrate expires after approx. 8 - 12 weeks when stored at room temperature (20 - 25 °C / 68 - 77 °F).

10.2. Test procedure

Accurate washing is very important. Do not allow microwells to dry up completely and avoid prolonged intervals between the working steps. Therefore, the processing should be done quickly. Carefully follow the recommended washing procedure as outlined in the test procedure. Reproducibility in any enzyme immunoassay is largely dependent upon the consistency with which the microwells are washed. It may not be processed more than 3 strips per test run.

Avoid direct sunlight during all incubations. Therefore, covering the microtiter plates is recommended.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Pipette 50 µl of standard or prepared sample into separate wells; use a new pipette tip for each standard or sample.
3. Add 50 µl of conjugate to each well. Mix gently by shaking the plate manually and incubate for 5 min at room temperature (20 – 25 °C / 68 – 77 °F) in the dark.
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat the washing procedure two times.
5. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 3 min at room temperature (20 - 25 °C / 68 – 77 °F) in the dark.
6. Add 100 µl of stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 15 minutes after addition of stop solution.

11. Results

A specific software, the RIDASOFT® Win.NET (Art. No. Z9996), is available to evaluate the RIDASCREEN® enzyme immunoassays.

Please note that the following method Ridawin.NET \ FOOD \ Mycotoxins FAST \ R5402 FAST Ochratoxin A ECO.met must be used for the evaluation.

The evaluation is done with logit/log.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a

system of coordinates on semilogarithmic graph paper against the ochratoxin A concentration [$\mu\text{g}/\text{kg}$].

The ochratoxin A concentration in $\mu\text{g}/\text{kg}$ corresponding to the extinction of each sample can be read from the calibration curve.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG