RIDASCREEN[®] FAST Mandel / Almond

Art. No. R6901

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RIDASCREEN®FAST Mandel / Almond

Brief information

RIDASCREEN®FAST Mandel / Almond (Art. No. R6901) is a sandwich enzyme immunoassay for the quantitative analysis of almond or parts of almond in food. As representative samples for the food groups pastries and sweets the following matrices have been tested during the assay's development: cereals, cookies, ice cream and chocolate. The assay can be likely used for analysis of other food samples too; this has to be checked by the user.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit. The test kit is sufficient for 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

- Sample preparation: homogenization, extraction and centrifugation
- Time requirement: sample preparation (for 10 samples) approx. 20 min
test implementation (incubation time) 30 min
- Limit of detection: 0.1 (0 – 0.23) mg/kg (ppm) almond
- Limit of quantification: 2.5 mg/kg (ppm) almond
- Specificity: The antibodies specifically detect proteins from almonds.
Cross reaction to apricot stone, mulberries, mahaleb cherry (cross reactions to other plants from genus Prunus are possible).

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments. Further information is described in the updated validation report.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in

the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded from the website <http://www.r-biopharm.com/products/food-feed-analysis>.

Related products

Bioavid Lateral Flow Mandel/Almond (BL601-10 und -25)
SureFood® ALLERGEN Almond (S3604)

1. Intended use

RIDASCREEN®FAST Mandel / Almond is a sandwich enzyme immunoassay for the quantitative analysis of almond or parts of almond in food. As representative samples for the food groups pastries and sweets the following matrices have been tested during the assay's development: cereals, cookies, ice cream and chocolate. The assay can be likely used for analysis of other food samples too; this has to be checked by the user. Limit of detection and limit of quantification depend on the tested food matrix, the level of food processing, and the extraction method. Detailed results and further validation information about other sample matrices can be found in the brochure Product Information. Additionally, results from proficiency tests and inter laboratory trials may exist for different food matrices. Further applications are frequently validated in our laboratory and provided in the Application Notes.

2. General

The allergen can be present as an ingredient or as a contamination in raw and cooked products. According to the **regulation (EU) No. 1169/2011**, almond must be declared on food labels as it can induce allergic reactions. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

3. Test principle

The wells of the microtiter strips are coated with specific antibodies to almond proteins. By adding standards and samples to the wells, almond protein present will bind to the specific antibodies. In a washing step components not bound are removed. Then antibody conjugated to peroxidase is added. This antibody conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. The detection of almond protein takes place by adding substrate/chromogen. The enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the almond proteins concentration of the sample. The result is expressed in mg/kg almond.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard measurements).

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use	-	48 wells
Allergen extraction buffer	Green	Concentrate	10x	100 ml
Standard 1*	Transparent	Ready to use	0.0 mg/kg	2.6 ml
Standard 2*	Transparent	Ready to use	2.5 mg/kg	2.6 ml
Standard 3*	Transparent	Ready to use	5.0 mg/kg	1.3 ml
Standard 4*	Transparent	Ready to use	10.0 mg/kg	1.3 ml
Standard 5*	Transparent	Ready to use	20.0 mg/kg	1.3 ml
Wash buffer	Brown	Concentrate	10x	100 ml
Conjugate	Red	Concentrate	11x	0.7 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use	-	10 ml
Stop solution	Yellow	Ready to use	-	14 ml

*) The dilution factor **20** for the sample preparation has already been considered when labeling. Therefore, the almond concentrations of samples can directly be read from the standard curve.

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal vials
- shaker
- water bath
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenisator
- graduated pipettes
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes

5.2. Reagents:

- distilled or deionized water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instructions for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen prior to test implementation
- a value of less than 1.2 absorbance units ($A_{450\text{ nm}} < 1.2$) for standard 5

9. Preparation of Samples

Working devices such as a mill, glass vials or spatulas must be cleaned before and after each sample preparation to remove any remains of almond and to avoid contamination.

The **Allergen Extraction buffer** is provided as a **10fold concentrate**. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99°F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted buffer is stable at 20 – 25 °C (68 -77 °F) for approx. four weeks.

9.1. Sample preparation (**without** skim milk powder)

- grind at least 5 g of the sample carefully and mix thoroughly
- weigh 1 g of the sample and add 20 ml diluted Allergen Extraction buffer (the extraction buffer should already have a temperature of approx. 60 °C (140 °F))
- or add 19 ml of the diluted Allergen Extraction buffer to 1 ml liquid sample (the extraction buffer should already have a temperature of approx. 60 °C (140 °F))
- mix intensively and incubate for 10 min at 60 °C (140 °F), afterwards cool down
- centrifuge: 10 min, at least 2500 g, if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl of the supernatant or filtrate per well in the assay

9.2. Sample preparation (**with** skim milk powder) for spices, sorbet and cacao

- at least 5 g of the sample should be ground well and thoroughly mixed
- weigh 1 g of sample and add 1 g skim milk powder
- add 20 ml extraction buffer (the extraction buffer should have already been heated to approx. 60 °C (140 °F))
- mix intensively and extract for 10 min at 60 °C (140 °F) by shaking casually
- centrifuge: 10 min / 2500 g / if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged for 10 minutes at high speed in reaction caps by using a microcentrifuge)
- use 100 µl of the supernatant or filtrate per well in the assay

Remark:

The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 3 days.

Further dilutions should be prepared with the finally diluted Allergen Extraction buffer.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **Conjugate** (bottle with red cap) is provided as an 11fold concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is

needed should be reconstituted. Before diluting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) with distilled water (e. g. 200 µl concentrate + 2 ml distilled water, sufficient for 2 microtiter strips).

The **Wash Buffer** is provided as a 10fold concentrate. Before dilution of the buffer concentrate, dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted wash buffer is stable at 20 - 25 °C (68 - 77 °F) for approx. four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than 3 strips (24 wells) at a time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101 or Mikrotiter Assembly Breakable Strip 1x8, Thermo Scientific) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Add 100 µl of each standard or prepared sample to the wells, mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell plate upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
4. Add 100 µl of the diluted conjugate (see 10.1.) to each well, mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell plate upside down vigorously (three times in a row) against absorbent paper to ensure complete

removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.

6. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 min after addition of stop solution.

11. Results

A special software, the RIDASOFT® Win.NET (Art. Nr. Z9996) is available for evaluation of the RIDASCREEN® enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit. The document 'Compliance Criteria' provides criteria for evaluating standard curves. For quality assurance assay controls should be used.

In comparison with the certificate, higher values of the absorbance ($A_{450\text{ nm}}$) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or allergen contamination.

Please note:

When working according to the described sample preparation, the dilution factor is 20. The allergen concentration can be read directly from the standard curve (the sample dilution factor of 20 is already taken into account for the standard concentrations (see 4. *)).

For sample dilutions of more than 1:20, the further dilution factor must be considered for the calculation of the almond concentration.

In general:

Samples tested negative could still contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

Allergen containing samples that have been heat treated show a reduced recovery because the proteins denature and are no longer recognized by the

antibody. The reduction in recovery depends strongly on the temperature and the duration of the heat treatment. If samples are heat treated at high temperature the recovery can be significantly reduced.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross-reactivities and exemplary analysed matrices are described in the updated validation report.

The protein content and the protein composition may vary considerably between different almond species. Therefore, different varieties may produce different results, since exemplary varieties were used for calibration. The protein content of the standard material is described in the validation report and can be used for the calculation to protein.

Recommendation:

In order to ensure a high analytical performance we recommend:

- in case of extremely acid or alkaline samples, the pH should be adjusted to a neutral pH
- to use allergen-free and allergen containing (spiked) samples as test controls
- due to the multitude of food types, matrix effects cannot be excluded
- to carry out spiking experiments for an accurate and correct procedure
- to analyze each sample material in duplicates
- to perform SureFood[®] PCR to confirm results
- to contact sales@r-biopharm.de if automates (e.g. ThunderBolt[®] / Bolt) are used

Further product information and application notes, please contact sales@r-biopharm.de.

The product information folder with further information is available on request from your local distributor or R-Biopharm AG.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

R-Biopharm AG

Postanschrift / Postal Address:

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64297 Darmstadt, Germany

Sitz / Corporate Seat: Pfungstadt

Tel.: +49 (0) 61 51 - 81 02-0

Fax: +49 (0) 61 51 - 81 02-40

E-mail: info@r-biopharm.de

www.r-biopharm.com

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are registered trademarks of R-Biopharm AG