RIDASCREEN®FAST Aflatoxin SC

REF R9002

Enzymimmunoassay zur quantitativen Bestimmung
von Aflatoxin

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RIDASCREEN® FAST Aflatoxin SC

Brief information

RIDASCREEN®FAST Aflatoxin SC (Art. No.: R9002) is a competitive enzyme immunoassay for the quantitative analysis of aflatoxin in corn, brown rice and feed.

All reagents required for the enzyme immunoassay, including standard, are contained in the test kit. The test kit is sufficient for 48 determinations (including standard). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: extraction, filtration, dilution

Time requirement: sample preparation (for 10 samples)..... approx. 10 min
test implementation (incubation time).....15 min

Limit of detection: Corn1.5 µg/kg (ppb)
Feed5.3 µg/kg (ppb)

Recovery rate: Corn 70 - 120 %
(naturally contaminated samples - Trilogy®) Brown rice70 - 120 %
Feed 70 - 100 %

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice. Such lists minimum standards concerning the framework conditions when using test kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Product offer

RIDASCREEN® Aflatoxin B1 30/15 (R1211)
RIDASCREEN® Aflatoxin Total (R4701)
RIDASCREEN®FAST Aflatoxin (R5202)
EASI-EXTRACT® Aflatoxin (RBRRP71/RBRRP70N)
RIDA® Aflatoxin column (R5001/R5002)

1. Intended use

RIDASCREEN®FAST Aflatoxin SC is a competitive enzyme immunoassay for the quantitative analysis of aflatoxin in corn, brown rice and feed.

2. General

Aflatoxins are secondary metabolites of the fungi species *Aspergillus flavus*, *parasiticus* and *nomius*. These fungi occur in humid tropical areas and the contamination of vegetable food takes place in the cultivable countries. Aflatoxins are highly toxic and belong to the strongest natural occurring carcinogenic substances.

Aflatoxin B₁, which is generally present together with aflatoxin B₂, G₁ and G₂, can be detected in cereals, corn, cotton-seed and some nuts.

3. Test principle

The basis of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with capture antibodies directed against anti-aflatoxin antibodies. Aflatoxin standard or sample solutions, aflatoxin enzyme conjugate and anti-aflatoxin antibodies are added. Free aflatoxin and aflatoxin enzyme conjugate compete for the aflatoxin antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-aflatoxin antibodies are also bound by the immobilized capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the aflatoxin concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for max. 48 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format	Volume
Microtiter plate M	-	ready to use	48 wells
Standard 1^{*)}	white	ready to use	0 µg/L 1.3 ml
Conjugate	red	ready to use	3 ml
Antibody	black	ready to use	3 ml
Substrate/Chromogen Red Chromogen Pro	brown	ready to use	10 ml
Wash buffer salt Tween		Dissolve the salt	
Stop solution	yellow	ready to use	14 ml

*) Only standard 1 (0 µg/L) is included in the test kit. The standard curve (B/B₀) is provided with the QC certificate of the test kit. For the calculation of results see 11. Results.

5. Reagents required but not provided

5.1 Equipment

- Microtiter plate spectrophotometer (450 nm)
- Graduated cylinder (plastic or glass) 100 ml
- Glassware for preparing sample extract: filter funnel and 50 ml flask
- Shaker
- grinder (mill)
- filter paper: Whatman No. 1 or equivalent
- 50 µl, 100 µl and 1000 µl micropipettes

for USDA/GIPSA extraction method:

- ultra-turrax (or equivalent)
- filter syringe (JM 1000 or equivalent)

5.2 Reagents

- Distilled or deionized water
- Methanol
- 70 % methanol solution: prepare 70 % methanol solution by mixing 70 ml methanol (100 %) with 30 ml distilled or deionized water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Ensure the proper and responsible disposal of all reagents and materials after their use. For disposal, please adhere to national regulations.

Decontamination of the glassware and aflatoxin solutions is best carried out using a sodium hypochlorite (bleach) solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Dilution or adulteration of these reagents may result in loss of sensitivity

The test kit can be regularly used at least up to the expiry date (indicated on the kit label), if stored correctly.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the red stained substrate/chromogen prior to test implementation
- a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for the zero standard

9. Preparation of Samples

The samples should be stored in a cool place, protected from light. A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

Bring all reagents and samples to room temperature (20 - 25 °C) before use and perform sample preparation at room temperature.

- weigh 5 g of ground sample into a suitable container and add 25 ml of 70 % methanol *)
- shake vigorously for three minutes (manually or with shaker)
we recommend the extraction by using a shaker or blender
- filter the extract through Whatman No. 1 filter (or equivalent)
- dilute 1 ml of the obtained filtrate with 1 ml of distilled or deionized water
- use 50 µl of the diluted filtrate per well in the test

*) sample size may be increased if required, but the volume of methanol/water must be adapted accordingly, e.g.: 10 g in 50 ml of 70 % methanol

USDA/GIPSA extraction method:

- weigh 50 g of ground sample into a suitable container and add 250 ml of 70 % methanol
- blend the sample by ultra-turrax (or equivalent) for two minutes or shake vigorously for three minutes (manually or with shaker)
- filter the extract using a filter syringe (JM 1000 or equivalent): pour the sample extract into the syringe within 1 cm from the top; using the plunger, push approx. 1.5 ml of the extract through the syringe filter
- collect the filtered extract into a clean tube
- dilute 1 ml of the filtrate with 1 ml of distilled or deionized water
- use 50 µl of the diluted filtrate per well in the test

Remark:

Samples found out of range > 100 µg/kg (ppb) have to be further diluted, e.g. 1:3 (1+2) with methanol/water (35/65), in order to get a correct result. The resulting dilution factor 3 has to be taken into account when calculating the result (the results need to be multiplied by 3).

10. Test implementation

10.1 Test preparation

1. Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.
2. The specific reaction starts with the addition of the specific antibody. Nevertheless, not more than three strips should be applied to the test, when a single step pipette is used. More strips (up to 6) could be applied, when a multistep pipette is used.
3. Return all reagents to 2 - 8 °C (35 - 46 °F) immediately after use.

The **aflatoxin standard 1 (0 µg/L)** is provided ready to use. B/B₀-values of Aflatoxin standards 2 - 7 (2, 4, 10, 20, 50 and 100 ppb) are reported on the QC certificate of the test. The standard curve is calculated with the **RIDASOFT® Win.NET** (Art. No. Z9996FF) (see 11. Results) according to those values. The dilution factor 10 for the sample has been considered in this calculation.

As **washing buffer** a PBS-Tween buffer is needed, please use the washing buffer salt (pouch) contained in the kit (see 4.). Please use the buffer salt contained in the kit (see 4.). Dissolve the entire buffer salt in one liter of distilled water. The ready-to-use wash buffer expires after approx. 4 - 6 weeks when stored at 2 - 8 °C (36 - 46 °F).

Alternatively: Dissolve the contents of the envelope in 100 ml of distilled water to obtain a 10fold concentrated washing buffer. The 10fold concentrate expires after approx. 8 - 12 weeks when stored at room temperature (20 - 25 °C / 68 - 77 °F). Use 1 part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use wash buffer.

10.2 Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry up completely and avoid prolonged intervals between the working steps. Therefore, the processing should be done quickly. Carefully follow the recommended washing procedure as outlined in the test procedure. Reproducibility in any EIA is largely dependent upon the consistency with which the microwells are washed.

Avoid direct sunlight during all incubations. Therefore, covering the microtiter plates is recommended. 1

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for the standard and samples to be run. Record standard and sample positions.
2. Pipet 50 µl of standard or prepared sample into separate wells; use a new pipette tip for the standard or each sample.
3. Add 50 µl of conjugate to the bottom of each well.
4. Add 50 µl of antibody to each well. Mix gently by shaking the plate manually and incubate for 10 min (+/- 1) at room temperature (20 - 25 °C / 68 - 77 °F).
5. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Using a wash bottle or multichannel pipette, fill the wells (250 µl per well) with washing buffer (see 10.1.). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.
6. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 5 min (+/- 0.5) at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes.

11. Results

For calculation of results you need to transfer the measured absorbance for standard 1 (0 ppb) as well as the B/B₀ values for standards 2 - 7 (2, 4, 10, 20, 50, 100 ppb) provided with the test kit QC certificate, into the **RIDASOFT® Win.NET** (Art. No. Z9996FF). From this the software program calculates the standard curve and the content of aflatoxin in the samples.

Please note:

The RIDASOFT® Win.NET (Art. No. Z9996FF) can be ordered from your local distributor.

Further product information, please contact your local distributor or R- Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2005-05-04	General revision
2007-05-18	General revision
2010-01-27	General revision

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.