CONGEN

**SureFood® GMO QUANT
T25 Corn**

Art. No. S2017

2 x 50 rxn

User Manual

December 2018

 Inhalt /  Content

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1 General Information

1.1 Description

The test detects the relative quantitative T25-Corn DNA amount. Therefore the kit contains two PCR systems, one specific for the T25-Corn (OECD unique identifier ACS-ZMØØ3-2) and the other one specific for Corn (reference system). The event-specific detection of T25-Corn is according to the validated method of the European Commission. The real-time PCR assay can be used with established real-time PCR instruments. The technical validation of instruments was performed on Agilent Mx3005P, Agilent Aria Dx, BioRad CFX 96, Roche LightCycler® 480 II, Roche cobas Z 480 Analyzer, Applied Biosystems 7500 and Qiagen Rotor-Gene Q.

1.2 Limit of Detection

The T25 corn has a limit of detection of ≤ 5 DNA-copies. The assay limit of detection depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The limit of quantitation depends on the concentration of the sample DNA used in the analysis. For example, if 50,000 target-sequence copies of Corn specific reference gene are present, the relative quantitation limit for T25 DNA is 0.1 %.

1.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended.

1.4 Kit components and storage

| Kit Code | Reagent | Amount | Lid Color |
|----------|-------------------|----------------------|------------|
| 1 | Corn Reaction Mix | 1 x 1100 µl | Orange |
| 2 | T25 Reaction Mix | 1 x 1100 µl | Yellow |
| 3 | Taq Polymerase | 1 x 11 µl | Red |
| 4 | Dilution Buffer | 1 x 1300 µl | White |
| 5 | Standard DNA | 1 x 50 µl | Dark Blue |
| 6 | Positive Control | 1 x 100 µl / 1 % T25 | Light Blue |

Store all reagents at -20°C and protected from light.

1.5 Additionally required equipment and materials

- real-time PCR instrument
- real-time PCR consumables (plates, tubes, foils, capillaries, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Setup

| | Blockcycler/Bio Molecular Systems MIC | Rotorcycler |
|--|--|---|
| Initial Denaturation (HOLD) | 5 min, 95°C | 1 min, 95°C |
| Cycles | 45 | 45 |
| Denaturation | 15 sec, 95°C | 10 sec, 95°C |
| Annealing/Extension (CYCLE) | 30 sec, 60°C | 15 sec, 60°C |
| Temperature Transition Rate/ Ramp Rate | Maximum | Maximum |
| Fluorescence Detection Setup (exemplary) | Detection: End of extension phase Reporter Dye: FAM Quencher Dye: BHQ Passive Reference: none | Rotorgene Q Reporter Dye: FAM (Green) |
| Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: http://www.congen.de/en/company/downloads | | |

2 Quantitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples, controls and standards).

Reactions needed for the corn detection:

For each run: 5 reactions for the standard curve
3 reactions for controls (1x no-template and 2x Positive Control)

For each sample: at least 1 reaction with each sample DNA

Reactions needed for the T25 corn detection:

For each run: 5 reactions for the standard curve
3 reactions for controls (1x no-template and 2x Positive Control)

For each sample: at least 1 reaction with each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

| Components for master-mix | Amount per reaction | 10 reactions (with 10% excess) |
|--|----------------------------|---------------------------------------|
| Corn Reaction Mix or T25 Reaction Mix | 19.9 µl | 218.9 µl |
| Taq Polymerase | 0.1 µl | 1.1 µl |
| Total volume | 20 µl | 220.0 µl |

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the master-mix

| Components for master-mix | Amount per reaction | 10 reactions (with 10% excess) |
|---------------------------|---------------------|--------------------------------|
| Reaction Mix | 19.9 µl | 218.9 µl |
| Taq Polymerase | 0.1 µl | 1.1 µl |
| Total volume | 20 µl | 220 µl |

Mix each master-mix well and centrifuge shortly before use.

2.1.3 Preparation of the standard DNA dilutions

Dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different DNA concentrations for the standard curves of the Corn reference gene and the T25 Corn detection gene. Prepare 5 dilutions of the supplied Standard DNA with the supplied Dilution Buffer. Prepare 5 reaction tubes (labeled S1 to S5) and add 45 µl Dilution Buffer (**Code 4**) each. The following procedure is recommended:

| standard | dilution | copy number per µl | final copy number per reaction* |
|----------|---|--------------------|---------------------------------|
| S1 | 45 µl Dilution Buffer + 5 µl Standard DNA | 100,000 copies | 500,000 copies |
| S2 | 45 µl Dilution Buffer + 5 µl DNA of S1 | 10,000 copies | 50,000 copies |
| S3 | 45 µl Dilution Buffer + 5 µl DNA of S2 | 1000 copies | 5000 copies |
| S4 | 45 µl Dilution Buffer + 5 µl DNA of S3 | 100 copies | 500 copies |
| S5 | 45 µl Dilution Buffer + 5 µl DNA of S4 | 10 copies | 50 copies |

***Note:** 5 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

If the diluted DNA standards (S1 to S5) are not immediately used, store them at -20°C. Before use allow the reagents to thaw, mix them on a vortex and centrifuge carefully before opening and use.

2.1.4 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries and close them.
- Pipette 5 µl of the Positive Control and the standard dilutions into the designated tubes/wells or capillaries and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates or capillaries into the PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The calculation for both reaction systems (**Corn**, **T25 Corn**) has to be made separately. Mark the standards, the controls and the samples for the specific system (**T25 Corn**) and make the evaluation according to the usual analysis program recommended by the real-time PCR instrument manufacturer. Repeat the same procedure for the **Corn** reference gene system. The value for the slope of the standard curve has to be between -3.1 and -3.6 and the correlation coefficient $R^2 > 0.98$. In case of different values for the standard curve, it should not be used for calculation.

By using the calculated copy numbers for **T25 Corn** and **Corn** the relative GMO content of the sample DNA and the Positive Control can be determined in the following way (example):

| | | | |
|------------------------|---------------|----------------------------------|---------------|
| Sample T25 Corn | 1350 copies | Positive Control T25 Corn | 400 copies |
| Sample Corn | 45,000 copies | Positive Control Corn | 28,000 copies |

Multiply the copy number of the specific system by 100 and divide by the copy number of the reference gene system to obtain the percentage.

T25 Corn content = **T25 Corn** copy number * 100 / **Corn** copy number

sample DNA **T25 Corn** content = $1350 * 100 / 45,000$ sample DNA **T25 Corn** content = 3 %

For the given example the numbers lead to a **T25 Corn** content of 3.0 % for the sample, and 1.4 % for the Positive Control with the same calculation.

For a final calculation the use of a correction factor K for the correction of run-to-run fluctuations is necessary. The correction factor is the relation of the true percentage value of the Positive Control (1 % GMO content) and the measured GMO percentage of the Positive Control. The factor is calculated in the following way:

$K = \text{true GMO percentage of Positive Control} / \text{measured GMO percentage of Positive Control}$

$K (\text{example}) = 1 \% / 1.4 \% = 0.7$

The calculated value for the sample is multiplied with K to obtain a corrected value.

$\text{GMO percentage sample} = \text{measured GMO percentage of sample} * K$

$\text{sample (example)} = 3.0 \% * 0.7 = 2.1 \%$

For that example the **T25 Corn** content is 2.1 %.

December 2018

3 Further Information

3.1 Product Information

- Validation Report
- Microsoft Excel template of calculation
(Download: www.congen.de/en/company/downloads)

3.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

3.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com