**SureFood® GMO SCREEN 4plex 35S/NOS/FMV+IAC (100 Reakt.)**

Art. Nr. S2126

Dezember 2016

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Description

This kit can be used for screening of genetically modified organisms (GMOs) in food, feed and seeds. For this purpose, PCR systems for detection of the 35S Cauliflower Mosaic Virus (CaMV) promoter DNA sequence, *A. tumefaciens* NOS terminator DNA sequence and for detection of the 34S FMV promoter DNA sequence are applied. The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of four fluorescence emissions at 522 nm, 553 nm, 610 nm and 670 nm (FAM, VIC, ROX and Cy5) at the same time (Rotor-Gene Q, Agilent MxSeries, BioRad CFX96, ABI7500, Roche LightCycler® 480² etc.). The kit is according to German Food Law § 64 for the detection of genetically modified DNA sequences.

Additionally the kit contains an internal amplification control (inhibition control, IAC) to examine DNA for inhibiting substances which interfere with the PCR mix. The Reaction Mix contains all reagents necessary for PCR as well as an artificial target (concentration: 100 copies / reaction). The target is amplified during the PCR run and detected in the VIC/HEX-channel. When the DNA contains PCR inhibiting substances, amplification of the target will be suppressed and the signal will be affected.

Some examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). Also spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

Limit of Detection

The SureFood® GMO SCREEN 4plex 35S/NOS/FMV+IAC assay has a limit of detection of ≤ 5 DNA copies. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended.

Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 11 µl	Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at -20°C and protected from light.

Additionally required equipment and materials

- real-time PCR instrument, equipped with four detection channels (522 nm, 553 nm, 610 nm, 670 nm)
- real-time PCR consumables (plates, tubes, foils, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

² Note: For the use of the Roche LightCycler® 480 instrument a Color Compensation is necessary. The SureCC Color Compensation Kit I (Art. No. F4009) must be used for the color compensation of such devices.

Protocol

1. Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and extraction control. Beside the 35S/NOS/FMV detection assays the Reaction Mix includes an internal amplification control (inhibition control) for each reaction.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20.0 µl	220.0 µl

Mix each master-mix well and centrifuge shortly before use.

2. Setup

	Blockcycler / LightCycler® 480	Rotor-Gene Q
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Detection 35S: in the FAM-Channel (LC480 I 483 nm, 533 nm; LC480 II 465 nm, 510 nm) Detection NOS: in the Cy5-Channel (LC480 I 615 nm, 670 nm; LC480 II 618 nm, 660 nm) Detection FMV: in the ROX-Channel (LC480 I 483 nm, 610 nm; LC480 II 533 nm, 610 nm) Detection internal amplification control: in the VIC/HEX-Channel Stratagene - HEX; CFX96 - VIC) (LC480 I 523 nm, 568 nm; LC480 II 533 nm, 580 nm;)	Detection: End of extension phase Detection 35S: in the FAM-Channel (Green) Detection NOS: in the Cy5-Channel (Red) Detection FMV: in the ROX-Channel (Orange) Detection internal amplification control: in the VIC-Channel (Yellow)
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: http://www.congen.de/en/company/downloads		

3. Preparation of the PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes.
- Pipette 5 µl of the Positive Control into the designated tubes.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the PCR instrument and start the run according to the setup.

Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer. The control reactions need to give the correct results.

The 35S promoter DNA is detected in the FAM-channel, the NOS terminator DNA is detected in the Cy5-channel, the FMV promoter DNA is detected in the ROX-channel and the internal amplification control is detected in the VIC/HEX-channel.

A sample is stated **positive** for the respective parameter (35S/NOS/FMV), if the sample DNA shows amplification in the respective channel and the internal signal of the sample shows amplification in the VIC/HEX-channel. A sample is stated **negative** for the respective parameter (35S/NOS/FMV), if the sample DNA shows no amplification in the respective channel and the internal signal (VIC/HEX) of the sample is **positive** (see also table below).

result in the respective channel				
FAM channel 35S promoter	Cy5 channel NOS terminator	ROX channel FMV promoter	VIC/HEX channel amplification control	result
positive	negative	negative	positive	35S promoter DNA detected ¹
negative	positive	negative	positive	NOS terminator DNA detected
negative	negative	positive	positive	FMV promoter DNA detected
negative	negative	negative	negative	invalid ²

¹ The 35S promoter is originally derived from the Cauliflower Mosaic Virus. Brassicaceae species are the host of the CaMV. Therefore, the analysis for the 35S promoter also detects the naturally occurring virus. It is possible that a sample containing the virus yields a positive result, although no GMO is present. As a supplementary kit, CONGEN offers the SureFood® GMO SCREEN CaMV kit, which helps to assure the absence of the natural Cauliflower Mosaic Virus in the sample in case of a single positive 35S result.

² In the case of a negative result in all channels including the VIC/HEX-channel the sample contains PCR inhibiting substances. Under these circumstances an evaluation of the sample is not possible. DNA isolation and purification for the sample need to be improved.

It may appear in some cases that only one of the two DNAs prepared from a test sample is 35S and/or NOS and/or FMV **positive**. This indicates that the amount of genetically modified DNA is very low and at the limit of detection. If such results are obtained in at least two repetitions of the analysis, the sample is stated negative.

Product Information

- Validation Report

Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

Distribution and ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com