**CONGEN**

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## 1 General Information

### 1.1 Description

The test detects pork (*Sus scrofa*), chicken (*Gallus gallus*) and turkey (*Meleagris gallopavo*) DNA. Each reaction contains an internal amplification control and an internal detection assay for vertebrates DNA (IAAC). The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of four fluorescence emissions at 522 nm, 553 nm, 610 nm and 670 nm (FAM, VIC, ROX and Cy5) at the same time. The technical validation of instruments was performed on Roche LightCycler® 480 II/ Roche cobas z 480 Analyzer<sup>2</sup>, Qiagen Rotor-Gene Q, Applied Biosystems 7500, Bio-Rad CFX96, Agilent Mx3005P and Agilent AriaDx.

### 1.2 Limit of Detection

The SureFood® ANIMAL ID 4plex Pork/Chicken/Turkey+IAAC real-time PCR is developed for the detection of pork, chicken and turkey DNA in muscle meat mixture at a relative amount of 0.1 %. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

**Note:** In mixed samples, inconsistent mixing ratios (e.g. 99.9 % pork and 0.1 % chicken) may cause a loss of sensitivity in the low concentration channel.

### 1.3 DNA-preparation

For DNA-preparation the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended.

### 1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 11 µl	Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at –20°C and protected from light.

### 1.5 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFood® PREP Basic or SureFood® PREP Advanced)
- real-time PCR instrument with four detection channels (522 nm, 553 nm, 610 nm and 670 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps) pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

<sup>2</sup> note: For the use of the Roche LightCycler® 480 II a Color Compensation is necessary. The SureCC Color Compensation Kit I (Art. No. F4009) must be used for the color compensation of such devices.

January 2019

## 1.6 Setup

	<b>Blockcycler/MIC/LightCycler® 480 II</b>	<b>Rotorgene Q</b>
Initial Denaturation (HOLD) Cycles Denaturation Annealing/Extension (CYCLE)	5 min, 95°C 35 15 sec, 95°C 30 sec, 60°C	1 min, 95°C 35 10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup (exemplary)	Detection: End of Extension Phase  Detection System pork: Various devices <b>FAM</b> -channel, Quencher: BHQ LightCycler 480 II 465 nm - 510 nm  Detection System vertebrates and internal amplification control (IAAC): Various devices <b>VIC/HEX</b> -channel, Quencher: BHQ LightCycler 480 II 533 nm - 580 nm  Detection System turkey: Various devices <b>ROX</b> -channel, Quencher: BHQ LightCycler 480 II 533 nm - 610 nm  Detection System chicken: Various devices <b>Cy5</b> -channel, Quencher: BHQ LightCycler 480 II 618 nm - 660 nm	
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: <a href="http://www.congen.de/en/company/downloads">http://www.congen.de/en/company/downloads</a>		

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and extraction control. The master-mix includes an internal amplification control (inhibition control) for each reaction.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

<b>Components for master-mix</b>	<b>Amount per reaction</b>	<b>10 reactions (with 10% excess)</b>
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix each master-mix well and centrifuge shortly before use.

**2.1.2 Preparation of the real-time PCR-mix**

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries.
- Pipette 5 µl of the Positive Control into the designated tubes/wells or capillaries.
- Centrifuge all tubes/wells or capillaries shortly at low speed.
- Place tubes/wells or capillaries into the real-time PCR instrument and start the run according to the setup.

**2.2 Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The results of the control reactions have to be correct.

Pork DNA is detected in the FAM-channel, turkey DNA is detected in the ROX-channel and chicken DNA is detected in the Cy5-channel. In the VIC-channel it is possible to detect vertebrates DNA in the sample as well as the amplification control (IAC) in a sample with no vertebrates DNA inside.

A sample is stated **positive** for the respective parameter (pork/turkey/chicken), if the sample DNA shows amplification in the respective channel. A sample is stated **negative** for the respective parameter (pork/turkey/chicken), if the sample DNA shows no amplification in the respective channel and the internal signal (VIC-channel) of the sample is **positive** (see also table below).

result in the respective channel				result
FAM - channel pork	VIC/HEX - channel <sup>1</sup> vertebrates + IAC	ROX - channel turkey	Cy5 - channel chicken	
<b>positive</b>	<b>positive</b>	negative	negative	<b>pork DNA detected</b>
negative	<b>positive</b>	negative	negative	<b>vertebrates DNA<sup>2</sup> detected</b>
negative	<b>positive</b>	<b>positive</b>	negative	<b>turkey DNA detected</b>
negative	<b>positive</b>	negative	<b>positive</b>	<b>chicken DNA detected</b>
negative	negative	negative	negative	Invalid <sup>3</sup>

<sup>1</sup> If the signal of the internal control (VIC-channel) of the sample DNA is detected significantly before the signal of the negative control (master-mix without DNA) the sample contains vertebrates DNA.

<sup>2</sup> Is the Cp-value of the internal control (VIC-channel) in the range of the negative control the sample contains no PCR-inhibiting substances but only a low amount or no vertebrates DNA.

<sup>3</sup> In the case of a negative result in all channels including the VIC-channel the sample contains PCR inhibiting substances. Under these circumstances an evaluation of the sample is not possible. DNA isolation and purification for the sample need to be improved.

**Note:** If the sample contains more than one animal species the DNA mixture can have a competitive influence on the absolute fluorescence. The lower the concentration of the determinant DNA is in a mixture of animal DNAs the lower is the fluorescence level of the amplification curve.

January 2019

### **3 Further Information**

#### **3.1 Product Information**

- Validation Report

#### **3.2 Technical Support**

For further questions please contact your distributor or send an e-mail to [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

#### **3.3 Distribution and ordering**

R-Biopharm AG  
An der neuen Bergstrasse 17,  
64297 Darmstadt, Germany  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-20  
E-Mail: [orders@r-biopharm.de](mailto:orders@r-biopharm.de)  
[www.r-biopharm.com](http://www.r-biopharm.com)