RIDASOFT®Win.NET Food & Feed

REF **Z9996FF**

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RIDASOFT® Win.NET Food & Feed

1. Intended use

RIDASOFT® Win.NET Food & Feed is a software for evaluating immunoassays for food and feed analytics. The software processes data, which is transmitted by a serial interface.

2. Minimum System requirements

- OS Win XP SP3 or higher
- .NET Framework from Version 4.0
- For the installation of the software: administration rights

3. Installation

- Insert the installation disk into the disk drive
- Open *setupDotNet.exe*
- The setup program will install RIDASOFT® Win.NET Food & Feed on your computer
- Start RIDASOFT® Win.NET Food & Feed
- Select your preferred language
- You will be asked if you would like step-by-step instructions if you confirm you will be introduced with the program functionality step by step
- If you do not wish to be asked again each time the program starts please click the check button
- All help topics including the test implementation can be accessed by using the *help index* command in the help menu

4. Use of the software

4.1 Starting the program

4.1.1 General

Start the RIDASOFT®Win.NET Food & Feed, the *select category* window opens, select the category and confirm your choice by clicking on the *OK* button. All the subfolders stored for the selected item are listed in the window that follows.

Select the subfolder and then the method you require and confirm your choice with *OK*. See point 4.2. for changing the settings of the methods.

4.1.3 Most-recently accessed method file

Start the RIDASOFT®Win.NET Food & Feed and the *Select category* window opens, select the category *most-recently accessed method file* and confirm your choice by clicking on the *OK* button. See point 4.2. for changing the settings of the methods.

4.2 Changing a method

Change a method by clicking on this item or on the menu item *File/Open/Method*. Select the method you want to change from the list and confirm your choice with *OK*.

In order to change this method, click on this icon or on the menu item *Edit/Method*. You can choose the specifications for the method in the window that then opens:

4.2.1 Standards

Quantity: Enter the desired number of standards. You can now switch between the standards one by one with the cursor buttons and define their properties.

ID (optional): Use a maximum of 20 characters to name the standard. Initially standards are named Std1...Stdn. The standards' names are visible in the method view and in the print output.

Replicates: Enter the desired number of replicates for each standard. If you would like to use the same number of replicates for all standards, please press the auto button after putting in the desired number.

Concentration: Enter the desired concentration for each standard. Zero concentrations are only allowed for the first standard, which is the zero standard.

Serial dilution: If you put in a dilution factor greater than zero, the concentration of each standard beginning with the highest is calculated. Enter the concentration for the highest standard in the concentration field. The concentration of the remaining standards is calculated using following formula: $\text{Conc. Std}_i = \text{Conc. Std}_{i+1} / \text{dilution factor}$.

4.2.2 Properties

Name: The name states the name of the method to be changed. When saving edited methods you have to use new filenames. The supplied method files are write protected and not changeable.

Dimension: Input the unit in which you want the results displayed, usually ppm, ppb, etc. for Food and Feed analysis.

Evaluation: The evaluation indicates the calculation for the method. Depending on the selected item, there are different calculation methods available. You can choose between linear regression, cubic spline and logit/log, for example.

4.2.3 Samples

Amounts: Reflects the remaining maximum amount of available samples. Switch the samples one by one with the cursor buttons and define their properties.

Id (optional): Use a maximum of 20 characters to name the samples. Initially samples are named P1 ... Pn. The sample names are visible in the method view and in the print output.

Replicates: Enter the desired number of replicates for each sample. If you would like to use the same number of replicates for all samples, please press the auto button after inputting the desired number.

Factor (optional): Enter the desired dilution factor for each sample. If you would like to use the same dilution factor for each sample, please press the auto button after inputting the desired factor. If you would like to change only the samples' dilution factor for the current method, it is more suitable to use the Edit Samples command.

4.2.4 Controls

Amount: Enter the desired number of controls. You can now switch between the samples one by one with the cursor buttons and define their properties.

Id (optional): Use a maximum of 20 characters to name the controls. Two controls are normally used and are named QC1 and QC2. The control names are visible in the method view and in the print output.

Replicates: Enter the desired number of replicates for each control. If you would like to use the same number of replicates for all controls, please press the auto button after inputting the desired number.

4.3. Preparation of the reader for your assay measurements

Before performing your first assay measurement, you have to adjust RIDASOFT® Win.NET Food & Feed to your photometer.

To determine the photometer and the serial port use the *Reader Setup* command. Select your reader in the dialog box. The following dialog box determines the serial port settings:

The dialog controls are initialized with the selected device specifications. These specifications refer to the reader manual. Please compare the settings. It is possible that your reader uses different specifications. In this case, you have to change the settings:

Port: Select the serial port connected to your reader.

Baud Rate: Select the speed the reader transfers its data.

Data Bits: Select the number of data bits.

Parity: Select the kind of parity bit.

Stop Bits: Select the number of stop bits.

Flow (optional): Select the kind of data flow control. Normally readers use a hardware handshake for communication with the computer. Please change this parameter only if you know that your reader uses a different kind of flow control than the one already checked.

If the device is remotely controlled by RIDASOFT® Win.NET Food & Feed, please select the filter wavelengths in the following step.

These settings are saved as a file and are used every time you start a new measurement.

Important:

Be sure your reader works in computer mode. Readers obtained from the manufacturer are in some cases set to printer mode, thus the data is transferred to a printer connected to the reader. In this case you have to switch to computer mode. Please refer to the reader manual. Be sure the serial cable is RS-232C standard. In most cases you must obtain the serial cable from the reader manufacturer.

4.4 Plate layout

4.4.1 Methods

The layout of the plate with standards and controls is already stipulated. The layout shifts if you change the characteristics of the method (number of standards, number of repeats, etc., see point 4.2.). If you want to alter the layout made automatically, mark the *BLK* entry in the small right-hand window. Then position the mouse over all entries (standards and controls) one after the other and press the left-hand mouse button each time.

In order to now occupy the plate with the standards and controls according to the stipulations of the selected method, mark the desired entry in the list and then position the mouse over the desired position in the layout window and then press the left-hand mouse button.

Select *P* in the list to occupy your samples and click on as many fields in the layout window as the number of samples you want to specify. In order to allocate an ID to the samples, select the menu item *Edit/Samples*. A window then opens in which you can allocate an ID to each sample, which then appears in the printout of results. You can delete an existing entry with *BLK*.

You can occupy the plate in one step with *Auto*. Mark any assignment of fields and the standards/sample are occupied according to their number and replicates in the desired allocation in chronological sequence on the plate.

4.5 Measuring

Before you can measure the occupied and processed plate, you must ensure that the photometer is connected and switched on and that the correct reader is set in the software (see point 4.3.).

Place the plate in your photometer and start the measurement either by clicking on this icon or by clicking on the menu item *File/New/Assay*.

4.6 Evaluation

Click on this icon or click on the menu item *Results/Default* to start the evaluation. You can view the measured values, the plate occupancy, the method used, the standard curve (quantitative evaluations) and the sample results on the lower edge of the window by clicking on the respective tab.

Save the test by clicking on this icon or click on the menu item *File/Save/Assay* and enter a name under which you want to save the test.

Stored tests can be reopened by clicking on this icon or by clicking on the menu item *File/Open/Assay* and then on the name under which the test was stored. As only the raw data were stored, the evaluation must be made again to enable the results to be visualised.

4.7 Printing and exporting data

Click on this icon or click on the menu item *File/Print* to print the results. Select Results if you want to print out all results, including the standard curve of the plate layout etc., in list form. Click on the desired printouts if you do not want to print out all data. In order to make the results available to other database programs, click on the menu item *File/Save Sample Results*. The data is then stored in the RIDADB.TXT file and can be read in by other database programs.

5. Explanation of the command line

5.1. File menu

New Method command (same as Edit Method command)

Use this command to create a new method based on the current method. You can open an existing method with the *Open Method* command (Shortcut Keys: CTRL+N).

New Assay command

Use this command to start a new measurement. You need to start the reader after activating this command for manual reader operation. The data can be seen in the plate values view. You can open an existing assay with the Open Assay command (Shortcut Keys: CTRL+M).

Open Method command

Use this command to open an existing method. You can create a new method with the New Method command (Shortcut Keys: CTRL+O).

Open Assay command

Use this command to open an existing assay. You can start a new assay measurement with the New Assay command (Shortcut Keys: CTRL+A).

Save Method command

Use this command to save and name the active method along with the current plate layout. RIDASOFT® Win.NET Food & Feed displays the Save As dialog box so you can name your method. The supplied method files are write protected (Shortcut Keys: CTRL+S).

Save Assay command

Use this command to save and name the active assay along with the current plate layout. RIDASOFT® Win.NET Food & Feed displays the Save As dialog box so you can name your assay. Existing assay files are write protected (Shortcut Keys: CTRL+R).

Save Sample Results command

Use this command to save the current samples' data. The file RIDADB.txt is automatically updated. You can import this file to a database application to create your own database.

Classify Samples command

Use this command to classify selected samples (for example suppliers, customers.).

Add (a category):

Press this button after filling the edit box and the list will be updated.

Remove (a category):

Press this button after selecting a list item and the item will be removed from the list.

Classify:

Press this button after selecting the desired samples in the layout view and after selecting the desired list item. The samples will be classified with the selected list item. Successful classification will be confirmed with an 'ok' message.

This command is unavailable if the layout view is not active.

The classification of samples is displayed in the method view of the method layout. Use the Save Sample Result command to add the data to the samples list.

Print command

Use this command to print a selection of method specific data. This command presents the Selection dialog box, where you may select the desired information to be printed (Shortcut Keys: CTRL+P).

5.2. Edit menu

Cut command

Use this command to remove the currently selected data from the document and copy it to the clipboard. This command is unavailable if there is no data currently selected. Cutting data to the clipboard replaces the contents previously stored there (Shortcut Keys: CTRL+X).

Copy command

Use this command to copy selected plate values onto the clipboard. This command is unavailable if there is no plate value currently selected. Copying data to the clipboard replaces the contents previously stored there (Shortcut Keys: CTRL+C).

Paste command

Use this command to insert a copy of the clipboard contents at the insertion point. This command is unavailable if the clipboard is empty (Shortcut Keys: CTRL+V).

Method command (same as New Method command)

Use this command to edit the current method. You can open an existing method with the *Open Method* command (Shortcut Keys: CTRL+N).

Samples command

Use this command to change the sample dilution factor and/or the sample ID for the current method.

Factor:

Put in the desired dilution factor for each sample. If you like to use the same dilution factor for each sample, please press the auto button after putting in the desired factor. Use the cursor buttons to switch between samples.

ID (optional):

Use a maximum of twenty characters to name the samples. Initially samples are named P1 ... Pn. The sample names are visible in the method view and in the print output.

5.3. Results menu

Default command

Use this command to start the method specific calculation. The view will be switched to method layout, and you will have access to the results in the graph view (quantitative evaluations) and in the samples view.

4P-Method command

Use this command to calculate the 4P-curve from the standard OD values. The view will be switched to method layout, and you will have access to the results in the graph view and in the samples view. This command is unavailable if the current method uses a qualitative or point-to-point calculation.

Linear Regression command

Use this command to calculate the regression line from the standard OD values. The view will be switched to method layout, and you will have access to the results in the graph view and in the samples view. This command is unavailable if the current method uses a qualitative or point-to-point calculation.

Logit/Log command

Use this command to calculate the Logit/log from the standard OD values. The view will be switched to method layout, you have access to the results in the graph view and in the samples view. This command is unavailable if the current method uses a qualitative or point-to-point calculation.

Cubic Spline command

Use this command to calculate the cubic spline from the standard OD values. The view will be switched to method layout, and you will have access to the results in the graph view and in the samples view. This command is unavailable if the current method uses a qualitative or point-to-point calculation.

5.4. View menu

Toolbar command

Use this command to display and hide the Toolbar, which includes buttons for some of the most common commands in RIDASOFT® Win.NET Food & Feed, such as *File Open*. A checkmark appears next to the menu item when the Toolbar is displayed.

Status Bar command

Use this command to display and hide the Status Bar, which describes the action to be executed by the selected menu item or depressed toolbar button, and keyboard latch state. A checkmark appears next to the menu item when the Status Bar is displayed.

5.5. Windows menu

Plate Layout command

 Use this command to switch to the Plate Layout mode. This command is unavailable if the mode is already active (Shortcut Keys: CTRL+L).

Method Layout command

 Use this command to switch to the Method Layout mode. In this mode you can check the template of the current method in tabbed views. When performing a calculation you can check the results in the updated view. This command is unavailable if the mode is already active (Shortcut Keys: CTRL+H).

5.6. Help menu

Index command

 Use this command to display the opening screen of *Help*. From the opening screen, you can jump to step-by-step instructions for using RIDA®SOFT Win.NET Food & Feed and various types of reference information. Once you open Help, you can click the Contents button whenever you want to return to the opening screen.

Using Help command

Use this command for instructions about using Help.

About Ridasoft Win command

Use this command to display the copyright notice, version number and the serial number of your copy of RIDA[®]SOFT Win.NET Food & Feed.

6. Troubleshooting

| Error Description | Possible Source | Solution |
|---|--|---|
| Error message during measurement “Reader not found or reader not started!” | Reader not correctly connected or not turned on | Connect and/or turn on the reader correctly (see 4.3.) |
| Error message during measurement „Data transfer failed“. | Reader is not working properly | Check functionality of the reader |
| Error message in reader properties: „Connection failed”. | Reader not correctly connected or not turned on | Connect and/or turn on the reader correctly (see 4.3.) |
| Error during reading/initilizing/writing of files | No reading/writing permission to the program directory | Grant reading/writing permissions |
| Connection issues with Biotek ELX800 | Wrong baudrate set and/or data flow control DTR/DSR not active | Set baudrate DTR/DSR 2400 and/or activate data flow control DTR/DSR (see 4.3) |
| Connection issues with Multiskan EX | Data flow control active | Disable all data flow control options (see 4.3) |
| Connection issues with Metertech Accureader | Data is not provided to the PC | Reader: SETUP->COMPUTER ->RAW-DATA Measurement: MODE->END POINT(S) ->CONTINUOUS |

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

| Version number | Chapter and title |
|----------------|---|
| 2015-02-20 | Release version |
| 2020-07-10 | Adaptation to new master template, change of the article number, all references to clinical diagnostics have been removed |

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.