CONGEN

SureFood® GMO QUANT MON810 Corn

Art. No. S2019

2 x 50 rxn

User Manual

May 2018

 **Inhalt** /  **Content**

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Mai 2018

3.1	Product Information	12
3.2	Technical Support	12
3.3	Distribution and Ordering	12

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	Blockcycler/LightCycler® 480/Bio Molecular Systems MIC	Rotorcycler / LTF MyGo Pro
Initial Denaturation (HOLD) Cycles	5 min, 95°C 45	1 min, 95°C 45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye: FAM Quencher Dye: BHQ Passive Reference: none	LightCycler® 2.0 Channel: 530 Acquisition mode: Single in extension phase Rotor-Gene Q Reporter Dye: FAM (Green)
Detaillierte Informationen zur Einstellung bestimmter real-time PCR Geräte stehen auf der CONGEN-Homepage zur Verfügung: http://www.congen.de/unternehmen/download		

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1 General Information

1.1 Description

The test detects the relative quantitative MON810 Corn DNA amount. Therefore the kit contains two PCR systems, one specific for the MON810 Corn (OECD unique identifier MON-00810-6) and the other one specific for Corn (the reference gene). The event specific detection of MON810 Corn is according to EN ISO 21570 and the official method of the European Commission. The real-time PCR assay can be used with established real-time PCR instruments. The technical validation of instruments was performed on Agilent Mx3005P, Agilent Aria Dx, BioRad CFX 96, Roche LightCycler® 480 II, Roche cobas Z 480 Analyzer, Roche LightCycler® 2.0, Applied Biosystems 7500, Qiagen Rotor-Gene Q and LTF MyGo Pro.

1.2 Limit of Detection

The MON810 corn has a limit of detection of ≤ 5 DNA-copies. The assay limit of detection depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The limit of quantitation depends on the concentration of the sample DNA used in the analysis. For example, if 50,000 target-sequence copies of Corn specific reference gene are present, the relative quantitation limit for MON810 DNA is 0.1 %.

1.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Corn Reaction Mix	1 x 1100 µl	Orange
2	MON810 Reaction Mix	1 x 1100 µl	Yellow
3	Taq Polymerase	1 x 11 µl	Red
4	Dilution Buffer	1 x 1300 µl	White
5	Standard DNA	1 x 50 µl	Dark Blue
6	Positive Control	1 x 100 µl / 1 % MON810	Light Blue

Store all reagents at -20°C and protected from light.

1.5 Additionally required equipment and materials

- real-time PCR instrument
- real-time PCR consumables (plates, tubes, foils, capillaries, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Setup

	Blockcycler/LightCycler® 480/Bio Molecular Systems MIC	Rotorcycler / LTF MyGo Pro
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup (exemplary)	Detection: End of extension phase Reporter Dye: FAM Quencher Dye: BHQ Passive Reference: none	LightCycler® 2.0 Channel: 530 Acquisition mode: Single in extension phase Rotor-Gene Q Reporter Dye: FAM (Green)
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: http://www.congen.de/en/company/downloads		

2 Quantitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples, controls and standards).

Reactions needed for the corn detection:

- For each run: 5 reactions for the standard curve
3 reactions for controls (1x no-template and 2x Positive Control)
 - For each sample: at least 1 reaction with each sample DNA
- Reactions needed for the MON810 corn detection:

- For each run: 5 reactions for the standard curve
3 reactions for controls (1x no-template and 2x Positive Control)
- For each sample: at least 1 reaction with each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Corn Reaction Mix or MON810 Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220.0 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the master-mix

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.3 Preparation of the standard DNA dilutions

Dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different DNA concentrations for the standard curves of the Corn reference gene and the MON810 Corn detection gene. Prepare 5 dilutions of the supplied Standard DNA with the supplied Dilution Buffer. Prepare 5 reaction tubes (labeled S1 to S5) and add 45 µl Dilution Buffer (**Code 4**) each. The following procedure is recommended:

standard	dilution	copy number per µl	final copy number per reaction*
S1	45 µl Dilution Buffer + 5 µl Standard DNA	100,000 copies	500,000 copies
S2	45 µl Dilution Buffer + 5 µl DNA of S1	10,000 copies	50,000 copies
S3	45 µl Dilution Buffer + 5 µl DNA of S2	1000 copies	5000 copies
S4	45 µl Dilution Buffer + 5 µl DNA of S3	100 copies	500 copies
S5	45 µl Dilution Buffer + 5 µl DNA of S4	10 copies	50 copies

***Note:** 5 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

If the diluted DNA standards (S1 to S5) are not immediately used, store them at -20°C. The dilutions are stable up to two months at -20°C. Before use allow the reagents to thaw, mix them on a vortex and centrifuge carefully before opening and use.

2.1.4 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries and close them.
- Pipette 5 µl of the Positive Control and the standard dilutions into the designated tubes/wells or capillaries and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates or capillaries into the PCR instrument and start the run according to the setup.

May 2018

2.2 Interpretation of results

The calculation for both reaction systems (**Corn**, **MON810 Corn**) has to be made separately. Mark the standards, the controls and the samples for the specific system (**MON810 Corn**) and make the evaluation according to the usual analysis program recommended by the real-time PCR instrument manufacturer. Repeat the same procedure for the **Corn** reference gene system. The value for the slope of the standard curve has to be between -3.1 and -3.6 and the correlation coefficient $R^2 > 0.98$. In case of different values for the standard curve, it should not be used for calculation.

By using the calculated copy numbers for **MON810 Corn** and **Corn** the relative GMO content of the sample DNA and the Positive Control can be determined in the following way (example):

Sample MON810 Corn	1350 copies	Positive Control MON810 Corn	400 copies
Sample Corn	45,000 copies	Positive Control Corn	28,000 copies

Multiply the copy number of the specific system by 100 and divide by the copy number of the reference gene system to obtain the percentage.

MON810 Corn content = $\text{MON810 Corn copy number} * 100 / \text{Corn copy number}$

sample DNA **MON810 Corn** content = $1350 * 100 / 45,000$ sample DNA **MON810 Corn** content = 3 %

For the given example the numbers lead to a **MON810 Corn** content of 3.0 % for the sample, and 1.4 % for the Positive Control with the same calculation.

For a final calculation the use of a correction factor K for the correction of run-to-run fluctuations is necessary. The correction factor is the relation of the true percentage value of the Positive Control (1 % GMO content) and the measured GMO percentage of the Positive Control. The factor is calculated in the following way:

$K = \text{true GMO percentage of Positive Control} / \text{measured GMO percentage of Positive Control}$

K (example) = $1 \% / 1.4 \% = 0.7$

The calculated value for the sample is multiplied with K to obtain a corrected value.

GMO percentage sample = $\text{measured GMO percentage of sample} * K$

sample (example) = $3.0 \% * 0.7 = 2.1 \%$

For that example the **MON810 Corn** content is 2.1 %.

May 2018

3 Further Information

3.1 Product Information

- Validation Report
- Microsoft Excel template of calculation
(Download: www.congen.de/en/company/downloads)

3.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

3.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com