CONGEN

SureFood[®] GMO ID Bt63 Rice

Art. No. S2024

2 x 50 rxn

User Manual

July 2018

 Inhalt /  Content

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	Blockcycler / LightCycler® 480	Rotorcycler
Initial Denaturation (HOLD) Cycles Denaturation Annealing/Extension (CYCLE)	5 min, 95°C 45 15 sec, 95°C 30 sec, 60°C	1 min, 95°C 45 10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye: FAM Quencher Dye: BHQ Passive Reference: none	LightCycler® Channel: 530 oder F1 Acquisition mode: Single in extension phase Rotor-Gene Q Reporter Dye: FAM (Green)
Detaillierte Informationen zur Einstellung bestimmter real-time PCR Geräte stehen auf der CONGEN-Homepage zur Verfügung: http://www.congen.de/unternehmen/download		

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1 General Information

1.1 Description

The test detects Bt63 Rice DNA. The detection of Bt63 Rice is a construct-specific real-time PCR assay. Additionally the SureFood® GMO ID Bt63 Rice kit contains a reference-PCR system for the detection of rice. The real time PCR system used in the SureFood® GMO ID Bt63 rice test is according to the method published by: D. Maede et al.; European Food Research and Technology, 224:271-278 (2006) and therefore according to Commission decision 2008/289/EC. The real-time PCR assay can be used with established real-time PCR instruments (Roche LightCycler®, Rotor-Gene Q, ABI PRISM, Eppendorf realplex, BioRad CFX96, Agilent MxSeries etc.).

1.2 Limit of Detection

The P35S:BAR PCR has a limit of detection of ≤ 5 DNA copies. This is equivalent to approx. 0.01 % for unprocessed Rice grain. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

1.3 DNA-preparation

The real-time PCR kit can be used in combination with the SureFood® PREP Basic or SureFood® PREP Advanced Kit. The use of SureFood® PREP Advanced kit protocol 1 (preparation of 100 mg of the homogenized sample with 1.3 ml Lysis Buffer) increases the DNA yield. This ensures the required limit of detection of 0.01 % in unprocessed Rice grain.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Rice Reaction Mix	1 x 1100 µl	Orange
2	Bt63 Rice Reaction Mix	1 x 1100 µl	Yellow
3	Taq Polymerase	1 x 11 µl	Red
4	Positive Control	1 x 200 µl	Light Blue

Store all reagents at -20°C and protected from light.

1.5 Additionally required equipment and materials

- real-time PCR instrument
- real-time PCR consumables (plates, tubes, foils, capillaries, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Setup

	Blockcycler / LightCycler® 480	Rotorcycler
Initial Denaturation (HOLD) Cycles Denaturation Annealing/Extension (CYCLE)	5 min, 95°C 45 15 sec, 95°C 30 sec, 60°C	1 min, 95°C 45 10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye: FAM Quencher Dye: BHQ Passive Reference: none	LightCycler® Channel: 530 or F1 Acquisition mode: Single in extension phase Rotor-Gene Q Reporter Dye: FAM (Green)
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: http://www.congen.de/en/company/downloads		

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control, extraction control and inhibition control for each sample. For the preparation of the inhibition control the use of the SureFood® GMO Plant PLUS (Art. No. S2049) is recommended.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Rice Reaction Mix or P35S:BAR Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20.0 µl	220.0 µl

Mix the master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries and close them.
- Pipette 5 µl of the Positive Control into the designated tubes/wells or capillaries and close them.
- Centrifuge all tubes/wells or capillaries shortly at low speed.
- Place tubes/wells or capillaries into the PCR instrument and start the run according to the setup.

July 2018

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer. The control reactions need to give the correct results.

A sample is stated **positive**, if the sample DNA shows amplifications in the Bt63 Rice system and the Rice reference system. A sample is stated **negative**, if the sample DNA shows no amplification in the Bt63 Rice system. In case of a **negative** result in the Bt63 Rice system the inhibition control of the sample must be **positive**. Is this not the case the sample contains PCR-inhibiting substances. Under these circumstances an evaluation of the samples is not possible. DNA isolation and purification for the sample need to be improved. If the sample DNA is negative for Rice, then either the sample contains no rice material or the DNA preparation was not successfully.

3 Further Information

3.1 Product Information

- Validation Report

3.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

3.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com