CONGEN

SureFood® GMO SCREEN CaMV

Art. No. S2027

100 rxn

User Manual

November 2018

 Inhalt /  Content

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	Blockcycler/LightCycler® 480/RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD) Cycles Denaturation Annealing/Extension (CYCLE)	5 min, 95°C 45 15 sec, 95°C 30 sec, 60°C	1 min, 95°C 45 10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye: FAM Quencher Dye: BHQ Passive Reference: none	LightCycler® Channel: 530 oder F1 Acquisition mode: Single in extension phase Rotor-Gene Q Reporter Dye: FAM (Green)
Detaillierte Informationen zur Einstellung bestimmter real-time PCR Geräte stehen auf der CONGEN-Homepage zur Verfügung: http://www.congen.de/unternehmen/download		

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1 General Information

1.1 Description

The test detects CaMV DNA. The presence of CaMV DNA indicates an infestation with the cauliflower mosaic virus. Host of the CaMV are *Brassicaceae* plants. With the virus specific detection system SureFood® GMO CaMV it is possible to verify the detection of a genetic modification. The real-time PCR assay can be used with established real-time PCR instruments (Roche LightCycler®, Rotor-Gene Q, Agilent Mx3005P, BioRad CFX 96, Applied Biosystems 7500, Agilent Aria Dx, R-Biopharm RIDA®CYCLER etc.).

1.2 Limit of Detection

The SureFood® GMO SCREEN CaMV PCR has a limit of detection of ≤ 5 DNA copies. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

1.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 μ l	Yellow
2	Taq Polymerase	1 x 11 μ l	Red
3	Positive Control	1 x 200 μ l	Light Blue

Store all reagents at -20°C and protected from light.

1.5 Additionally required equipment and materials

- real-time PCR instrument
- real-time PCR consumables (plates, tubes, foils, capillaries, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Setup

	Blockcyler/LightCycler® 480/RIDA®CYCLER	Rotorcyclcr
Initial Denaturation (HOLD) Cycles Denaturation Annealing/Extension (CYCLE)	5 min, 95°C 45 15 sec, 95°C 30 sec, 60°C	1 min, 95°C 45 10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye: FAM Quencher Dye: BHQ Passive Reference: none	LightCycler® Channel: 530 or F1 Acquisition mode: Single in extension phase Rotor-Gene Q Reporter Dye: FAM (Green)
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: http://www.congen.de/en/company/downloads		

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control, extraction control and an inhibition control. For the preparation of the inhibition control the use of the SureFood® GMO SCREEN 4plex 35S/NOS/FMV+IAC kit (Art. No. S2126) or the SureFood® GMO SCREEN 35S/NOS/FMV kit (Art. No. S2026), respectively, is recommended.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix the master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries and close them.
- Pipette 5 µl of the Positive Control into the designated tubes/wells or capillaries and close them.
- Centrifuge all tubes/wells or capillaries shortly at low speed.
- Place tubes/wells or capillaries into the PCR instrument and start the run according to the setup.

November 2018

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions need to give the correct results.

A sample is stated **positive**, if the sample DNA shows an amplification signal in the CaMV system.

A sample is stated **negative**, if the sample DNA shows no amplification in the CaMV system and a performed amplification control (inhibition control) with S2026 or S2126 of the sample is **positive**. If the sample DNA and a performed amplification control are **negative** the sample contains PCR-inhibiting substances. Under these circumstances an evaluation of the samples is not possible. DNA isolation and purification for the sample need to be improved.

3 Further Information

3.1 Product Information

- Validation Report
- Microsoft Excel template of calculation
(Download: www.congen.de/en/company/downloads)

3.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

3.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com