# **RIDA<sup>®</sup>CHECK**

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## **1. Intended use**

RIDA<sup>®</sup>CHECK is a rapid swab test for monitoring the efficiency of surface cleaning procedures in the production environment. The test is intended to be used as a support for hygiene monitoring as well as general protein screening in allergen management.

## **2. General**

Within the scope of a reasonable HACCP concept it is not sufficient to only carry out microbiological status within a production environment after the food is manufactured, as this control would only monitor the success of sanitation steps within the hygiene program of the food factory. Residues of foods left remaining on surfaces may additionally provide the ideal nutrient base for airborne microbes which can re-contaminate disinfected surfaces. Hence an efficient cleaning of surfaces and equipment is needed and therefore it is expedient to regularly monitor the effectiveness of these procedures.

## **3. Test principle**

During production processes protein residues from raw materials or food stuffs may remain on surfaces. These should be removed during the normal cleaning process. However cleaning actions are not always successful and consequently contamination may still remain and possibly lead to the cause of future health issues.

RIDA<sup>®</sup>CHECK's methodology is based on the detection of protein residues and is independent of ATP. During the assay the swabbed protein residues react with the colour indicator. This interaction induces a pH decrease which leads to a rapid colorimetric reaction from yellow through to green. According to the intensity of the colour change (yellow → light green → green → dark green) the guideline determination of the contamination level is possible. The operator then determines whether the level of contamination in the tested area is still acceptable (according to routinely observed and recorded colour levels) to decide if corrective action has to be initiated.

#### **4. Reagents provided**

Each package of RIDA<sup>®</sup>CHECK Art.No. R1091 contains 50 swabs for the performance of 100 tests. Each package of RIDA<sup>®</sup>CHECK Art.No. R1092 contains 20 swabs for the performance of 40 tests.

#### **5. Materials required but not provided**

For the testing with RIDA<sup>®</sup>CHECK indicators no additional reagents are needed.

Laboratory equipment is not required.

#### **6. Warnings and precautions for users**

Used RIDA<sup>®</sup>CHECK indicators can be disposed with the regular household refuse. Even though the chemicals of the reaction liquid are not harmful to humans, animals or the environment it is recommended to rinse the sampled surface with alcohol or water subsequently to the application of the indicator. A relevant material safety data sheet is available on request.

#### **7. Storage instructions**

It is recommended to store the RIDA<sup>®</sup>CHECK indicators at temperatures between 4 °C (39.2 °F) and 24 °C (75.2 °F).

#### **8. Indication of instability or deterioration of reagents**

Do not use swabs which are actually dried out or where the colour already changed to green.

#### **9. Preparation of samples**

For the swab test on surfaces with RIDA<sup>®</sup>CHECK-Indicators no sample preparation is needed.

#### **10. Test implementation**

RIDA<sup>®</sup>CHECK indicators are ready-to-use and applicable without any preparation of reagents. Touching the cotton heads which are moisturised with reagent liquid should be avoided.

Performance:

1. Rip the packaging of the swab at its perforation
2. Hold the swab on one side with one hand and remove the other side of the packaging with the other hand
3. Firmly wipe the surface with the unpackaged part of the RIDA<sup>®</sup>CHECK indicator (Test 1)
4. Return the used part of the swab into its packaging
5. Remove the packaging of the second indicator and firmly wipe the surface of another sampling area (Test 2)\*

The size of the sampling area should equal to approximately 20 cm<sup>2</sup> (4.5 cm x 4.5 cm).

The developing time of the colour reaction depends on the level of contamination with protein residues on the sampled surface. In case of high pollution the colour change will be visible after only a few seconds. In case of only low pollution the time until the reaction ends could last 1 up to 2 minutes.

\* After opening the indicator package, both indicators have to be used within 1 hour.

## 11. Results

The interpretation of the results has to be done after the colour development is ended.

In case of the contamination of the surface with protein residues the colour of the indicator changes from yellow to green. The relative degree of surface contamination can be determined according to the following colour scale:

Colour	Contamination value
Yellow	Clean
Light green	Borderline
Green	Lightly contaminated
Dark green	Highly contaminated

A yellow 'clean' result indicates that the amount of protein is below the detection limit of 20 µg. The threshold has been determined with serum albumin. Due to the fact that food borne proteins are normally mixed combinations of different protein residues, a direct quantitative correlation between the intensity of indicator dye and amount of proteins isn't possible. The result of each test is only relative to Albumin and for this reason the interpretation of the degree of contamination is individual and cannot be compared directly with any individual protein standard. By routinely comparing test results to historical cleaning data the person in charge for the decontamination has to decide personally whether the cleansing of surfaces is acceptable or not and if corrective actions have to be undertaken.

Remark:

Possibly the RIDA<sup>®</sup>CHECK indicator may also detect residues of disinfectants. The strength of the colorimetric reaction depends on the amount and the kind of the disinfectant residue. For further information on the detection of disinfectants please refer to R-Biopharm.

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG