# **RIDASCREEN<sup>®</sup> FAST Aflatoxin**

Enzymimmunoassay zur quantitativen Bestimmung von  
Aflatoxin

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# RIDASCREEN<sup>®</sup>FAST Aflatoxin

## Brief information

RIDASCREEN<sup>®</sup>FAST Aflatoxin (Art. No.: R5202) is a competitive enzyme immunoassay for the quantitative analysis of aflatoxin in cereals and feed.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Technicians need no specialized training to perform the RIDASCREEN<sup>®</sup>FAST Aflatoxin test, however, free support is offered by the distributor on request, if necessary.

Sample preparation: extraction, filtration, dilution

Time requirement: sample preparation (for 10 samples)  
cereals and feed .....approx. 10 min  
test implementation (incubation time) ..... 15 min

Limit of detection: < 1.7 µg/kg (ppb)  
(corresponding to the standard substance)

## 1. Intended use

RIDASCREEN®FAST Aflatoxin assay is a competitive enzyme immunoassay for the quantitative analysis of aflatoxin in cereals and feed.

## 2. General

Aflatoxins are carcinogenic, highly toxic metabolites of the mould fungi *Aspergillus flavus* and *Aspergillus parasiticus*.

Aflatoxin B<sub>1</sub>, generally present together with aflatoxin B<sub>2</sub>, G<sub>1</sub> und G<sub>2</sub>, is occurring mainly in cereals, corn, cotton-seed and some nuts.

## 3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-aflatoxin antibodies.

Aflatoxin standards or sample solutions, aflatoxin enzyme conjugate and anti-aflatoxin antibodies are added. Free aflatoxin and aflatoxin enzyme conjugate compete for the aflatoxin antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-aflatoxin antibodies are also bound by the immobilized capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the aflatoxin concentration in the sample.

## 4. Reagents provided

Each kit contains sufficient materials for as many as 43 analyses (plus 5 standard analyses). Each test kit contains:

- 1 x Microtiter plate with 48 wells (6 strips with 8 removable wells each)  
coated with capture antibodies
- 5 x Aflatoxin standard solutions <sup>\*)</sup>, 1.3 ml each  
0 ppb (zero standard), 1.7 ppb, 5 ppb, 15 ppb, 45 ppb  
aflatoxin in methanol/water, ready to use
- 1 x Conjugate (3 ml) .....red cap  
peroxidase conjugated aflatoxin  
ready to use
- 1 x Anti-aflatoxin antibody (3 ml) ..... black cap  
ready to use
- 1 x Substrate/chromogen (10 ml) ..... brown cap  
stained red
- 1 x Stop solution (14 ml) .....yellow cap  
contains 1 N sulfuric acid
- 1 x Washing buffer (salt)  
for preparation of a 10 mM Phosphate Buffer, pH 7.4  
contains 0.05 % Tween 20

\*) The dilution factor 10 for the sample has already been considered. Therefore, the aflatoxin concentrations of samples can be read directly from the standard curve.

## 5. Materials required but not provided

### 5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- graduated cylinder (plastic or glass), 100 ml
- glassware for preparing sample extract: filter funnel and 50 ml flask
- grinder (mill)
- optional: shaker
- filter paper: Whatman No. 1 or equivalent
- 50 µl, 100 µl and 1000 µl micropipettes

### 5.2. Reagents:

- methanol
- 70 % methanol solution: prepare 70 % methanol solution by mixing 70 ml methanol (100 %) with 30 ml distilled or deionized water
- distilled or deionized water

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The standards contain aflatoxin, particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

Decontamination of the glassware and aflatoxin solutions is best carried out using a sodium hypochlorite (bleach) solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

The stop solution contains 1 N sulfuric acid (R36/38, S2-26).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen solution is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

The test kit can be regularly used at least up to the expiry date (indicated on the kit package), if stored correctly.

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- any bluish coloration of the reddish substrate/chromogen solution prior to test implementation
- a value of less than 0.6 absorbance units ( $A_{450\text{ nm}} < 0.6$ ) for the zero standard

## 9. Preparation of Samples

The samples should be stored in a cool place, protected from light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- weigh 5 g of ground sample into a suitable container and add 25 ml of methanol (70 %)\*
- shake vigorously for three minutes (manually or with shaker)
- filter the extract through Whatman No. 1 filter
- dilute 1 ml of the obtained filtrate with 1 ml of distilled or deionized water
- use 50 µl of the diluted filtrate per well in the test

\*) sample size may be increased if required, but the volume of methanol/water must be adapted accordingly, e.g. 10 g in 50 ml of methanol (70 %)

## 10. Test implementation

### 10.1. Preliminary comments

1. Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.
2. As **washing buffer** a PBS tween buffer is needed. Please use the buffer salt contained in the kit (see 4.). Dissolve the entire buffer salt in one liter of distilled water. The ready to use washing buffer expires after approx. 4 - 6 weeks 2 - 8 °C (36 - 46 °F).

Alternative: Dissolve the contents of the envelope in 100 ml of distilled water to obtain a 10fold concentrated washing buffer. The 10fold concentrate expires after approx. 8 - 12 weeks, store at room temperature (20 - 25 °C / 68 - 77 °F).

Use 1 part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use washing buffer.

3. The specific reaction starts with the addition of the specific antibody. Nevertheless, not more than three strips should be applied to the test, when a single step pipette is used. More strips (up to 6) could be applied, when a multistep pipette is used.
4. Return all reagents to 2 - 8 °C (35 - 46 °F) immediately after use.

The **aflatoxin standards** are provided ready to use. The dilution factor 10 for the sample has been considered when labelling. Therefore, the aflatoxin concentration of samples can be read directly from the standard curve.

## 10.2. Test procedure

Accurate washing is very important. Do not allow microwells to dry up totally and avoid prolonged intervals between the working steps. Reproducibility in any enzyme immunoassay is largely dependent upon the consistency with which the microwells are washed. Carefully follow the recommended washing sequence as outlined in the test procedure.

Avoid direct sunlight during all incubations. Covering the microtiter plates is recommended.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Pipet 50 µl of standard or prepared sample into separate wells; use a new pipette tip for each standard or sample.
3. Add 50 µl of enzyme conjugate (red cap) to each well.
4. Add 50 µl of anti-aflatoxin antibody solution (black cap) to each well. Mix gently by shaking the plate manually and incubate for 10 min (+/- 1) at room temperature (20 - 25 °C / 68 - 77 °F).
5. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Using a wash bottle or multichannel pipette, fill the wells (250 µl per well) with washing buffer (see 10.1.). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.
6. Add 100 µl of substrate/chromogen (brown cap) to each well. Mix gently by shaking the plate manually and incubate for 5 min (+/- 0.5) at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of stop solution (yellow cap) to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

## 11. Results

A special software, the RIDA<sup>®</sup>SOFT Win (Art. No. Z9999), is available to evaluate the RIDASCREEN<sup>®</sup> enzyme immunoassays.

For single determinations we recommend logit/log evaluation and for double or multiple determinations cubic spline should be used.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the aflatoxin concentration [ $\mu\text{g}/\text{kg}$ ].

The aflatoxin concentration in  $\mu\text{g}/\text{kg}$  corresponding to the extinction of each sample can be read from the calibration curve.

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG