

Colorimetric assay for the determination of iron in wine and food analysis
4 x 100 ml (200 assays)

For in vitro use only
Store between +2 and +8°C

Principle

The Iron is dissociated from proteins in particular conditions of ionic strength, and reduced to a bivalent state through ascorbic acid. It will then react with the chromogenic reagent FERENE-S and build a stable complex, giving a blue colour which is proportional to the concentration of Iron in the sample.

Assay specifications

Wavelength: 582 nm (575 – 582 nm)
Light path: 1.00 cm (glass; plastic)
Temperature: 25 to 37°C
Method: end point
Reaction: 15 minutes
Measurement: against air or against water
Linearity: 2 – 40 mg/L

Reagents

- # 1: Buffer, 4 vials with approx. 84 ml each (Buffer > 0.1 mol/L)
- # 2: Ferene-S, 4 vials with approx. 16 ml each (Ferene-S > 0.1 mmol/L; Ascorbic acid > 0.1 mol/L).
- # 3: Calibrator, 1 x 5 ml (Iron = 20 mg/l).

All reagents are ready for use. They are stable at 2-8 °C up to the expiry date shown on the package if not contaminated during handling.

Let the reagents reach the laboratory temperature (20 – 25 °C) before use. Mix kindly before pipeting. Close immediately after handling.

The reagents are not hazardous. The general safety rules for the work in chemical laboratories should be applied. After use, the reagents can be disposed of with the laboratory waste. Packaging materials may be recycled.

Sample preparation

- Wine can be used directly.
- Use colorless, clear and neutral liquid samples directly if Iron conc. is between 2 – 40 mg/L; otherwise, dilute with water to reduce it in this range.
- Strongly coloured samples have to be treated with PVPP (polyvinylpyrrolidone e.g. 1 g/100 mL sample)
- For application on biochemistry analysers, it is recommended to add PVP (polyvinylpyrrolidone) at a final concentration of 5 g/l into R1 (2.1 ml of a mother solution 200 g/l in each vial)
- Turbid solutions have to be filtered or centrifuged
- Samples containing carbon dioxide have to be degassed

Procedure

Pipette into cuvettes:	Reagent blank (RB)	Calibrator	Samples
Reagent 1 (buffer)	1680 µl	1680 µl	1680 µl
Distilled water	100 µl	-	-
Standard (vial 3)	-	100 µl	-
Sample	-	-	100 µl
Mix carefully. Read the absorbance A ₁ after 5 minutes at 25 to 37 °C, then add:			
Reagent 2 (Ferene-S)	320 µl	320 µl	320 µl
Mix carefully. Read the absorbance A ₂ after 10 min. at 25 to 37 °C. The color is stable 30 min. at room temperature.			

Calculation

$$\Delta A = (A_2 - df \times A_1)_{\text{sample or standard}} - (A_2 - df \times A_1)_{\text{RB}}$$

with df = dilution factor of the optical densities by reagent volumes:
df = (sample volume + R1) / (sample volume + R1 + R2) = 0.848

$$\text{and } C_{\text{sample}} [\text{mg/L}] = \frac{C_{\text{standard}} [\text{mg/L}]}{\Delta A_{\text{standard}}} \times \Delta A_{\text{sample}}$$

Since the concentration of the standard is fixed at 20 mg/L, this gives the following calculation formula depending on the unit that is preferred:

$$C_{\text{sample}} [\text{mg/L}] = (\Delta A_{\text{sample}} / \Delta A_{\text{standard}}) \times 20$$

$$\text{or } C_{\text{sample}} [\mu\text{mol/L}] = (\Delta A_{\text{sample}} / \Delta A_{\text{standard}}) \times 358$$

Notes

1. A proportional variation of the reaction volumes does not change the result.
2. Use one way cuvettes or very clean tubes washed with diluted HCl and distilled water.
3. Specificity: this test is specific for Iron, no interferences were detected.
4. The sample volume is low compared to other tests, because the Ferene-S chromogen is very sensitive; the reduced sample volume has the advantage of reducing interferences by the sample matrix
5. The test takes place under acidic pH conditions, so in many cases no pH adjustment of the sample is needed (for instance for wine samples or after perchloric acid treatment)

Additional application examples

- Determination of iron in flour:
 - weigh precisely 0.8 grams of flour
 - resuspend in 6 – 7 ml of 7% (w/v) perchloric acid
 - shake and mix for 10 min.
 - fill to the mark at 10 ml with 7 % (w/v) perchloric acid
 - mix a few times by inversion
 - centrifuge at 1,700 - 2,000 x g for 20 min; alternatively, filtration is also possible (e.g. Macherey Nagel ref. 531012, MN 615 1/4, diameter 125 mm)
 - use the supernatant solution as sample for Iron determination

The content in iron is calculated with the formula:

$$\text{Content [g/100 g]} = \frac{C_{\text{iron}} [\text{g/l}]}{\text{weight}_{\text{sample}} [\text{g/l}]} \times 100$$

(weight_{sample} = 80 g/L if initial weight is 0.8 g in 10 ml as in our example)

Literature

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Duffy J. R., Gaudin J., Clin. Biochem. 10, 122 (1977).
3. Higgins T., Clin. Chem. 27, 1619 (1981).

Applications on biochemistry analysers

General application features (all analyzers)	
Temperature	37°C (25°C is also possible)
Wavelength	575 - 600 nm (primary) / 700 nm (secondary)
Test type and sequence	End-point test with 2 reagents test sequence: - pipeting R1 (5 Vol.) + sample (20 – 50 µl) - pre-incubation 2 – 3 min - measurement A1 (before addition of R2) - pipeting R2 (1 Vol.) - incubation >= 5 min at 37 °C - measurement A2 - calculation A2-A1 against calibration curve
Calibration	2 - 5 calibrators from 0 to 10 mg/l Linear calibration curve
Reagent blank	Yes
Reagents	R1 and R2 in the ratio approx. 5 to 1 20 µl for low sensitivity (0 – 20 mg/l)
Sample volume	50 µl for high sensitivity (0 – 10 mg/l)

Example Konelab / Arena			
Test Definition:			
Full name	Enzytec Iron		
Online name	#		
Test used	yes		
Test type	Photometric		
Result unit	mg/l		
Number of decimals	1		
Temperature	37 °C		
Test limit	Low 0	High 20	Unit mg/l
Initial absorbance	0.000	2.000	E
Dilution limit	0	20	mg/l
Secondary dilution 1 +	#	#	
Acceptance	Automatic		
Dilution 1 +	0.0		
Sample type	Wine (Other)		
Correction factor	1.00		
Correction bias	0.00		
Calibration parameter			
Calibration type	Linear		
Curve direction	Ascending		
Repeat time (d)	1	Err. (mA)	#
Points/calibrator	simple	Err. (%)	#
Acceptance	Automatic		
Response limit min/max	# / #		
Bias correction in use	No		
Type of calibrators	Series (separate is also possible)		
Calibrator identification	Iron calibrator (20 mg/l)		
Concentrations	0 to 10 mg/l		
Distribution			
Additional blank	YES (or fixed time)	Antigen excess	No
Reagent 1 volume (µl)	ENZYT Iron R1		
dispense with wash reagent	150	Volume (µl)	10
Sample volume (µl)	Water	(adjust to sensitivity needs)	
dispense with wash reagent	None	Volume (µl)	10
Incubation (sec)	180		
Measurement BLANK response min./max	End-point		
Reagent 2 volume (µl)	# / #		
wash reagent dispense with	ENZYT Iron R2		
Incubation (sec)	25	Volume (µl)	10
Measurement	None		
λ 1 / λ 2 (nm)	Water		
measurement type	300		
	End-point		
	600 / 700		575 / 700 is also possible
	Fixed timing		

Example Lisa 200		
Full Test name	Enzytec Iron Ferene	
Short name	Iron	
Units	mg/l	
Type of test	End-point calibration	
Filter (nm)	580	
1st reading = zero	NO (A1 must be measured for each sample)	
Waiting time 1 (tours before addition of R2)	4 (4 tours = 4 x 24 sec.)	
Number of measurements (= incubation time)	15 (15 tours x 24 sec = 6 min)	
Reagent 1	VOL (µl)	250
	DIL (µl)	0
	POS	#
Reagent 2	VOL (µl)	45
	DIL (µl)	0
	POS	#
Sample.	VOL (µl)	20 - 50 µl (adjust to sensitivity needs)
	DIL	0
Starter (of reaction)	REAGENT 2	
Waiting time 2 (tours after addition of R2)	0 (measurements start immediately after addition of R2)	
Calibration	1 degree (= linear calibration; 2 degrees = non-linear calibration is also possible)	
Blank = Standard	YES (water used for blank is used as 1 st point of the calibration curve)	
Number of standards	3	
Standard 1	VAL (mg/l)	2 (= example with 2 mg/l)
	POS	# (position on the instrument)
Standard 2	VAL (mg/l)	5 (= example with 5 mg/l)
	POS	#
Standard 3	VAL (mg/l)	10 (= example with 10 mg/l)
	POS	#
Nb repeat Sample/Control	1 (number of repetitions ; it is possible to choose 2)	
Control	VAL (mg/l)	# (concentration of QC sample)
	POS	# (position on the instrument)
	DEV	# (deviation accepted)
Pre-dilution	1	
Post-dilution	4	
Sample diluents	Physiologic solution	
Type rincing	3	
Normal values	UPPER	20 (in mg/l)
	LOWER	0
Limit linearity	20 (in mg/l)	
Blanc OD limit	LOWER	0
	UPPER	500

value entered by the operator

Example Hitachi 717	
TEST	[IRON]
ASSAY CODE	[2POINTS] - [24] - [50]
SAMPLE VOLUME (µl)	[30] - [5] (adjust to sensitivity needs)
REAGENT VOLUME R1	[250] - [50] - [NO]
REAGENT VOLUME R2	[40] - [20] - [NO] (can be lower)
WAVELENGTH	[700] [600]
CALIBRATION TYPE	[LIN] [1] [4]
STD.(1) CONC. POS.	[0] - [1] (example)
STD.(2) CONC. POS.	[2] - [2] (2 mg/l =example)
STD.(3) CONC. POS.	[5] - [3] (5 mg/l =example)
STD.(4) CONC. POS.	[10] - [4] (10 mg/l =example)
SD LIMIT	[0.1]
DUPLICATE LIMIT	[#]
SENSITIVITY LIMIT	0
LIMITE D.O (CINET)	[0] - [INCR.]
PROZONE CHECK	[32000] - [1]
EXPECTED VALUE	[#] - [#]
PANIC VALUE	[#] - [#]
INSTRUMENT FACTOR	[1.0]