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Esempio Konelab / Arena (software inglese)			
Test Definition:			
Full name	Enzytec Iron		
Online name	#		
Test used	yes		
Test type	Photometric		
Result unit	mg/l		
Number of decimals	1		
Temperature	37 °C		
	Low	High	Unit
Test limit	0	20	mg/l
Initial absorbance	0.000	2.000	E
Dilution limit	0	20	mg/l
Secondary dilution 1 +	#	#	
Acceptance	Automatic		
Dilution 1 +	0.0		
Sample type	Wine (Other)		
Correction factor	1.00		
Correction bias	0.00		
Calibration parameter			
Calibration type	Linear		
Curve direction	Ascending		
Repeat time (d)	1	Err. (mA)	#
Points/calibrator	simple	Err. (%)	#
Acceptance	Automatic		
Response limit min/max	# / #		
Bias correction in use	No		
Type of calibrators	Series	(separate is also possible)	
Calibrator identification	Iron calibrator (20 mg/l)		
Concentrations	0 to 10 mg/l		
Distribution			
Additional blank	YES (or fixed time)	Antigen excess	No
Reagent 1 volume (µl)	ENZYT Iron R1		
dispense with wash reagent	Water	Volume (µl)	10
Sample volume (µl)	15	(adjust to sensitivity needs)	
dispense with wash reagent	Water	Volume (µl)	10
Incubation (sec)	None		
Measurement BLANK response min./max	180		
Reagent 2 volume (µl)	End-point		
wash reagent	# / #		
dispense with	ENZYT Iron R2		
Incubation (sec)	25		
Measurement	None		
λ 1 / λ 2 (nm)	Water	Volume (µl)	10
measurement type	300		
	End-point		
	600 / 700	575 / 700 is also possible	
	Fixed timing		

Esempio Lisa 200 (software inglese)	
Full Test name	Enzytec Iron Ferene
Short name	Iron
Units	mg/l
Type of test	End-point calibration
Filter (nm)	580
1st reading = zero	NO (A1 must be measured for each sample)
Waiting time 1 (tours before addition of R2)	4 (4 tours = 4 x 24 sec.)
Number of measurements (= incubation time)	15 (15 tours x 24 sec = 6 min)
Reagent 1 VOL (µl)	250
DIL (µl)	0
POS	#
Reagent 2 VOL (µl)	45
DIL (µl)	0
POS	#
Sample. VOL (µl)	20 - 50 µl (adjust to sensitivity needs)
DIL	0
Starter (of reaction)	REAGENT 2
Waiting time 2 (tours after addition of R2)	0 (measurements start immediately after addition of R2)
Calibration	1 degree (= linear calibration; 2 degrees = non-linear calibration is also possible)
Blank = Standard	YES (water used for blank is used as 1 st point of the calibration curve)
Number of standards	3
Standard 1 VAL (mg/l)	2 (= example with 2 mg/l)
POS	# (position on the instrument)
Standard 2 VAL (mg/l)	5 (= example with 5 mg/l)
POS	#
Standard 3 VAL (mg/l)	10 (= example with 10 mg/l)
POS	#
Nb repeat Sample/Control	1 (number of repetitions ; it is possible to choose 2)
Control VAL (mg/l)	# (concentration of QC sample)
POS	# (position on the instrument)
DEV	# (deviation accepted)
Pre-dilution	1
Post-dilution	4
Sample diluents	Physiologic solution
Type rincing	3
Normal values UPPER	20 (in mg/l)
LOWER	0
Limit linearity	20 (in mg/l)
Blanc OD limit LOWER	0
UPPER	500

value entered by the operator

Esempio Hitachi 717 (software inglese)	
TEST	[IRON]
ASSAY CODE	[2POINTS] - [24] - [50]
SAMPLE VOLUME (µl)	[30] - [5] (adjust to sensitivity needs)
REAGENT VOLUME R1	[250] - [50] - [NO]
REAGENT VOLUME R2	[40] - [20] - [NO] (can be lower)
WAVELENGTH	[700] [600]
CALIBRATION TYPE	[LIN] [1] [4]
STD.(1) CONC. POS.	[0] - [1] (example)
STD.(2) CONC. POS.	[2] - [2] (2 mg/l =example)
STD.(3) CONC. POS.	[5] - [3] (5 mg/l =example)
STD.(4) CONC. POS.	[10] - [4] (10 mg/l =example)
SD LIMIT	[0.1]
DUPLICATE LIMIT	[#]
SENSITIVITY LIMIT	0
LIMITE D.O (CINET)	[0] - [INCR.]
PROZONE CHECK	[32000] - [1]
EXPECTED VALUE	[#] - [#]
PANIC VALUE	[#] - [#]
INSTRUMENT FACTOR	[1.0]