RIDASCREEN[®] FAST Aflatoxin M₁

Enzymimmunoassay zur quantitativen Bestimmung von
Aflatoxin M₁

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RIDASCREEN®FAST Aflatoxin M₁

Brief information

RIDASCREEN®FAST Aflatoxin M₁ (Art. No.: R5812) is a competitive enzyme immunoassay for the quantitative analysis of aflatoxin M₁ in milk and milk powder. All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Technicians need no specialized training to perform the RIDASCREEN®FAST Aflatoxin M₁ test, however, free support is offered by the distributor on request, if necessary.

Sample preparation: centrifugation

Time requirement: sample preparation (for 10 samples)
milk and milk powder approx. 15 min
test implementation (incubation time) 15 min

Limit of detection: milk: < 125 ppt
(corresponding to the milk powder (referring to reconstituted milk): < 125 ppt
standard substance)

1. Intended use

RIDASCREEN®FAST Aflatoxin M₁ is a competitive enzyme immunoassay for the quantitative analysis of aflatoxin M₁ in milk and milk powder.

2. General

Aflatoxins are secondary metabolites of the fungi species *Aspergillus flavus*, *parasiticus* and *nomius*, which are carcinogenic and highly toxic.

Aflatoxin M₁ is produced as a metabolite of aflatoxin B₁. It is secreted with the milk after the feeding of aflatoxin B₁ content feedstuffs to lactating cows.

As aflatoxin M₁ is relatively stable towards the pasteurizing process, not only a comprehensive routine check of the raw materials to be processed is required, but also of the final products.

Since the first of January 1999 EU-wide uniform residue limits for aflatoxins exist. For aflatoxin M₁ the limit has been fixed at 0.05 µg/l (50 ppt).

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-aflatoxin M₁ antibodies.

Aflatoxin M₁ standards or sample solutions, aflatoxin M₁ enzyme conjugate and anti-aflatoxin M₁ antibodies are added. Free aflatoxin M₁ and aflatoxin M₁ enzyme conjugate compete for the aflatoxin M₁ antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-aflatoxin M₁ antibodies are also bound by the immobilized capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the aflatoxin M₁ concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for as many as 42 analyses (plus 6 standard analyses). Each test kit contains:

- 1 x Microtiter plate with 48 wells (6 strips with 8 removable wells each)
coated with capture antibodies
- 6 x Standard solutions, 1.3 ml each
0 ppt (zero standard), 125 ppt, 250 ppt, 500 ppt, 1000 ppt, 2000 ppt
aflatoxin M₁ in milk buffer
ready to use
- 1 x Conjugate (3 ml)red cap
peroxidase conjugated aflatoxin
ready to use
- 1 x Anti-aflatoxin M₁ antibody (3 ml).....black cap
ready to use
- 1 x Red Chromogen Pro (10 ml) brown cap
Substrate/chromogen solution, stained red
contains tetramethylbenzidine
- 1 x Stop solution (14 ml)..... yellow cap
contains 1 N sulfuric acid
- 1 x Washing Buffer (Salt)
for preparation of a 10 mM phosphate buffer (pH 7.4)
contains 0.05 % Tween 20

5. Materials required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge
- graduated pipettes
- variable 20 - 200 µl and 200 - 1000 µl micropipettes

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The standards contain aflatoxin M₁, particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

Decontamination of the glassware and aflatoxin solutions is best carried out using a sodium hypochlorite (bleach) solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

The stop solution contains 1 N sulfuric acid (R36/38, S2-26).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen solution is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

The test kit can be regularly used at least up to the expiry date (indicated on the kit package), if stored correctly.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the reddish substrate/chromogen solution prior to test implementation
- a value of less than 0.6 absorbance units ($A_{450\text{ nm}} < 0.6$) for the zero standard

9. Preparation of Samples

The samples should be stored in a cool place, protected from light.

9.1. Milk

- centrifuge milk samples for degreasing: 10 min / 3500 g / 10 °C (50 °F) (if a refrigerated centrifuge is not available, cool sample to 10 °C (50 °F) prior to centrifugation)
- after centrifugation, remove upper cream layer completely by aspirating through a Pasteur pipette
- use the skimmed milk (= defatted supernatant) directly in the test (50 µl per well)

9.2. Milk powder

- weigh 10 g of milk powder in a flask and fill up to 100 ml with deionized water *)
- dissolve by stirring for 5 min
- continue with the preparation of milk as described in 9.1.

*) sample size may be increased if required, but the volume of water must be adapted accordingly

10. Test implementation

10.1. Preliminary comments

1. Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.
2. The specific reaction starts with the addition of the specific antibody. Nevertheless, not more than three strips should be applied to the test, when a single step pipette is used. More strips (up to 6) could be applied, when a multistep pipette is used.
3. Return all reagents to 2 - 8 °C (35 - 46 °F) immediately after use.

As **washing buffer** a PBS tween buffer is needed. Please use the washing buffer salt (see 4.) contained in the kit. Dissolve the entire buffer salt in one liter of distilled water. The ready to use washing buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternatively: Dissolve the contents of the pouch in 100 ml of distilled water to obtain a 10fold concentrated washing buffer. Use 1 part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use washing buffer.
This 10fold concentrate expires after approx. 8 - 12 weeks, stored at room temperature (20 - 25 °C / 68 - 77 °F).

10.2. Test procedure

Accurate washing is very important. Do not allow microwells to dry up totally and avoid prolonged intervals between the working steps. Reproducibility in any EIA is largely dependent upon the consistency with which the microwells are washed. Carefully follow the recommended washing sequence as outlined in the EIA test procedure.

Avoid direct sunlight during all incubations. Therefore covering the microtiter plates is recommended.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Pipet 50 µl of standard or prepared sample into separate wells; use a new pipette tip for each standard or sample.
3. Add 50 µl of enzyme conjugate (red cap) to each well.
4. Add 50 µl of anti-aflatoxin M₁ antibody solution (black cap) to each well. Mix gently by shaking the plate manually and incubate for 10 min (+/- 1) at room temperature (20 - 25 °C / 68 - 77 °F).
5. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Using a wash bottle or multichannel pipette, fill the wells (250 µl per well) with washing buffer (see 10.1.). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.

6. Add 100 µl of substrate/chromogen (brown cap) to each well. Mix gently by shaking the plate manually and incubate for 5 min (+/- 0.5) at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of stop solution (yellow cap) to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes.

11. Results

A special software, the RIDA®SOFT Win (Art. No. Z9999), is available to evaluate the RIDASCREEN® enzyme immunoassays.

For single determinations we recommend logit/log evaluation and for double or multiple determinations cubic spline should be used.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the aflatoxin M₁ concentration [ng/l].

In order to obtain the aflatoxin M₁ concentration in ng/l (ppt) actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factors are as follows:

Milk.....	1
Milk powder (referring to reconstituted milk)	1
Milk powder (referring to g-weight)	10

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG