RIDASCREEN[®] FAST T-2 Toxin

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RIDASCREEN®FAST T-2 Toxin

Brief information

RIDASCREEN®FAST T-2 Toxin (Art. No.: R5302) is a competitive enzyme immunoassay for the quantitative analysis of T-2 Toxin in cereals and feed.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Technicians need no specialized training to perform the RIDASCREEN®FAST T-2 Toxin test, however, free support is offered by the distributor on request, if necessary.

The test was validated with the matrices corn, corn seedlings and mixed feed for chicken breeding and for fattened pigs.

Sample preparation: extraction, filtration and dilution

Time requirement: sample preparation (for 10 samples)
cereals and feed.....approx. 10 min
test implementation (incubation time).....15 min

Limit of detection: < 20 µg/kg (ppb)
(corresponding to the standard substance)

Limit of quantification: 50 µg/kg (ppb)
(corresponding to the standard substance)

1. Intended use

RIDASCREEN®FAST T-2 Toxin is a competitive enzyme immunoassay for the quantitative analysis of T-2 toxin in cereals and feed.

2. General

T-2 toxin belongs to the trichothecene group of mycotoxins and is produced by fungi of the genus *Fusarium*. T-2 toxin is often found in agricultural commodities, although the incidence and the concentrations found show a broad regional variation. Due to its cytotoxic and immunosuppressive activity, T-2 toxin is a threat for human and animal health.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-T-2 toxin antibodies. T-2 toxin standards or sample solution, T-2 toxin enzyme conjugate and anti-T-2 toxin antibodies are added. Free T-2 toxin and T-2 toxin enzyme conjugate compete for the T-2 toxin antibody binding sites (competitive enzyme immunoassay). At the same time, the T-2 toxin antibodies are also bound by the immobilized capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the T-2 toxin concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for as many as 43 analyses (plus 5 standard analyses). Each test kit contains:

- 1 x Microtiter plate with 48 wells (6 strips with 8 removable wells each)
coated with capture antibodies
- 5 x T-2 Toxin standard solutions ^{*}), 1.3 ml each
0 ppb (zero standard), 50 ppb, 100 ppb, 200 ppb, 400 ppb
T-2 Toxin in methanol/water, ready to use
- 1 x Conjugate (3 ml)red cap
peroxidase conjugated T-2 toxin
ready to use
- 1 x Anti-T-2 Toxin antibody (3 ml)black cap
monoclonal, ready to use
- 1 x Substrate/Chromogen solution (10 ml) brown cap
stained red
- 1 x Stop solution (14 ml)..... yellow cap
contains 1 N sulfuric acid

^{*}) The dilution factor 10 for the sample has already been considered. Therefore, the T-2 Toxin concentrations of samples can be read directly from the standard curve.

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- graduated cylinder (plastic or glass) 100 ml
- glassware for preparing sample extract: filter funnel and 50 ml flask
- grinder (mill)
- optional: shaker
- filter paper: Whatman No. 1 or equivalent
- 50 µl, 100 µl and 1000 µl micropipettes

5.2. Reagents:

- methanol
- 70 % methanol solution: prepare 70 % methanol solution by mixing 70 ml pure methanol with 30 ml distilled or deionized water
- distilled or deionized water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The standards contain T-2 Toxin, particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

Decontamination of the glassware and T-2 Toxin solutions is best carried out using a sodium hypochlorite solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

The stop solution contains 1 N sulfuric acid (R36/38, S2-26).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen solution is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

The test kit can be regularly used at least up to the expiry date (indicated on the kit package), if stored correctly.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the reddish substrate/chromogen solution prior to test implementation
- a value of less than 0.6 absorbance units ($A_{450\text{ nm}} < 0.6$) for the zero standard

9. Preparation of Samples

The samples should be stored in a cool place, protected from light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- weigh 5 g of ground sample and add it to a suitable container with 25 ml of methanol (70 %) *)
- shake vigorously for 3 min (manually or with shaker)
- filter the extract through Whatman No. 1 filter
- dilute 1 ml of the obtained filtrate with 1 ml of distilled or deionized water
- use 50 µl of the filtrate per well in the test

*) sample size may be increased if required, but the volume of methanol/water must be adapted accordingly

10. Test implementation

10.1. Preliminary comments

1. Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.
2. The specific reaction starts with the addition of the specific antibody. Nevertheless, not more than three strips should be applied to the test, when a single step pipette is used. More strips (up to 6) could be applied, when a multistep pipette is used.
3. Return all reagents to 2 - 8 °C (36 - 46 °F) immediately after use.

The **T-2 Toxin standards** are provided ready to use. The dilution factor 10 for the sample has been considered when labeling. Therefore, the T-2 Toxin concentration of samples can be read directly from the standard curve.

10.2. Test procedure

Accurate washing is very important. Do not allow microwells to dry up totally and avoid prolonged intervals between the working steps. Reproducibility in any EIA is largely dependent upon the consistency with which the microwells are washed. Carefully follow the recommended washing sequence as outlined in the EIA test procedure.

Avoid direct sunlight during all incubations. Covering the microtiter plates is recommended.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Pipet 50 µl of standard or prepared sample to separate wells; use a new pipette tip for each standard or sample.
3. Add 50 µl of enzyme conjugate (red cap) to the bottom of each well.
4. Add 50 µl of the anti-T-2 toxin antibody solution (black cap) to each well. Mix gently by shaking the plate manually and incubate for 10 min (+/- 1) at room temperature (20 - 25 °C / 68 - 77 °F).
5. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Using a wash bottle or multichannel pipette, fill the wells with distilled or deionized water (250 µl per well). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.
6. Add 100 µl of substrate/chromogen (brown cap) to each well. Mix gently by shaking the plate manually and incubate for 5 min (+/- 0.5) at room temperature (20 - 25 °C; 68 - 77 °F) in the dark.
7. Add 100 µl of stop solution (yellow cap) to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win / RIDA[®]SOFT Win.net (Art. Nr. R9996), is available to evaluate the RIDASCREEN[®] enzyme immunoassays.

For single determinations we recommend logit/log evaluation and for double or multiple determinations cubic spline should be used.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the T-2 Toxin concentration [$\mu\text{g}/\text{kg}$]. The T-2 Toxin concentration in $\mu\text{g}/\text{kg}$ corresponding to the extinction of each sample can be read from the calibration curve.

12. Sensitivity

12.1. Limit of detection

The limit of detection (LOD) was investigated with cereal samples and mixed feed by repeated measurement of the zero-matrix. LOD was defined as the concentration that corresponds to the mean absorbance of the measurement - 3 times the standard deviation, and was determined by extrapolation from the standard curve. All results were found between 10 and 20 ppb.

12.2. Limit of quantification

The limit of quantification (LOQ) was investigated by repeated spiking experiments at the lower end of the standard curve; LOQ was found at 50 $\mu\text{g}/\text{kg}$ (ppb), which concentration is well above LOD, corresponding recoveries were found between 70 and 110 %, with CV-values below 20 %.

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG