**User Manual**

# **SureFast<sup>®</sup> PREP Aqua**

## **Art. No. F1023**

**100 extractions**

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## Description

This kit is intended to be used for the isolation of bacteria DNA from water samples.

## Principle

1. Membrane filtration
2. Lysis of bacterial cells
3. Setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

## Kit components (per box) and storage

50x Membrane Filter (white with grid; pore size 0.45 µm) <sup>1)</sup>	<b>(Code F)</b>	1x Elution Buffer (10 ml)	<b>(Code E)</b>
		1x Receiver Tubes 2.0 ml (50x)	<b>(Code R)</b>
1x Lysis Buffer A (50 ml)	<b>(Code L)</b>	1x Receiver Tubes 1.5 ml (50x)	<b>(Code T)</b>
1x Binding Buffer (10 ml)	<b>(Code B)</b>	1x Spin Filter (50x)	<b>(Code S)</b>
1x Wash Buffer (60 ml) <sup>2)</sup>	<b>(Code W)</b>		

<sup>1)</sup> The membrane included in the kit is made of mixed cellulose esters. As an alternative of the membrane filter (included in the kit) polycarbonate membrane filter can be used for the detection of *Legionella* spp. according to ISO/TS 12869.

<sup>2)</sup> After adding ethanol (purity ≥ 96 %; not supplied with the kit)

All reagents of the kit should be stored dry and at room temperature (14 – 25°C).

## Additionally required equipment and materials

- membrane filtering system
- reaction tubes free from DNA and DNase 1.5 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- sterile tweezers and spatula
- Vortex mixer
- Thermomixer/ heating block (up to 95°C)
- micro centrifuge (up to 12,000 rpm)
- ethanol (purity ≥ 96 %) for preparation of Wash Buffer
- disposal bags or waste bin
- autoclave
- biological safety cabinet

## Preparations

### General

- Add 42 ml ethanol to the Wash Buffer **(Code W)** and mix thoroughly.

### Before each preparation

- Preheating the Elution Buffer **(Code E)** - Transfer the needed amount of Elution Buffer **(Code E)** under calculation of a reserve volume into a reaction tube (not supplied with the kit) and equilibrate to 60°C (the Elution Buffer is necessary for step 7).

**Extraction protocol****1. Preparation of the basic material**

Filtrate 500 ml (0.1 – 1.0 l depending on the degree of contamination) of the water sample using the Membrane Filter (**Code F**). Place the filter (white with grid) with the grid showing upwards onto the filter holder.

**Note:** If the SureFast® PREP Aqua kit is used together with the SureFast® Legionella Screen PLUS real time PCR kit for the detection of *Legionella* spp. in combination with the German drinking water ordinance 2001 limit of 100 CFU/100 ml (Positive Control), 100 ml of the sample have to be filtered.

Fold the filter with a sterile tweezer (because the filter needs to be folded more than once, a sterile spatula to fix the filter might be helpful) and transfer it into a 1.5 ml reaction tube (not supplied with the kit).

**2. Thermal lysis**

Add 1.0 ml Lysis Buffer (**Code L**) to the filter, mix briefly to ensure that the whole filter is completely surrounded by Lysis Buffer. Close the reaction tube with a cap and incubate under continuously shaking for 10 min at 95°C.

**3. Setting of optimal binding conditions**

Transfer 130 µl of the lysate into a 1.5 ml reaction tube (not supplied with the kit) and add 200 µl Binding Buffer (**Code B**).

**4. Binding of the nucleic acids on a Spin Filter**

Place a Spin Filter (**Code S**) into a yellow 2.0 ml Receiver Tube (**Code R**). Mix the sample solution, transfer it directly onto the filter and incubate at room temperature for 1 min.

Centrifuge for 1 min at 12,000 rpm. Discard the filtrate and place the Spin Filter back into the yellow Receiver Tube.

**5. Purification of the bound nucleic acids**

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge for 1 min at 12,000 rpm. Discard the filtrate and place the Spin Filter back into the yellow Receiver Tube.

Once more add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge for 1 min at 12,000 rpm. Discard the filtrate and place the Spin Filter back into the yellow Receiver Tube.

**6. Drying of the Spin Filter**

Remove the residual ethanol by final centrifugation for 2 min at 12,000 rpm.

**7. Elution of nucleic acids from the Spin Filter**

Place the Spin Filter into a clear 1.5 ml Receiver Tube (**Code T**) and add 50 µl of the preheated Elution Buffer (**Code E**) directly onto the Spin Filter.

Incubate 3 min at room temperature and centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

The eluted nucleic acids are ready-to-use for the PCR. The DNA can be stored for up to 24 hours at 4°C. For a storage time of more than 24 hours it should be kept at -20°C.

**Product Information**

- Material Safety Data Sheet

## Technical Support

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).

## Safety information

### Binding Buffer

Danger

H225-319-336 P210-233-305-351-338

- H225:** Highly flammable liquid and vapour.  
**H319:** Causes serious eye irritation.  
**H336:** May cause drowsiness or dizziness.  
**P210:** Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
**P233:** Keep container tightly closed.  
**P305+P351+P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For further information we offer a Material Safety Data Sheet. Please contact your distributor or send an e-mail to [info@congen.de](mailto:info@congen.de).