RIDASCREEN® Sulfamethazin

REF R3011

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RIDASCREEN® Sulfamethazin

Brief information

RIDASCREEN® Sulfamethazin (R3011) is a competitive enzyme immunoassay for the quantitative analysis of sulfamethazine in milk, meat, honey, liver, kidney, fish, shrimp and egg.

All reagents required for performing the enzyme immunoassay, including the standards, are contained in the test kit. The test kit is sufficient for a maximum of 48 duplicate determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: Milk: direct use
 Meat: homogenization, centrifugation
 Honey: extraction, clean-up with RIDA® C18 column, evaporation, reconstitution, precipitation
 Liver, kidney, fish: homogenization, centrifugation, precipitation
 Shrimp, egg: homogenization, extraction, evaporation, reconstitution

Time requirement: sample preparation (for 10 samples)
 according to sample matrix 5 min to 2 h
 test implementation (incubation time) 45 min

Limit of detection: Milk approx. 4 µg/l
(corresponding to the Meat (bovine, porcine) approx. 5 µg/kg
standard substance) Meat (poultry) approx. ca. 10 µg/kg
 Honey* approx. ca. 10 µg/kg
 Liver approx. ca. 6 µg/kg
 Kidney approx. ca. 10 µg/kg
 Fish approx. ca. 7 µg/kg
 Shrimp approx. ca. 15 µg/kg
 Egg approx. ca. 16 µg/kg

* For honey the LOD can deviate due to matrix effects.

| | | |
|---|------------------------|-----------------------|
| Detection capability (CC β): (corresponding to the standard substance) | Milk | approx. 10 μ g/l |
| | Meat (bovine, porcine) | approx. 10 μ g/kg |
| | Meat (poultry) | approx. 15 μ g/kg |
| | Honey | approx. 20 μ g/kg |
| | Liver | approx. 15 μ g/kg |
| | Kidney | approx. 15 μ g/kg |
| | Fish | approx. 20 μ g/kg |
| | Shrimp | approx. 20 μ g/kg |
| | Egg | approx. 15 μ g/kg |

| | | |
|--|------------------------|---------------|
| Recovery rate: (corresponding to the standard substance) | Milk | approx. 110 % |
| | Meat (bovine, porcine) | approx. 94 % |
| | Meat (poultry) | approx. 115 % |
| | Honey | approx. 88 % |
| | Liver | approx. 89 % |
| | Kidney | approx. 96 % |
| | Fish | approx. 88 % |
| | Shrimp | approx. 111 % |
| | Egg | approx. 107 % |

The specificity of the RIDASCREEN[®] Sulfamethazine assay was determined by analyzing the cross-reactivity to the respective substances in the buffer system. In samples, the specificity may differ from the values determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the limit of detection and the recovery rate in the relevant matrix. The test cannot distinguish between analytes and cross-reactive substances.

| | | |
|----------------------------------|---|-------------|
| Specificity: in buffer system | Sulfamethazine (standard substance) | 100 % |
| | Sulfamerazine | approx. 4 % |
| | Sulfacetamide, Sulfachloropyridazine, Sulfachloropyridazine-(phenyl- ¹³ C6), Sulfadiazine, Sulfadimethoxine, Sulfadoxine, Sulfaguandine, Sulfamethoxazole, Sulfamethoxypyridazine, Sulfamethizole, Sulfamoxole, Sulfanilamide, Sulphaphenazole, Sulfapyridine, Sulfaquinoxaline, Sulfisoxazole | < 1 % |

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice manual. It lists minimum standards concerning the framework conditions when using test kits of R-Biopharm AG and performing ELISA analyses with them. The manual can be retrieved, printed and downloaded from www.r-biopharm.com/products/food-feed-analysis.

Product offer

RIDA® Sulfamethazin Spiking Solution (R3098)

RIDASCREEN® Sulfonamide (R3004)

RIDA® Sulfonamide-Sulfamethoxypyridazin Spiking Solution (R3099)

1. Intended use

RIDASCREEN® Sulfamethazin (R3011) is a competitive enzyme immunoassay for the quantitative analysis of sulfamethazine in milk, meat, honey, liver, kidney, fish, shrimp and egg.

2. General information

Sulfonamides are widely used as feed additives, mainly for fattening of calves and pigs. Combined with inhibitors of dihydrofolate reductase such as trimethoprim, tetroxoprim, or pyrimethamine sulfonamides are also used in veterinary medicine for the treatment of intestinal infections, mastitis, pneumonitis and other (systemic) diseases. Sulfonamide residues may therefore occur in food of animal origin such as meat and milk. Particularly the transmission of the carcinogenic sulfamethazine represents a threat to human health. The EU established a maximum residue limit (MRL) for all substances of the sulfonamide-group of 100 µg/kg (ppb) in muscle, fat, liver, kidney and milk.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-sulfamethazine antibodies. The microtiter plate is already functionalized with anti-sulfamethazine antibodies, which are bound by the immobilized capture antibodies. Standards or sample and sulfamethazine conjugate are added into the wells. Free and enzyme conjugated sulfamethazine competes for the antibody binding sites (competitive enzyme immunoassay). Any unbound conjugate is then removed in a washing step. The conjugate converts the chromogen into a blue product. The addition of the stop

solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the sulfamethazine concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for a max. 96 duplicate determinations (including standard analyses). Each test kit contains:

| Component | Cap color | Format | Volume |
|---|-----------|---------------------|--------------------|
| Microtiter plate | - | Ready-to-use | 96 wells |
| Standard 1 | White | Ready-to-use | 0 µg/l 1.3 ml |
| Standard 2 | White | Ready-to-use | 3 µg/l 1.3 ml |
| Standard 3 | White | Ready-to-use | 9 µg/l 1.3 ml |
| Standard 4 | White | Ready-to-use | 27 µg/l 1.3 ml |
| Standard 5 | White | Ready-to-use | 81 µg/l 1.3 ml |
| Standard 6 | White | Ready-to-use | 162 µg/l 1.3 ml |
| Wash buffer salt Tween | - | Salt for dissolving | |
| Conjugate | Red | Ready-to-use | 6 ml |
| Substrate/Chromogen Red Chromogen Pro | Brown | Ready-to-use | 10 ml |
| Stop solution | Yellow | Ready-to-use | 14 ml |

5. Reagents required but not provided

5.1 Equipment

| Equipment | milk | meat | honey | liver | kidney | fish | shrimp egg |
|---|------|------|-------|-------|--------|------|------------|
| microtiter plate spectrophotometer (450 nm) | • | • | • | • | • | • | • |
| graduated pipettes | | • | • | • | • | • | • |
| variable 5 - 10 µl micropipettes | | | • | | • | • | |
| variable 20 - 200 µl and 200 - 1000 µl micropipettes | • | • | • | • | • | • | • |
| mixer | | • | | • | • | • | • |
| shaker | | • | | • | • | • | • |
| centrifuge | | • | • | • | • | • | • |
| evaporator | | | • | | | | • |
| vortex | • | • | • | • | • | • | • |
| RIDA® C-18 column | | | • | | | | |

5.2 Reagents

| Reagent | milk | meat | honey | liver | kidney | fish | shrimp egg |
|-------------------|------|------|-------|-------|--------|------|------------|
| extraction buffer | | | • | | | | |
| methanol | | | • | | | | • |
| Carrez I | | | • | • | • | • | |
| Carrez II | | | • | • | • | • | |

Extraction buffer for honey: 50 mM sodium acetate buffer, pH 5

- Dissolve 4.1 g sodium acetate (CH_3COONa) in 900 ml distilled water, adjust pH to 5 with HCl and fill up to 1000 ml with distilled water

Carrez I:

- Carrez I: 1.52 g Potassium ferrocyanide(II) x 3 H_2O ; fill up to 10 ml with distilled water

Carrez II:

- Carrez II: 2.99 g Zinc sulfate x 7 H_2O ; fill up to 10 ml with distilled water

6. Warnings and precautions for the users

Only trained laboratory employees should carry out this test. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Ensure the proper and responsible disposal of all reagents and materials after their use. For disposal, please adhere to national regulations.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The red-colored substrate/chromogen is photosensitive, therefore avoid direct exposure to light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the red-stained substrate/chromogen prior to test implementation
- A value of less than 0.6 absorbance units ($A_{450\text{ nm}} < 0.6$) for the zero standard

9. Sample preparation

The samples should be stored in a cool place, protected against light.

9.1 Milk

- Vortex a representative sample amount for homogenization
- Use 50 µl milk per well in the assay

9.2 Meat (beef, pork, poultry)

- Remove fat from a representative amount of meat and homogenize the meat with an Ultra-Turrax or a mixer
- Add 2 ml washing buffer (see 10.1.) to 1 g of homogenized sample
- Vortex for 10 seconds and then shake vigorously upside down for 5 minutes
- Centrifuge: 10 min / 4000 g / room temperature (20 - 25 °C)
- If there is an upper fat layer, remove it completely (e.g. with a pasteur pipette)
- Use 50 µl of the supernatant per well in the assay

9.3 Honey

- Dissolve 1 g of honey in 6 ml of extraction buffer (see 5.2)
- Centrifuge: 10 min / 3,000 g / room temperature (20 - 25 °C)

Purify all of the sample extract with RIDA® C-18 column (Art. No. R2002):

- Wash the column with 2 ml methanol
- Equilibrate the column with 2 ml of extraction buffer
- Apply all of the sample extract and wash with 8 ml extraction buffer
- Completely remove any fluid from the column by flushing the column with N₂ or air
- Elute sample slowly with 1 ml of 100 % methanol (approx. 15 drops per second)

- Evaporate the eluate to dryness under an nitrogen or air stream at 60 °C
- Reconstitute the dry residue in 1 ml washing buffer
- Add 10 µl of Carrez I, vortex
- Add 10 µl of Carrez II, vortex
- Centrifuge: 10 min / 4,000 g / room temperature
- Use 50 µl of the supernatant per well in the assay

9.4 Liver

- Remove fat from a representative amount of sample and homogenize the liver with an Ultra-Turrax or a mixer
- Add 2 ml washing buffer (see 10.1.) to 1 g of homogenized sample
- Vortex vigorously for 10 seconds and then shake upside down for 5 minutes
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- To 1 ml of the supernatant add 20 µl Carrez I, vortex
- Add 20 µl Carrez II, vortex
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- Use 50 µl of the supernatant per well in the assay

9.5 Kidney

- Remove fat from the sample and homogenize the kidney with an Ultra-Turrax or a mixer
- Add 2 ml washing buffer (see 10.1.) to 1 g of homogenized sample
- Vortex for 10 seconds and then shake vigorously upside down for 5 minutes
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- To 1 ml of the supernatant add 10 µl Carrez I, vortex
- Add 10 µl Carrez II, vortex
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- Use 50 µl of the supernatant per well in the assay

9.6 Fish

- Remove fat and skin from a representative amount of the sample and homogenize the fish with an Ultra-Turrax or a mixer
- Add 2 ml washing buffer (see 5.1.) to 1 g of homogenized sample
- Vortex for 10 seconds and then shake vigorously upside down for 5 minutes
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- To 1 ml of the supernatant add 5 µl Carrez I, vortex
- Add 5 µl Carrez II, vortex
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- Use 50 µl of the supernatant per well in the assay

9.7 Shrimp/ egg

- Homogenize the sample with an Ultra-Turrax or a mixer
- Add 2 ml methanol to 1 g of homogenized sample
- Vortex vigorously for 30 seconds
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- Transfer 1 ml of the supernatant into a new vial and evaporate to complete dryness under nitrogen or air stream at 60 °C
- Reconstitute the dried residue in 1 ml washing buffer (see 10.1.)
- Use 50 µl of the supernatant per well in the assay

10. Test implementation

10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

As wash buffer PBS Tween buffer is needed; please use the enclosed buffer salt (see section 4. Reagents provided). To prepare the buffer, dissolve the entire contents of the pouch in 1 l distilled water. The dissolved wash buffer can be stored for approximately 4 to 6 weeks at 2 - 8 °C.

Alternative: Dissolve the contents of the pouch in 100 ml distilled water (10-fold concentrate). The solution can be stored for approximately 8 - 12 weeks at room temperature (20 - 25 °C). To prepare the ready-to-use solution, mix 1 part of the 10-fold concentrate with 9 parts of distilled water.

Return all reagents to 2 - 8 °C immediately after use.

10.2 Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples for run in duplicate. Record standard and sample positions.
2. Pipette 50 µl each of the standard and the samples (see section 9. Preparation of Samples) as duplicates into the respective wells. Use new pipette tips for each standard or sample.
3. Pipette 50 µl each of conjugate into the respective wells, mix carefully by shaking the plate per hand, and incubate in the dark for 30 minutes at room temperature (20 - 25 °C).
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times) against absorbent paper to ensure complete removal of liquid from the wells. Wash each well using 250 µl diluted wash buffer (see 10.1 Test preparation). Repeat the washing step another two times.
5. Pipette 100 µl each of substrate/chromogen into the wells. Mix carefully by hand and incubate in the dark for 15 minutes at room temperature (20 - 25 °C).
6. Pipette 100 µl stop solution into each well and mix carefully by hand. Measure the absorbance at 450 nm within 30 minutes after adding the stop solution.

11. Results

A special software, the **RIDASOFT® Win.NET Food & Feed** (Art. No. Z9996FF) is optionally available for evaluation of the RIDASCREEN® enzyme immunoassays. The course of the standard curve is shown in the Quality Assurance Certificate (certificate of analysis) enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance standard (or sample)}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the sulfamethazine concentration [µg/l].

Interpretation of results

In order to obtain the sulfamethazine concentration in µg/l (µg/kg) actually contained in a sample, the concentration read from the standard curve must be further multiplied by the corresponding dilution factor. When working according to the instructions, the dilution factors are as follows:

| | |
|--------------|---|
| Milk | 1 |
| Meat | 3 |
| Honey | 1 |
| Liver | 3 |
| Kidney | 3 |
| Fish | 3 |
| Shrimp | 3 |
| Egg | 3 |

Recommendation

In order to ensure a high analytical performance we recommend:

- To do spike experiments to ensure an accurate and correct test procedure

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

| Version number | Chapter and title |
|----------------|--|
| 2017-02-10 | Release version |
| 2021-01-11 | Brief information, Product offer, 4. Reagents provided, 6. Warnings and precautions for the users, 7. Storage instructions |

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.