## RIDASCREEN® Fumonisin

**REF** R3401

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# RIDASCREEN® Fumonisin

## Brief information

RIDASCREEN® Fumonisin (Art. No. R3401) is a competitive enzyme immunoassay for the quantitative determination of fumonisin in corn and corn products (corn flour, corn semolina, polenta, corn for popcorn).

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: extraction, filtration and dilution

Time requirement: sample preparation (for 10 samples)..... approx. 30 min  
test implementation (incubation time)..... 45 min

Detection limit: .....25 µg/kg (ppb)  
(corresponding to the standard substance)

Recovery rate: Corn .....80 – 120 %  
(spiked samples)

Specificity: The specificity of the RIDASCREEN® Fumonisin test was established by analyzing the cross-reactivity to corresponding mycotoxins in buffer system.

Fumonisin B1 ..... 100 %  
Fumonisin B2.....approx. 40 %  
Fumonisin B3.....approx. 100 %

The specificity of the RIDASCREEN® Fumonisin test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the limit of detection and the recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our good ELISA practice (GEP) – manual in the respective version. Such lists minimum standards concerning the framework conditions when using test kits of R-Biopharm AG and performing ELISA analysis. The manual can be retrieved, printed and downloaded from the website [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

## **Product offer**

RIDASCREEN®FAST Fumonisin (R5602)  
Trilogy®Dried Standard Fumonisin B1, B2 (TS-202-2)  
Trilogy®Liquid Standard Fumonisin B1, B2 (TSL-202-2)  
Trilogy®Liquid Standard Fumonisin B1 (TSL-204-2)  
Trilogy®Liquid Standard Fumonisin B2 (TSL-205-2)

## **1. Intended use**

RIDASCREEN® Fumonisin is a competitive enzyme immunoassay for the quantitative determination of fumonisin in corn and corn products (corn flour, corn semolina, polenta, corn for popcorn).

## **2. General**

Fumonisin are carcinogenic, neuro, hepato and pneumo toxic metabolites of *Fusarium moniliforme*, a mould fungi which grows hostspecific on corn.

The dose of fumonisin for the release of toxic effects differs significantly depending on the animal species. A concentration of approx. 5 - 10 mg/kg fumonisin in feed induces neuro toxic effects in horses. In pigs the ingestion of 4 - 16 mg/kg body weight may result in liver cirrhosis and more than 16 mg/kg body weight may lead to pulmonary edema. Chickens tolerate higher concentrations of fumonisin in feed, up to 75 mg/kg. Cattle seem to be insensitive to high fumonisin concentrations.

## **3. Test principle**

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-fumonisin antibodies.

Fumonisin standards or sample solutions, fumonisin enzyme conjugate and anti-fumonisin antibodies are added. Free fumonisin and fumonisin enzyme conjugate compete for the fumonisin antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-fumonisin antibodies are also bound by the

immobilized capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the fumonisin concentration in the sample.

#### 4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
<b>Microtiter plate K</b>	-	Ready to use		96 wells
<b>Standard 1*</b>	White	Ready to use	0 mg/L	1.3 mL
<b>Standard 2*</b>	White	Ready to use	0,025 mg/L	1.3 mL
<b>Standard 3*</b>	White	Ready to use	0,074 mg/L	1.3 mL
<b>Standard 4*</b>	White	Ready to use	0,222 mg/L	1.3 mL
<b>Standard 5*</b>	White	Ready to use	0,666 mg/L	1.3 mL
<b>Standard 6*</b>	White	Ready to use	2 mg/L	1.3 mL
<b>Conjugate</b>	Red	Ready to use		6 mL
<b>Antibody</b>	Black	Ready to use		6 mL
<b>Substrate/Chromogen</b> Red Chromogen Pro	Brown	Ready to use		10 mL
<b>Stop solution</b>	Yellow	Ready to use		14 mL

\*) The dilution factor 70 for the sample has already been considered. Therefore, the fumonisin concentrations of samples can be read directly from the standard curve.

#### 5. Materials required but not provided

##### 5.1. Equipment

- grinder (mill)
- graduated cylinder: (plastic or glass) 250 ml
- for preparing sample extract: filter funnel and 300 ml flask
- filter paper: Whatman No. 1 or equivalent
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes
- microtiter plate spectrophotometer (450 nm)
- optional: shaker, Ultra-Turrax

## 5.2. Reagents

- methanol
- 70 % methanol solution: prepare 70 % methanol solution by mixing 70 ml methanol (100 %) with 30 ml distilled water
- distilled (or deionized) water

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The standard solutions contain fumonisin, particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

Decontamination of the glassware and toxin-content solutions is best carried out using a sodium hypochlorite (bleach) solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

Ensure the proper and responsible disposal of all reagents and materials after their use. For disposal, please adhere to national regulations.

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

Fumonisin is light sensitive, therefore, avoid exposure to direct light.

The red stained substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen solution prior to test implementation
- a value of less than 0.8 absorbance units ( $A_{450\text{ nm}} < 0.8$ ) for the zero standard

## 9. Preparation of Samples

The samples should be stored in a cool place, protected from light.

Bring all reagents and samples to room temperature (20 - 25 °C / 68 - 77 °F) before use and perform the sample preparation at room temperature.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- weigh 5 g of ground and homogenized sample into a suitable container and add 25 ml of 70 % methanol \*)
  - blend the sample by Ultra-Turrax for 2 minutes or shake vigorously for 3 minutes (manually or with shaker)
  - filter the extract through Whatman No. 1 filter or centrifuge (10 min / 3.500 g / room temperature)
  - dilute the filtered sample extract 1:14 (1+13) with distilled or deionized water (e. g. 100 µl extract + 1.3 ml water)
  - use 50 µl of the diluted filtrate per well in the test
- \*) sample size may be increased if required, but the volume of methanol/water must be adapted accordingly, e.g.: 50 g in 250 ml of 70 % methanol

## 10. Test implementation

### 10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use and perform the test at room temperature.

The fumonisin standards are provided ready to use. The dilution factor 70 for the sample has been considered when labeling. Therefore, the fumonisin concentration of samples can be read directly from the standard curve.

Return all reagents to 2 - 8 °C (35 - 46 °F) immediately after use.

## 10.2. Test procedure

Accurate washing is very important. Carefully follow the recommended washing procedure as outlined in the test procedure. Do not allow microwells to dry up totally and avoid prolonged intervals between the working steps.

Insert a sufficient number of microtiter wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.

1. Add 50 µl of standard or prepared sample to separate duplicate wells. Use a new pipette tip for each standard or sample.
2. Add 50 µl of conjugate to each well.
3. Add 50 µl of antibody to each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl distilled water and pour the liquid out again. Repeat the washing procedure two times.
5. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
6. Add 100 µl of stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of the stop solution.

## 11. Results

Specific software, the RIDASOFT® Win.NET (Art. No. Z9996), is available for evaluation of the RIDASCREEN® enzyme immunoassays.

For sample analysis we recommend to perform double determinations. Please use the RIDASOFT® Win.NET cubic spline function for test evaluation.

**For further information and applications please contact your local distributor of R-Biopharm AG ([sales@r-biopharm.de](mailto:sales@r-biopharm.de)).**

## Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG.