
r-biopharm®

# **RIDASCREEN® FAST Cashew**

**Art. No. R6872**

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# RIDASCREEN®FAST Cashew

## Brief information

The RIDASCREEN®FAST Cashew (Art. No.: R6872) is a sandwich enzyme immunoassay developed for the quantitative analysis of Cashew in food. Exemplarily for the food groups pastries and sweets have been tested during the assay's development: cookies, chocolate and ice cream. The assay can be likely used for analysis of other food samples too; this has to be checked by the user.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization, extraction and centrifugation

Time requirement: sample preparation (for 10 samples) .....approx. 20 min  
test implementation (incubation time) .....30 min

Limit of detection: 0.13 (0.1 – 0.19) mg/kg (ppm) cashew

Limit of quantification: 2.5 mg/kg (ppm) cashew

Specificity: The antibodies specifically detect proteins from cashew.  
There is a cross-reactivity to pistachio.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. corn bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

## **Related products**

Bioavid Cashew Kerne/Cashew Kernel (Art.-No. BL610-10 und 25)

SureFood® ALLERGEN Cashew (Art.-No. S3615)

## 1. Intended use

The RIDASCREEN®FAST Cashew (Art. No.: R6872) is a sandwich enzyme immunoassay developed for the quantitative analysis of Cashew in food. Exemplarily for the food groups pastries and sweets have been tested during the assay's development: cookies, chocolate and ice cream. The assay can be likely used for analysis of other food samples too; this has to be checked by the user. Limit of detection and limit of quantification depend on the tested food matrix, the level of food processing, and the extraction method. Detailed results and further validation information about other sample matrices can be found in the brochure Product Information. Additionally, results from proficiency tests and inter laboratory trials may exist for different food matrices. Further applications are frequently validated in our laboratory and provided in the Application Notes.

## 2. General

Cashew proteins can be present as an ingredient or as a contamination in raw and processed products. According to the **regulation (EU) No. 1169/2011**, cashew must be declared on food labels as it can induce allergic reactions. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand. The cashew tree (*Anacardium occidentale*) produces the cashew apple and cashew kernels. 60 % of cashew is consumed as snacks, remaining 40 % are used in candies. The presence of undeclared cashew is a potential risk for cashew allergic patients. The tree nut allergy is the major allergy worldwide. Beside walnut, cashew is the second most common tree nut allergy. Cashew is a very potent allergen and has to be avoided by allergic patients.

## 3. Test principle

The wells of the microtiter strips are coated with specific antibodies to cashew proteins. By adding standards and samples to the wells, cashew protein present will bind to the specific antibodies. In a washing step components not bound are removed. Then antibody conjugated to peroxidase is added. This antibody conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. The detection of cashew protein takes place by adding Substrate/Chromogen. The enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the cashew concentration of the sample. The result is expressed in mg/kg cashew.

## 4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses).

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
Allergen extraction buffer	Green	<b>Concentrate</b>	<b>10x</b>	100 ml
Standard 1*	Transparent	Ready to use	0.0 mg/kg	2.6 ml
Standard 2*	Transparent	Ready to use	2.5 mg/kg	2.6 ml
Standard 3*	Transparent	Ready to use	5.0 mg/kg	1.3 ml
Standard 4*	Transparent	Ready to use	10.0 mg/kg	1.3 ml
Standard 5*	Transparent	Ready to use	20.0 mg/kg	1.3 ml
Wash buffer	Brown	<b>Concentrate</b>	<b>10x</b>	100 ml
Conjugate	Red	<b>Concentrate</b>	<b>11x</b>	0.7 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 ml
Stop solution	Yellow	Ready to use		14 ml

\*) The dilution factor 20 for the sample has already been considered when labeling. Therefore, the cashew concentrations of samples can directly be read from the standard curve.

## 5. Reagents required but not provided

### 5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal glass vials
- water bath
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenisator
- graduated pipettes
- graduated cylinder
- paper filter
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes
- if possible multichannel pipette and/or multistep pipette

## 5.2. Reagents:

- distilled or deionized water
- fish gelatin (Fisch-Gelatine liquid, Serva, No.: 22156.02)
- skim milk powder (food quality)

## 6. Warnings and precautions

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at room temperature

Reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen prior to test implementation
- a value of less than 0.8 absorbance units ( $A_{450\text{ nm}} < 0.8$ ) for standard 5

## 9. Preparation of Samples

Working devices such as a mill, glass vials or spatulas must be cleaned before and after each sample preparation to remove any remains of allergen and to avoid contamination.

The **Allergen extraction buffer** is provided as a **10fold concentrate**. Before dilution, dissolve possibly formed crystals by heating (water bath 37 °C / 98.6 °F)

and mix well. Dilute the heated concentrate 1:10 (1+9) with distilled water (e.g. 100 ml concentrate + 900 ml dist. Water). The diluted Allergen Extraction buffer is either stable at 20 – 25 °C (68 - 77 °F) for four weeks or at 2 - 8 °C (35 - 46 °F) for 12 weeks.

#### 9.1. Sample preparation solid samples (**without** skim milk powder)

- homogenize a representative amount of sample (5 – 50 g)
- weigh 1 g of a representative sample and add 20 ml preheated (60 °C / 140 °F) and finally diluted Allergen Extraction buffer, close the vial, mix vigorously and incubate the solution for 10 min at 60 °C / 140 °F in a water bath
- let the sample cool down (e.g. in an ice bath) and centrifuge 10 min 2500 g if possible at 4 °C (39 °F) and / or filter the extract (alternatively, 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl extract per well in the assay

#### 9.2. Sample preparation liquid samples (**without** skim milk powder)

- take 1 ml of a representative sample and add 19 ml preheated (60 °C / 140 °F) and finally diluted Allergen Extraction buffer, close the vial, mix vigorously and incubate the solution for 10 min at 60 °C / 140 °F in a water bath
- let the sample cool down (e.g. in an ice bath) and centrifuge 10 min 2500 g if possible at 4 °C (39 °F) and / or filter the extract (alternatively, 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl extract per well in the assay

#### 9.3. Sample preparation chocolate with a cocoa fraction of >75 % and pure cocoa (**with** gelatin allergen extraction buffer)

For the extraction of chocolate with a cocoa fraction of >75 % and pure cocoa a gelatin Allergen extraction buffer is required: e.g. add 1.25 g liquid fish gelatin (e.g. Fish-gelatine liquid, Serva, No.: 22156.02, 45 % solid components) to a graduated cylinder and fill up to 100 ml with final diluted Allergen extraction buffer, mix well.

- homogenize a representative amount of sample (5 – 50 g)
- weigh 1 g of a representative sample and add 20 ml preheated (60 °C / 140 °F) gelatin Allergen extraction buffer, close the vial, mix vigorously and incubate the solution for 10 min at 60 °C / 140 °F in a water bath

- let the sample cool down (e.g. in an ice bath) and centrifuge 10 min 2500 g if possible at 4 °C (39 °F) and / or filter the extract (alternatively, 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl extract per well in the assay

#### 9.4. Sample preparation spices (**with** skim milk Allergen extraction buffer)

For the extraction of spices a skim milk Allergen extraction buffer is required: e.g. solve 5 g skim milk powder (food quality) in 100 ml final diluted allergen extraction buffer.

- homogenize a representative amount of sample (5 – 50 g)
- weigh 1 g of a representative sample and add 20 ml preheated (60 °C / 140 °F) skim milk Allergen extraction buffer, close the vial, mix vigorously and incubate the solution for 10 min at 60 °C / 140 °F in a water bath
- let the sample cool down (e.g. in an ice bath) and centrifuge 10 min 2500 g if possible at 4 °C (39 °F) and / or filter the extract (alternatively, 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl extract per well in the assay

#### **Remark:**

**The gelatin Allergen extraction buffer can be stored at room temperature (20 - 25 °C / 68 - 77 °F) for 1 week, the skim milk allergen extraction buffer can be stored at 2 - 8 °C (35 - 46 °F) for 1 day.**

**Further dilutions should be prepared with the finally diluted Allergen Extraction buffer.**

**The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 3 days.**

## **10. Test implementation**

### 10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **Conjugate** (bottle with red cap) is provided as a 11fold concentrate. Since the diluted conjugate solution has a limited stability, only the amount which actually is needed should be diluted. Before pipetting, the conjugate concentrate

should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) with distilled water (e. g. 200 µl concentrate + 2 ml distilled water, sufficient for 2 microtiter strips).

The **Wash Buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml dist. water). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at room temperature for approx. four weeks.

## 10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at the same time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one, Cat.-No. 655101 or Mikrotiter Assembly Breakable Strip 1x8, Thermo Scientific) should be used to avoid a time shift over the microtiter plate. All standards and samples (at least 150 µl for every well) are pipetted into the uncoated plate and then 100 µl are quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standard and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard or prepared sample to separate duplicate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1). Repeat two more times.
4. Add 100 µl of the conjugate to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 -77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.

7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

## 11. Results

A special software, the RIDA<sup>®</sup>SOFT Win / RIDA<sup>®</sup>SOFT Win.net (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN<sup>®</sup> enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

In comparison with the certificate, higher values of the absorbance ( $A_{450\text{ nm}}$ ) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or allergen contamination.

A further dilution and new detection of the samples is recommended for absorbance values ( $A_{450\text{ nm}}$ ) > standard 5.

### Please note

When working in accordance with the regulation stated, the dilution factor is 20. The cashew concentration can be read directly from the standard curve (see 4. \*) - the sample dilution factor of 20 is already taken into account.

For sample dilutions of more than 1:20, the further dilution factor must be considered for the calculation of the cashew concentration.

### In general

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

The reduction in recovery depends strongly on the temperature and the duration of the heat treatment. If samples are heat treated at high temperature the recovery can be significantly reduced.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross reactivities and exemplary analyzed matrices are described in the Validation report.

## **Recommendation**

In order to ensure a high analytical performance we recommend:

- to analyze each sample material in duplicates
- to use also allergen-free and allergen- containing (spiked) samples as test controls
- in case of extremely acid or alkaline samples, the pH should be adjusted to a neutral pH
- due to the multitude of food types, matrix effects cannot be excluded.
- to carry out spiking experiments for an accurate and correct procedure
- to perform SureFood<sup>®</sup> PCR for confirmation of the result
- For details using the ChemWell<sup>®</sup> or GEMINI automation please contact sales@r-biopharm.de.

For further product information and application notes, please contact sales@r-biopharm.de.

**The product information folder with further information is available on request from your local distributor or R-Biopharm AG.**

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG