# Premi<sup>®</sup>Test Incubator

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## 1. Introduction

Premi<sup>®</sup>Test is a microbial screening test for the detection of antibiotic residues in food products of animal origin.

Premi<sup>®</sup>Test is based on the sensitivity of the *Bacillus stearothermophilus var. calidolactis* to antibiotics. *Bacillus stearothermophilus* is a thermophilic bacterium with optimal growth conditions at 64 °C (147.2 °F). A standardized number of spores are imbedded in agar medium containing selected nutrients and the pH-indicator Bromocresol. If no antibiotic residues are present in a tested sample, the germinated spores will multiply during incubation at 64 °C (147.2 °F) with the production of acid. This will be visible by a color change of the pH-indicator from purple to yellow. When antibiotic residues are present in concentrations above the detection limit, the spores will be unable to germinate and therefore no color change will be observed.

While meat (beef, pork, poultry), fish, liver and serum/plasma samples have to be incubated continuously at 64 °C (147.2 °F), shrimp, kidney, egg and feed samples have to be heated additionally at 80 °C (176 °F) for 10 min before continuous incubation at 64 °C (147.2 °F).

Therefore, the Premi<sup>®</sup>Test incubator comes with 2 different programs, which allow the incubation of all samples:

- 'Program A' allows a continuous incubation of the Premi<sup>®</sup>Test ampoules at 64 °C (147.2 °F) with alarm time.
- 'Program b' allows incubating the Premi<sup>®</sup>Test ampoules automatically for 10 min at 80 °C (176 °F) followed by incubation at 64 °C (147.2 °F).

## 2. Incubator inventory

| Item                 | Quantity |
|----------------------|----------|
| Incubator            | 1        |
| Power supply         | 1        |
| Plug socket adaptors | 3        |
| Manual               | 1        |

## 3. Technical specifications


|                                  |   |
|----------------------------------|---|
| Model                            | MiniS   |
| Dimensions (WxDxH)               | 99 x 135 x 34 mm  |
| Weight                           | 0.4 kg  |
| Power supply                     | DC 12 V 18 W  |
| Power                            | 15 W  |
| Average power consumption        | 5 W   |
| Ambient temperature              | 5 to 35 °C (41 to 95 °F)                                |
| Temperature range                | (ambient temperature +5 °C / 9 °F)<br>to 80 °C / 176 °F |
| Accuracy of the temperature      | $\leq \pm 0.5$ °C / 0.9 °F                              |
| Display accuracy                 | $\pm 0.1$ °C / 0.18 °F                                  |
| Heating time from 20 °C to 80 °C | $\leq 20$ min   |
| Timing range                     | 1 - 999 h / 1 - 999 min / 1 - 999 s                     |

## 4. Safety precautions

- Read the manual carefully before use.
- Use the incubator indoor only.
- Use the incubator at an ambient temperature from 5 - 35 °C (41 - 95 °F) only.
- Do not place the incubator near an open window or close to air-conditioning or heating.
- Do not place the incubator in direct sunlight or strong lamplight.
- Keep the incubator away from strong magnetic fields or corrosive gas.
- Do not use the incubator in a wet or dusty environment.
- Avoid strong mechanical shocks.
- Switch the incubator off when you have finished your analysis and pull out the power supply.
- Store the incubator dry at an ambient temperature from 5 - 35 °C (41 - 95 °F). Avoid freezing.
- Do not open or repair the incubator. Contact your local R-Biopharm distributor for assistance.
- The incubator should not be used in a room where antibiotics are handled.
- **Do not fill the ampoule holes of the incubator with water e.g. for temperature control measurements. Danger to life due to electric shock!**
- In case of non-specified use, your safety is not guaranteed.

# 5. Schematic overview

## 6. Installing the incubator

- Place the incubator on a level surface.
- Insert the power connector of the power supply to the power input port of the incubator. Insert the connector of the power supply itself to the plug socket (100 - 240 V).
- The main switch is on the rear of the device. Push 'I' to switch on the incubator, and push 'O' to switch off the incubator.
- Switch on the incubator.
- The incubator will start up in the last used program and heats up to the last used temperature.
- During heating the 'Heating' indicator flashes once per second. After reaching the set temperature the 'Heating' indicator lights constantly.
- For best results it is recommended to wait further 5 min to allow complete temperature stabilization at the set temperature. In case of multiple test runs during one day, leave incubator switched on and don't put cold ampoules in the incubator during an on-going incubation.

## 7. Using the incubator

### 7.1 Check the current program

- Press 'Prog' and '▲' at the same time
- The display will show either 'Program A' or 'Program b'

### 7.2 Choosing the incubation program

- Press and hold 'Prog' and '▼' at the same time for 10 s
- Press '▲' or '▼' to switch between the two programs
- Press 'Start/Stop' to select a program


### 7.3 Program A

In 'Program A' incubation temperature and time can be set. After the set time is reached, the incubator sounds an alert signal and the temperature will stay at the set value.

'Program A' is intended to be used for the incubation of meat-, fish-, liver- and serum/plasma samples. Therefore, incubation temperature is pre-set to 64 °C (147.2 °F) and incubation time to 2 h 40 min.

#### Setting parameters of Program A

- Press and hold 'Prog' for 3 s
- The active parameter will flicker. To switch between parameters (temperature, time and time unit) press 'Prog'.
- To change parameters press '▲' or '▼'.
- If no button is pressed for 5 s, the incubator confirms the actual setting automatically and exits the setting mode.

## Starting 'Program A'

'Program A' cannot be started before the incubator has reached the set temperature (approx. 10 - 15 min to reach 64 °C (147.2 °F)).

- Press 'Start/Stop' to start of 'Program A'. The 'Running' indicator lights constantly.
- Insert the Premi<sup>®</sup>Test ampoules into the incubator and ensure that they are fully inserted, touching wall and bottom of the incubator wells.
- At the end of the incubation time the incubator sounds an alert signal and the 'Running' indicator turns off.
- To stop 'Program A' during incubation, press and hold 'Start/Stop' for 3 s.


## 7.4 Program b

In 'Program b' temperature and time can be set for two consecutive segments. After the incubation time of the first segment is reached, the second segment is started automatically. After the set time of the second segment is reached the incubator sounds an alert signal and the temperature will stay at the set value.

'Program b' is intended for the incubation of shrimp, kidney, egg and feed samples. Therefore, the first segment has a pre-set incubation temperature of 80 °C (176 °F) and incubation time of 10 min. The second segment has a pre-set incubation temperature of 64 °C (147.2 °F) and incubation time of 2 h 30 min.

### Setting parameters of 'Program b':

- Press and hold 'Prog' for 3 s
- The active parameter will flicker. To switch between parameter of the first segment (temperature, time and time unit) press 'Prog'. Press 'Prog' again to enter the parameters of the second segment.
- To change parameters in both segments press '▲' or '▼'.
- If no button is pressed for 5 s, the incubator confirms the actual setting automatically and exits the setting mode.

## Starting 'Program b'

'Program b' cannot be started before the incubator has reached the set temperature (approx. 15 - 20 min to reach 80 °C (176 °F)).

- Press 'Start/Stop' to start 'Program b'. The 'Running' indicator lights constantly. The incubator will automatically incubate the samples for 10 min at 80 °C (176 °F) followed by continuous incubation at 64 °C (147.2 °F) or other set temperature and time.
- Insert the Premi<sup>®</sup>Test ampoules into the incubator and ensure that they are fully inserted, touching wall and bottom of the incubator wells.
- At the end of the incubation time the incubator sounds an alert signal and the 'Running' indicator turns off.
- To stop 'Program b' during incubation, press and hold 'Start/Stop' for 3 s.

## 8. Cleaning

- Before cleaning, switch the incubator off, pull power plug and let cool down to room temperature.
- If the instrument gets dirty, clean it with a dry cloth.
- The wells in the block should be cleaned with a dry cloth or tissue. Do not use water or any other liquids.

## 9. Error messages

| <b>Error</b>             | <b>Cause</b>   | <b>Solution</b>  |
|--------------------------|--|--|
| No display               | No power connection<br>Power failure<br>Controller failure | Connect with power supply.<br>Check main fuse.<br>Contact your local R-Biopharm distributor. |
| Display “Er1” with alarm | Sensor broken circuit                                      | Contact your local R-Biopharm distributor.   |
| Display “Er2” with alarm | Sensor short circuit                                       | Contact your local R-Biopharm distributor.   |
| Display “Er4” with alarm | Over-heating temperature alarm                             | Contact your local R-Biopharm distributor.   |
| No heating of the block  | Heater failure   | Contact your local R-Biopharm distributor.   |
| Press invalid            | Keyboard failure   | Contact your local R-Biopharm distributor.   |

## **R-Biopharm AG**

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64297 Darmstadt, Germany

Telefon: +49 (0) 61 51 - 81 02-0

Fax: +49 (0) 61 51 - 81 02-40

E-Mail: [info@r-biopharm.de](mailto:info@r-biopharm.de)

[www.r-biopharm.com](http://www.r-biopharm.com)