RIDASCREEN[®] FAST Zearalenon

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RIDASCREEN[®]FAST Zearalenon

Brief information

RIDASCREEN[®]FAST Zearalenon (Art. No.: R5502) is a competitive enzyme immunoassay for the quantitative analysis of zearalenone in cereals and feed.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Technicians need no specialized training to perform the RIDASCREEN[®]FAST Zearalenon test, however, free support is offered by the distributor on request, if necessary.

The test was validated with different matrices of cereals and composed feed.

Sample preparation: extraction, filtration and dilution

Time requirement: sample preparation (for 10 samples)
cereals and feed.....approx. 10 min
test implementation (incubation time)15 min

Limit of detection: 17 - 41 µg/kg (ppb)

Limit of quantification: 50 µg/kg (ppb)

1. Intended use

RIDASCREEN®FAST Zearalenon is a competitive enzyme immunoassay for the quantitative analysis of zearalenone in cereals and feed.

2. General

The mycotoxin zearalenone is formed by fungi of the genus *Fusarium*. Zearalenone is a phytohormone which displays, apart from its anabolic properties, mainly estrogenic effects. Because of its estrogenic properties, zearalenone may induce fertility disorders in animals with clinical signs of hyperestrogenism - an aspect of a disease which although reported mainly in hogs, is described in other species such as cows, horses and sheep.

The potential health risk for man induced by this mycotoxin, which can be ingested with foods of vegetable or animal origin, has been discussed intensively.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-zearalenone antibodies.

Zearalenone standards or sample solutions, zearalenone enzyme conjugate and anti-zearalenone antibodies are added. Free zearalenone and zearalenone enzyme conjugate compete for the zearalenone antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-zearalenone antibodies are also bound by the immobilized capture antibodies. Any unbound zearalenone enzyme conjugate is then removed in a washing step. Substrate/Chromogen is added to the wells and bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the zearalenone concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for up to 43 analyses (plus 5 standard analyses). Each test kit contains:

- 1 x Microtiter plate with 48 wells (6 strips with 8 removable wells each)
coated with capture antibodies
- 5 x Zearalenone standard solutions ^{*)}, 1.3 ml each
0 ppb (zero standard), 50 ppb, 100 ppb, 200 ppb, 400 ppb
zearalenone in methanol/water, ready to use
- 1 x Conjugate (3 ml)red cap
peroxidase conjugated zearalenone, ready to use
- 1 x Anti-zearalenon antibody (3 ml).....black cap
monoclonal, ready to use
- 1 x Substrate/Chromogen (10 ml) brown cap
stained red
- 1 x Stop reagent (14 ml) yellow cap
contains 1 N sulfuric acid

*) The dilution factor 10 for the sample has been already considered. Therefore, the zearalenone concentrations of samples can be read directly from the standard curve.

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- graduated cylinder: 100 ml plastic or glass
- glassware for preparing sample extract: filter funnel and 50 ml flask
- grinder (mill)
- optional: shaker
- filter paper: Whatman No. 1 or equivalent
- 50 µl, 100 µl and 1000 µl micropipettes

5.2. Reagents:

- methanol
- 70 % methanol solution: prepare 70 % methanol solution by mixing 70 ml methanol (100 %) with 30 ml distilled or deionized water
- distilled or deionized water

6. Warnings and precautions for the users

The standards contain zearalenone, so avoid contact of the reagent with the skin (use gloves).

Decontamination of the glassware and zearalenone solutions is best carried out using a sodium hypochlorite solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

The stop solution contains 1 N sulfuric acid (R36/38, S2-26).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen solution is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

The test kit can be regularly used at least up to the expiry date (indicated on the kit package), if stored correctly.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the reddish substrate/chromogen solution prior to test implementation
- a value of less than 0.6 absorbance units ($A_{450\text{ nm}} < 0.6$) for the zero standard

9. Preparation of Samples

The samples should be stored in a cool place, protected from light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- weigh 5 g of ground sample and add it to a suitable container with 25 ml of methanol (70 %) *)
- shake vigorously for 3 min (manually or with shaker)
- filter the extract through Whatman No. 1 filter
- dilute 1 ml of the obtained filtrate with 1 ml of distilled or deionized water
- use 50 µl of the filtrate per well in the test

*) Sample size may be increased if required, but the volume of methanol/water must be adapted accordingly.

10. Test implementation

10.1. Preliminary comments

1. Bring all reagents to room temperature (20 - 25 °C, 68 - 77 °F) before use.
2. The specific reaction starts with the addition of the specific antibody. Nevertheless, not more than three strips should be applied to the test, when a single step pipette is utilized. More strips (up to 6) can be applied, when a multistep pipette is used.
3. Return all reagents to 2 - 8 °C (35 - 46 °F) immediately after use.

The **zearalenone standards** are provided ready to use. The dilution factor 10 for the sample has been considered when labeling. Therefore, the zearalenone concentration of samples can be read directly from the standard curve.

10.2. Test procedure

Accurate washing is very important. Do not allow microwells to dry up totally and avoid prolonged intervals between the working steps. Reproducibility in any EIA is largely dependent upon the consistency with which the microwells are washed. Carefully follow the recommended washing sequence as outlined in the EIA test procedure.

Avoid direct sunlight during all incubations. Covering the microtiter plates is recommended.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Pipet 50 µl of standard or prepared sample into separate wells; use a new pipette tip for each standard or sample.
3. Add 50 µl of enzyme conjugate (red cap) to each well.
4. Add 50 µl of the anti-zearalenone antibody solution (black cap) to each well. Mix gently by shaking the plate manually and incubate for 10 min (+/- 1) at room temperature (20 - 25 °C; 68 - 77 °F).
5. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Using a wash bottle or multichannel pipette and fill the wells with distilled or deionized water (250 µl per well). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.
6. Add 100 µl of substrate/chromogen (brown cap) to each well. Mix gently by shaking the plate manually and incubate for 5 min (+/- 0.5) at room temperature (20 - 25 °C; 68 - 77 °F) in the dark.
7. Add 100 µl of stop solution (yellow cap) to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win (Art. No. Z9999), is available to evaluate the RIDASCREEN[®] enzyme immunoassays.

For single determinations we recommend logit/log evaluation and for double or multiple determinations cubic spline should be used.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the zearalenone concentration [$\mu\text{g}/\text{kg}$]. The zearalenone concentration in $\mu\text{g}/\text{kg}$ corresponding to the extinction of each sample can be read from the calibration curve.

12. Sensitivity

12.1. Limit of detection

The limit of detection (LOD) was investigated with cereal samples and mixed feed by repeated measurement of the zero-matrix. LOD was defined as the concentration that corresponds to the mean absorbance of the measurement plus 3 times the standard deviation, and was determined by interpolation from the standard curve. All results were found between 17 and 41 $\mu\text{g}/\text{kg}$ (ppb).

12.2. Limit of quantification

The limit of quantification (LOQ) was investigated by repeated spiking experiments at the lower end of the standard curve; LOQ was found at 50 $\mu\text{g}/\text{kg}$ (ppb), which concentration is well above LOD, corresponding recoveries were found between 64 and 97 %, with CV-values below 20 %.

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG