---

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# Simplex<sup>®</sup> Easy Wine

## 1. Intended use

Fast and optimal DNA-extraction from microorganisms in wine samples and grape must

## 2. Content

The kit contains sufficient reagents for 100 or 10 extractions:

1 x Easy-Reagent (10 mL – 100 rxn; 1 mL – 10 rxn)

1 x Wine Solution

1 x Washing Solution

## 3. Materials required but not provided

### 3.1 Instruments

- Heating block or water bath 37 °C – 95 °C (98.6 - 203 °F)
- Microcentrifuge for 1.5 – 2.0 mL reaction tubes
- Pipettes
- “Vortex”

### 3.2 Reagents and plastic ware

- safe-lock reaction tubes 1.5 – 2.0 mL
- suitable filter tips
- single use gloves
- Lyticase (Sigma L2524)

#### 4. Warnings

Careful use with personal protection according to good laboratory practice is recommended. Do not incorporate. Avoid skin and eye contact with all solutions.

#### 5. Storage

Store the kit at room temperature.

#### 6. Indications of deterioration of reagents

In case of accurate handling deterioration unknown.

#### 7. Preliminary preparations

Prepare Lyticase stock solution:

5U /  $\mu$ L in lyticase buffer (50 % 1xTE + 50 % Glycerol, pH 7.5 - 8.0).

Store solution at  $-20\text{ }^{\circ}\text{C}$  ( $-4\text{ }^{\circ}\text{F}$ )

Prepare Easy-Reagent Wine solution:

**Before first use** add the total content of the Wine Solution (violet cap) to the Easy-Reagent and mix well. Please paste up the enclosed sticker "**Easy-Reagent Wine**". The resulting buffer is stable for 4 months at room temperature.


If the Washing Solution is turbid or frozen incubate the solution some minutes at  $45\text{ }^{\circ}\text{C}$  ( $113\text{ }^{\circ}\text{F}$ ) till it is clear.

## 8. DNA Extraction


### 8.1 DNA-Extraction from yeast

#### Protocol for 1 mL wine or grape must

1. mix the sample and transfer 1.0 mL to a 1.5 mL safe-lock reaction tube
2. place the tubes always in the same manner into the centrifuge (lid of the tube), so that even when no pellet will be visible, it will be clear where the cells/pellet will be

3. centrifuge 5 minutes for at least 6000 x g
4. discard carefully and completely the supernatant. Aspirate from the opposite side of the pellet


5. resuspend the pellet completely in **100 µL Washing Solution** (vortex or pipetting)
6. centrifuge 5 minutes at 6000 x g
7. discard carefully and completely the supernatant (**very important**)  
**Optional:** if the supernatant shows after the washing step strong red flavour, repeat the washing step once or twice (5.-7.)
8. add **100 µL Easy-Reagent Wine** (including Wine solution) to the pellet
9. resuspend the pellet completely (vortex or pipetting)
10. add **3 µL (15 U) Lyticase** and mix carefully
11. incubate 15 min at 37 °C +/- 1°C (98,6 °F +/- 1,8 °F)

12. incubate 10 min at 95 °C +/- 1°C (203 °F +/- 1,8 °F )
13. centrifuge 5 min at 12000 x g
14. transfer the supernatant (contains DNA) to a new 1.5 mL safe-lock tube  
Attention: transfer no pellet
15. use 2.5 µL DNA-suspension for PCR


## 8.2 DNA-Extraction from bacteria and yeast

### Protocol for 1 mL wine or grape must

1. mix the sample and transfer 1.0 mL sample to a 1.5 mL safe-lock reaction tube
2. place the tube always in the same manner into the centrifuge (lid of the tube), so that even when no pellet will be visible, it will be clear where the cells/pellet will be

3. centrifuge 5 minutes at 12000 x g
4. discard carefully and completely the supernatant. Aspirate from the opposite side of the pellet.

5. **resuspend** the pellet completely in **100 µL Washing Solution** (vortex or pipetting)
6. centrifuge 5 minutes at 12000 x g
7. discard carefully and completely the Washing Solution (**very important**)
8. add **100 µL Easy-Reagent Wine** (including Wine Solution)
9. resuspend completely (vortex) the pellet

10. add **3 µL (15 U) Lyticase** and mix carefully
11. incubate 15 min at 37 °C +/- 1°C (98,6 °F +/- 1,8 °F)
12. incubate 10 min at 95 °C +/- 1°C (203 °F +/- 1,8 °F )
13. centrifuge 5 min at 12000 x g
14. transfer the supernatant (contains DNA) to a new 1.5 mL safe-lock reaction tube  
**Attention:** transfer no pellet
15. use 2.5 µL DNA-suspension for PCR

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.