RIDASCREEN[®] Gliadin

Art. Nr. R7001

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RIDASCREEN[®] Gliadin

Brief information

RIDASCREEN[®] Gliadin (R7001) is a sandwich enzyme immunoassay for the quantitative analysis of prolamins from wheat (gliadin), rye (secalin) and barley (hordein) in as gluten-free declared food.

The R5 ELISA RIDASCREEN[®] Gliadin is
–accepted as AOAC-OMA (2012.01) and AACCI 38.50.01
–certified at AOAC-RI (120601)
–Codex Alimentarius Method (Type I)

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit. The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization and extraction
Standard material:	The RIDASCREEN [®] standard material is calibrated to the standard of the Prolamin Working Group.
Time requirement:	sample preparation (for 10 samples) approx. 2 h test implementation (incubation time) 1.5 h
Limit of detection:	0.5 mg/kg (ppm) gliadin or to 1 mg/kg (ppm) gluten (depending on matrix)
Limit of quantification:	2.5 mg/kg (ppm) gliadin or to 5 mg/kg (ppm) gluten
Specificity:	The monoclonal antibody R5 reacts with the gliadin-fractions from wheat and corresponding prolamins from rye and barley.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

RIDASCREEN®FAST Gliadin (R7002)

RIDASCREEN® Gliadin competitive (R7021)

Cocktail (patented) (R7006/R7016)

RIDA® Extraction Solution (colorless) (Art. No. R7098)

Set of 3 processed Gliadin Assay Controls (Art. No. R7012)

RIDA®QUICK Gliadin (R7003/R7004/R7005)

SureFood® Allergen QUANT Gluten (S3206)

SureFood® Allergen Gluten (S3106)

1. Intended use

RIDASCREEN® Gliadin is a sandwich enzyme immunoassay for the quantitative analysis of contaminations by prolamins from wheat (gliadin), rye (secalin), and barley (hordein) in raw products like flours (buckwheat, rice, corn, oats, teff) and spices as well as in processed food like noodles, ready-to-serve meals, bakery products, sausages, beverages and ice cream.

All samples should be extracted with the Cocktail (patented) (R7006/R7016, official R5-Mendez method).

2. General

The use of wheat flour and gluten in foodstuffs is extremely common because of their heat stability and useful effects on e.g. texture, moisture retention and flavour. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley.

Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

According to the Codex Alimentarius (CODEX STAN 118/1979) two categories for labeling of food according to the gluten content now exist:

1.) Food products which contain less than 20 mg/kg can be labeled as "**gluten-free**".

2.) Food products labeled as "**very low gluten**" can have a gluten content above 20 and up to 100 mg/kg.

The threshold of 20 mg/kg has been adopted by many national legislations in many countries.

3. Test principle

The basis of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific R5 antibodies against gliadins. By adding the standard or sample solution to the wells, present gliadin will bind to the specific capture antibodies. The result is an antibody-antigen-complex. Components not bound by the antibodies are then removed in a washing step. Then R5 antibody conjugated to peroxidase is added. This conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. Enzyme substrate and chromogen are added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the gliadin concentration of the sample.

4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 wells
Buffer	White	Concentrate	5x	60 ml
Standard 1	Transparent	Ready to use	0 ng / ml gliadin	1.3 ml
Standard 2	Transparent	Ready to use	5 ng / ml gliadin	1.3 ml
Standard 3	Transparent	Ready to use	10 ng / ml gliadin	1.3 ml
Standard 4	Transparent	Ready to use	20 ng / ml gliadin	1.3 ml
Standard 5	Transparent	Ready to use	40 ng / ml gliadin	1.3 ml
Standard 6	Transparent	Ready to use	80 ng / ml gliadin	1.3 ml
Wash buffer	Brown	Concentrate	10x	100 ml
Conjugate	Red	Concentrate	11x	1.2 ml
Substrate	Green	Ready to use		7 ml
Chromogen	Blue	Ready to use		7 ml
Stop Solution	Yellow	Ready to use		14 ml

5. Reagents required but not provided

5.1. Equipment

- microtiter plate spectrophotometer (450 nm)
- centrifuge, centrifugal vials (e.g. Brand 10742512)
- shaker
- laboratory mincer / grinder, pestle and mortar, ultra-turrax or homogenizer
- water bath (50 °C / 122 °F)
- graduated pipettes
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes

5.2. Reagents

- distilled or deionized water
- gluten-free **skim milk powder** (food quality)
- **Cocktail (patented)** (R7006 / R7016, 105 ml / 1000ml) and **ethanol solution (80 %)**: i.e. add 120 ml ethanol p.a. to 30 ml distilled water and shake well

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The colorless chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the chromogen solution prior to test implementation
- a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 6

9. Preparation of Samples

9.1. Preliminary comments

Airborne cereal dust and dirty laboratory equipment lead to gliadin contamination of the assay. Therefore, before starting and during the assay wear gloves.

- clean surfaces, glass vials, mincers and other equipment with 40 % ethanol or 2-propanol (see 5.2.)
- carry out the sample preparation in a room isolated from the ELISA procedure
- check for gliadin contamination of reagents and equipment with the test strips RIDA[®]QUICK Gliadin (R7003/R7004/R7005)
- it is recommended to work **under a chemical hood**, because of β -mercaptoethanol in the Cocktail (patented)
- β -mercaptoethanol can disturb the ELISA, therefore dilute the samples **at least 1:500** (recommendation 1:500 for samples with approx. 20 mg/kg gluten and 1:2500 for samples with approx. 100 mg/kg gluten).

9.2. Extraction with the Cocktail (patented) (R7006 / R7016, official AOAC method)

Homogenize well a sufficient amount (at least 5 g or 5 ml) of sample (grind it thoroughly to powder and mix well or mix well the solution respectively).

- liquid food samples:** use 0.25 ml of the homogenized sample and add 2.5 ml of the Cocktail (patented), close the vial and mix well
- other food samples (e.g. soy and quinoa containing samples):** weigh 0.25 g of the homogenized sample and add 2.5 ml of the Cocktail (patented), close the vial and mix well
- tannin and polyphenol containing food samples (e.g. chocolate, coffee, cocoa, chestnut flour, buckwheat, millet and spices):** weigh 0.25 g of the homogenized sample, add 0.25 g of skimmed milk powder and add 2.5 ml of the Cocktail (patented), close the vial and mix well
- meat and sausages:** in these matrices gliadin may be not distributed evenly, therefore, weigh 50 g sample and homogenize: weigh 0.25 g of the homogenized sample and add 2.5 ml of the Cocktail (patented), close the vial and mix well

–**oat samples:** gliadin may not be distributed evenly, furthermore the samples are difficult to homogenize. Therefore, homogenize 200g, then carry out the extraction with at least the fourfold amount of reagents: weigh 1 g of the homogenized sample and add 10 ml of the Cocktail (patented), close the vial and mix well

Please further extract all samples as described in the following:

- incubate for 40 min at 50 °C (122 °F)
- let the sample cool down and then mix it with 7.5 ml 80 % ethanol (see 5.2.) (for oat samples: 30 ml 80 % ethanol)
- close the vial and shake for 1 h up side down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F)
- centrifuge: 10 min, at least 2500 g, at room temperature (20 - 25 °C / 68 - 77 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged with high speed for 10 min in reaction caps by using a microcentrifuge)
- transfer the supernatant in a screw top vial
- dilute the sample 1:12.5 (1+11.5 / 80 µl + 920 µl) with diluted sample diluent (see 10.1.): the final dilution factor is 500
- use **immediately** 100 µl per well in the assay

Remark:

The supernatant obtained after the centrifugation or the filtrate can be stored in a tightly closed vial in the dark at room temperature (20 - 25 °C / 68 - 77 °F) up to eight weeks.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **buffer** is provided as a concentrate (5fold). Only the amount which actually is needed should be diluted 1:5 (1+4) with distilled water (e. g. 3 ml concentrate + 12 ml distilled water, sufficient for the dilution of 10 samples). Make sure that the buffer is not contaminated with gliadin.

The **conjugate** (bottle with red cap) is provided as a concentrate (11fold). Since the diluted conjugate has a limited stability, only the amount which actually is needed should be diluted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11

(1+10) with distilled water (e. g. 200 µl concentrate + 2 ml distilled water, sufficient for 2 microtiter strips). Take care that the water is not contaminated with gliadin.

The **washing buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml dist. water). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 -77 °F) for four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than 6 strips (48 wells) at a time. In the case of more than six strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard solution or sample to separate duplicate wells and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
4. Add 100 µl of the diluted conjugate (see 10.1.) to each well and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 50 µl of substrate and 50 µl of chromogen to each well. Mix gently by skaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 min after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win / RIDA[®]SOFT Win.net (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays. The calculation should be done by use of a cubic spline function.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit. For quality assurance Gliadin Assay Controls (R7012) should be used.

In comparison with the certificate, higher values of the absorbance ($A_{450\text{ nm}}$) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or gliadin contamination.

A further dilution and new detection of the samples is recommended for absorbance values ($A_{450\text{ nm}}$) > standard 6.

The gliadin concentration in ng/ml (ppb) is read from the RIDA[®]SOFT Win calibration curve and must be further multiplied by the dilution factor of at least 500. This result is then multiplied by 2 in order to obtain the gluten concentration (gliadin represents 50 % of the proteins present in gluten, Codex Definition). The RIDA[®]SOFT Win (Version 1.93 or newer) indicates the results in gliadin and gluten.

Example:

The absorbance value of a sample corresponds to 10 ng/ml gliadin in the calibration curve. Multiplying by the recommended dilution factor 500 leads to 5000 ng/ml, corresponding to 5 mg/kg (ppm) gliadin, respectively 0.0005 % gliadin. To calculate the gluten content, it is necessary to multiply by factor 2 which results in 10 mg/kg gluten, respectively 0.001 % gluten. This sample is considered to be gluten-free, because the gluten concentration is below 20 mg/kg.

In general:

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross reactivities and exemplary analyzed matrices are described in the Validation report.

Recommendation:

In order to ensure a high analytical performance:

- Each sample material should be analyzed in duplicates
- Use also gluten free and gluten containing (spiked) samples as test controls
- Due to the multitude of food types, matrix effects cannot be excluded. To ensure an accurate result spike experiments are recommended
- Confirm results with PCR (e.g. SureFood® Allergen QUANT Gluten, Art.No. S3201)
- In the production of foods such as beer or sourdough, proteins are fragmented. In sandwich ELISAs protein fragments lead to a reduced recovery, such samples should be analyzed with a competitive ELISA test systems like the RIDASCREEN® Gliadin competitive (R7021).
- For details using the ChemWell® or GEMINI automation please contact sales@r-biopharm.de.

Further application notes:

- Sample preparation for processed food with the RIDA® Extraction Solution (colorless) (Art. No. R7098) - **only after validation**
- Sample preparation for raw materials with ethanol.
- Sample preparation for polyphenol containing raw materials (e.g. chocolate, coffee, cacao, buckwheat) with fish gelatine and ethanol

The product information folder with further information is available on request from your local distributor or R-Biopharm AG.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG