# **RIDASCREEN<sup>®</sup> 17 $\beta$ -Östradiol**

**Art. No. R2301**

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# RIDASCREEN® 17β-Östradiol

## Brief information

RIDASCREEN® 17β-Östradiol (Art. No.: R 2301) is a competitive enzyme immunoassay for the quantitative analysis of 17β-estradiol in bovine plasma.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit. The test kit is sufficient for 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: extraction, freezing and evaporation

Time requirement: sample preparation (for 10 samples) ..... approx. 2 h  
test implementation (incubation time) ..... 2 h 30 min

Detection limit: approx. 20 ng/l (ppt)  
(corresponding to the standard substance)

Recovery rate: in spiked bovine plasma ..... approx. 85 % ± 12 %  
(corresponding to the standard substance)

The specificity of the RIDASCREEN® 17β-Östradiol test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Specificity:	17β-Estradiol.....	100 %
	Estradiol-3-benzoate .....	50 %
	17α-Estradiol.....	0.9 %
	Estron.....	0.7 %
	Trenbolone.....	1.0 %
	19-Nortestosterone.....	0.5 %
	Testosterone .....	< 0.25 %
	Zeranol.....	< 0.25 %
	DES.....	< 0.25 %
	Estriol.....	< 0.1 %

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

## Related products

RIDA® 17β-Östradiol Spiking Solution .....	(R2399)
RIDASCREEN® Testosteron .....	(R2401)
RIDA® Testosteron Spiking Solution .....	(R2499)

### 1. Intended use

RIDASCREEN® 17β-Östradiol is a competitive enzyme immunoassay for the quantitative analysis of 17β-estradiol in bovine plasma.

### 2. General

17β-estradiol is a natural sexual hormone. The levels of 17β-estradiol in bovine plasma from peripheral blood are very low. In calves not more than 1 - 2 pg/ml are present and in cycling females the concentration of 17β-estradiol is 1 - 5 pg/ml. Only in plasma of pregnant animals or in plasma of animals illegally treated with anabolic agents, levels in the range of 100 to 1000 pg/ml can be found.

### 3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-17β-estradiol antibodies. Standards or sample, 17β-estradiol conjugate and anti-17β-estradiol antibodies are added. Free and enzyme conjugated 17β-estradiol compete for the antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-17β-estradiol antibodies are also bound by the immobilized capture antibodies. Any unbound conjugate is then removed in a washing step. Substrate and chromogen are added to the wells and incubated. Bound conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the 17β-estradiol concentration in the sample.

## 4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses).

Komponente	Deckelfarbe	Zustand		Inhalt
Microtiter plate	-	Ready to use		96 wells
Standard 1	White	Ready to use	0 ng/ l	1.3 ml
Standard 2	White	Ready to use	50 ng/ l	1.3 ml
Standard 3	White	Ready to use	200 ng/ l	1.3 ml
Standard 4	White	Ready to use	800 ng/ l	1.3 ml
Standard 5	White	Ready to use	3200 ng/ l	1.3 ml
Standard 6	White	Ready to use	12800 ng/ l	1.3 ml
Conjugate	Red	Concentrate	11x	0.7 ml
Antibody	Black	Concentrate	11x	0.7 ml
Conjugate/antibody/sample buffer	White	Ready to use		50 ml
Substrate	Green	Ready to use		7 ml
Chromogen	Blue	Ready to use		7 ml
Stop solution	Yellow	Ready to use		14 ml

## 5. Materials required but not provided

### 5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- freezer -25 °C (-77 °F) or -60 °C (-140 °F)
- water bath 60 °C (140 °F)
- graduated pipettes
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes

### 5.2. Reagents:

- tert. butylmethylether
- petroleum ether

## **6. Warnings and precautions for the users**

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## **7. Storage instructions**

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The colorless chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

## **8. Indication of instability or deterioration of reagents**

- any coloration of the colorless chromogen prior to test implementation
- a value of less than 0.6 absorbance units ( $A_{450\text{ nm}} < 0.6$ ) for standard 1

## 9. Preparation of Samples

### 9.1. Bovine plasma

- extract 1 ml bovine plasma with 5 ml ether mixture (tert. butylmethylether / petroleum ether (30:70, v/v) in a glass tube by shaking thoroughly at room temperature 20 - 25 °C (68 - 77 °F) for 20 min
- freeze the extracted solution at -25 °C (-77 °F) for 60 min or at -60 °C (-140 °F) for 30 min
- decant the ether supernatant into a glass tube and evaporate it at 60 °C (140 °F), e. g. in a water bath
- redissolve the dried residue in 400 µl sample buffer and mix thoroughly
- employ 20 µl per well in the assay

#### Remark:

Plasma samples, as well as sample extracts obtained, can be stored at 2 - 8 °C (36 - 46 °F) for up to one week.

## 10. Test implementation

### 10.1. Preliminary comments

Bring all reagents to room temperature 20 - 25 °C (68 - 77 °F) before use.

The **17β-estradiol conjugate** (bottle with red cap) is provided as a concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is needed should be reconstituted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in conjugate buffer (e. g. 200 µl concentrate + 2 ml conjugate buffer, ready to use sufficient for 4 microtiter strips).

The **anti-17β-estradiol antibody solution** (bottle with black cap) is provided as a concentrate. Since the diluted antibody has a limited stability, only the amount which actually is needed should be reconstituted. Before pipetting, the antibody concentrate should be shaken thoroughly. For reconstitution, the antibody concentrate is diluted 1:11 (1+10) in antibody buffer (e. g. 200 µl concentrate + 2 ml antibody buffer, ready to use sufficient for 4 microtiter strips).

## 10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 20 µl of each standard or prepared sample to separate duplicate wells and add 50 µl of diluted conjugate to each well.
3. Add 50 µl of the diluted anti-17β-estradiol antibody to each well. Mix gently by shaking the plate manually and incubate 2 h at room temperature (20 - 25 °C / 68 - 77 °F).
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all wells with 250 µl of distilled water and pour out the liquid again. Repeat two more times.
5. Add 50 µl of substrate and 50 µl of chromogen to each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
6. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 minutes after addition of stop solution.

## 11. Results

A special software, the RIDA<sup>®</sup>SOFT Win / RIDA<sup>®</sup>SOFT Win.net (Art. No. Z9996), is available for evaluation of the RIDASCREEN<sup>®</sup> enzyme immunoassays.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = B/B_0 (\%)$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the 17 $\beta$ -estradiol concentration [ng/l].

In order to obtain the 17 $\beta$ -estradiol concentration in ng/l, which is actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factor is as follows:

bovine plasma .....0.4

**For further information or applications please contact your local distributor or R-Biopharm AG.**

**Further application notes:**

–Sample preparation for meat

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG