# **RIDASCREEN<sup>®</sup> Bacitracin**

**Art. No. R2901**

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# RIDASCREEN<sup>®</sup> Bacitracin

## Brief information

RIDASCREEN<sup>®</sup> Bacitracin (Art. No.: R2901) is a competitive enzyme immunoassay for the quantitative analysis of bacitracin in milk, meat, eggs, feed and urine.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:      milk: pH adjustment, defatting, dilution  
                                 meat, eggs and feed: homogenization, extraction,  
                                 centrifugation, dilution  
                                 urine: dilution

Time requirement:      sample preparation (for 10 samples)  
                                 milk .....approx. 30 min  
                                 meat.....approx. 45 min  
                                 eggs .....approx. 45 min  
                                 feed.....approx. 45 min  
                                 urine .....approx. 15 min  
                                 test implementation (incubation time) ..... 1.5 h

Detection limit:      milk ..... approx. 11 µg/kg (ppb)  
(corresponding to the      meat..... approx. 9 µg/kg (ppb)  
standard substance)      eggs ..... approx. 11 µg/kg (ppb)  
                                 feed..... approx. 82 µg/kg (ppb)  
                                 urine ..... approx. 23 µg/kg (ppb)

Recovery rate:      milk ..... approx. 99 %  
(corresponding to the      meat..... approx. 122 %  
standard substance)      eggs ..... approx. 110 %  
                                 feed..... approx. 87 %  
                                 urine ..... approx. 84 %

The specificity of the RIDASCREEN® Bacitracin test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Specificity:	Bacitracin (standard substance).....	100 %
	Virginiamycin .....	< 3 %
	Tylosin .....	< 2 %
	Spiramycin, Neomycin, Colistin.....	< 1 %

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

## 1. Intended use

RIDASCREEN® Bacitracin is a competitive enzyme immunoassay for the quantitative analysis of bacitracin in milk, meat, eggs, feed and urine.

## 2. General

Bacitracin belongs to the group of polypeptide antibiotics. It was named after the *Bacillus* strain producing it and the 7 year old girl 'Tracy' from whose open tibial fracture the strain was isolated in 1945 for the first time. Bacitracin is a mixture of different polypeptides, wherein bacitracin A as major component has the highest biological activity. Through inhibition of cell wall synthesis, it has a broad spectrum bactericidal effect against both gram-positive and gram-negative bacteria.

In addition to veterinary application, bacitracin can be used as antimicrobial growth promoter in animal husbandry. Thereby, residues can remain in foodstuff of animal origin and may pose a health risk to consumers. Improper use of antibiotics also promotes the formation of antibiotic resistant bacteria, which raise an increasing health problem for the population. Consequently, Bacitracin has been deleted by the regulation 2821/98/EC from the list of approved animal feed additives (Directive 70/524/EEG Appendix B) and maximum residue limits in foodstuffs have been set (Regulation 37/2010/EC).

### 3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with antibodies directed against bacitracin. Bacitracin standards or sample and bacitracin conjugate are added. Free bacitracin and bacitracin conjugate compete for the bacitracin antibody binding sites (competitive enzyme immunoassay). Any unbound conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the bacitracin concentration in the sample.

### 4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 wells
Sample buffer	White	Ready to use		40 ml
Standard 1	White	Ready to use	0 µg/L	2 ml
Standard 2	White	Ready to use	0.625 µg/L	1 ml
Standard 3	White	Ready to use	1.25 µg/L	1 ml
Standard 4	White	Ready to use	2.5 µg/L	1 ml
Standard 5	White	Ready to use	5 µg/L	1 ml
Standard 6	White	Ready to use	10 µg/L	1 ml
Standard 7	White	Ready to use	20 µg/L	1 ml
Wash buffer	White	Concentrate	20x	30 ml
Conjugate	Red	Concentrate	100x	0.1 ml
Conjugate buffer	Green	Ready to use		15 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		12 ml
Stop solution	Yellow	Ready to use		15 ml

## 5. Materials required but not provided

### 5.1. Equipment:

Equipment	milk	meat, eggs, feed	urine
microtiter plate spectrophotometer (450 nm)	•	•	•
graduated pipettes	•	•	•
variable 20 µl - 200 µl and 200 - 1000 µl micropipettes	•	•	•
pH-meter	•		
mixer		•	
shaker		•	
vortex	•		
centrifuge	•	•	

### 5.2. Reagents:

Reagent	milk	meat, eggs, feed
0.1 M NaOH for pH adjustment	•	
80 % (v/v) methanol in sample buffer		•

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen prior to test implementation
- a value of less than 0.6 absorbance units ( $A_{450\text{ nm}} < 0.6$ ) for the zero standard

## 9. Preparation of Samples

The samples should be stored in a cool place.

### 9.1. Milk

- adjust pH to  $7 \pm 0.5$  with 0.1 M NaOH
- centrifuge: 5 min / 2000 g / 4 °C (39.2 °F)
- remove upper fat layer
- dilute defatted milk sample 1:10 (1+9) in sample buffer and vortex (e.g. 20 µl defatted milk + 180 µl sample buffer)
- use 50 µl per well in the assay

### 9.2. Meat, eggs and feed

- homogenize sample completely
- add 2 ml of 80 % methanol in sample buffer to 1 g of homogenized sample
- shake: 15 min / over-head
- centrifuge: 10 min / 2000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- meat and eggs: dilute supernatant 1:5 (1+4) with sample buffer and vortex (e.g. 40 µl supernatant + 160 µl sample buffer)
- feed: dilute supernatant 1:20 (1+19) with sample buffer and vortex (e.g. 10 µl supernatant + 190 µl sample buffer)
- use 50 µl per well in the assay

### 9.3. Urine

- dilute urine 1:25 (1+24) in sample buffer (e.g. 20 µl urine + 480 µl sample buffer)
- use 50 µl per well in the assay

## 10. Test implementation

### 10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

Eventually occurring precipitates in the concentrates have to be dissolved by shaking at room temperature before dilution.

Dilute **wash buffer** concentrate 1:20 (1+19) with demineralized water (e.g. 2 ml concentrate + 38 ml demineralized water, sufficient for 4 microtiter strips). Diluted wash buffer can be stored at 2 - 8 °C (35 - 46 °F) until the shelflife of the concentrate expires.

The **bacitracin conjugate** is provided as a concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is needed should be reconstituted. Before dilution, centrifuge the conjugate concentrate in the vial shortly (1 min / 1000 g / room temperature (20 - 25 °C / 68 - 77 °F)). For reconstitution, the concentrate is diluted 1:100 (1+99) in conjugate buffer (e.g. 20 µl conjugate concentrate + 1980 µl conjugate buffer, sufficient for 4 microtiter strips).

### 10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µl of standard or prepared sample to separate duplicate wells and add 50 µl of diluted conjugate to each well. Mix gently by shaking the plate manually and incubate for 1 h at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all wells with 300 µl of wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
4. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
5. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm immediately after addition of stop solution.

## 11. Results

A special software, the RIDA<sup>®</sup>SOFT Win / RIDA<sup>®</sup>SOFT Win.net (Art. No. Z9996), is available to evaluate the RIDASCREEN<sup>®</sup> enzyme immunoassays.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the bacitracin concentration [ $\mu\text{g}/\text{kg}$ ].

In order to obtain the bacitracin concentration in  $\mu\text{g}/\text{kg}$  (ppb) actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factor is as follows:

milk .....	10
meat, eggs.....	15
feed.....	40
urine.....	25

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG