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## Brief information

RIDASCREEN<sup>®</sup>FAST Folsäure (Folic Acid) (Art. No. R3203) is a competitive enzyme immunoassay for the quantitative determination of added folic acid in milk, milk powders, food for special medical purpose\*, grain and cereals, fortified flour, vitamin powders, -mixtures, -tablets, and vitamin juices.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit. The test kit is sufficient for 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:           milk: can be applied directly, dilution if necessary  
  juices: can be applied directly, dilution if necessary  
  other matrices: extraction, dilution if necessary

Time requirement:           sample preparation (for 10 samples).....10 - 60 min  
  test implementation (incubation time).....25 min

Detection limit:           0.5 µg/L (ppb) (in measurement solution\*\*) (corresponding to standard substance)

The specificity of the RIDASCREEN<sup>®</sup>FAST Folic test was determined by analyzing the cross-reactivities to corresponding substances in a buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Specificity:

Folic Acid .....	100 %
Dihydrofolic Acid .....	approx. 14 %
Tetrahydrofolic Acid .....	< 0.5 %
5-Methyltetrahydrofolic Acid .....	< 0.5 %
Formyltetrahydrofolic Acid .....	< 0.5 %
Tetrahydrobiopterin.....	< 0.5 %
6-Biopterin.....	< 0.5 %

\* Commission Directive 1999/21/EC on dietary foods for special medical purposes

\*\* see product information

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) - Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

## **1. Intended use**

RIDASCREEN®FAST Folsäure (Folic Acid) is a competitive enzyme immunoassay for the quantitative determination of added folic acid in milk, milk powder, food for special medical purpose, grain and cereals, fortified flour, vitamin powder, -mixture, -tablets, and vitamin juice.

## **2. General**

Folic acid also known as vitamin B9 or vitamin B11 belongs to the group of water soluble vitamins. It can't be produced by the human body and therefore must be included in the diet in order to meet the daily requirements. Folic acid plays a key role in all growth and cell division processes. Adequate folic acid intake during pregnancy is important to avoid neural tube defects, for example. Some countries have agreed to fortify either bread, cereals, flour, pasta, rice or other grain products with folic acid in different concentrations as a preventive approach to improve micronutrient status of the population.

## **3. Test principle**

The basis of the test is the antigen-antibody reaction. The microtiter plate is coated with specific antibodies against folic acid. Folic acid standards, respectively sample and enzyme labelled folic acid (conjugate) are added. Free and enzyme conjugated folic acid compete for the antibody binding sites (competitive enzyme immunoassay). Any unbound conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The resulting absorbance values are inversely proportional to the folic acid concentration of the sample.

## 4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses). Each test kit contains:

Component	Cap Color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
Sample buffer	Transparent	2-fold concentrate	2x	125 ml
Standard 1	White	Ready to use	0 µg/L (ppb)	1.3 ml
Standard 2	White	Ready to use	0.75 µg/L (ppb)	1.3 ml
Standard 3	White	Ready to use	1.5 µg/L (ppb)	1.3 ml
Standard 4	White	Ready to use	3 µg/L (ppb)	1.3 ml
Standard 5	White	Ready to use	6 µg/L (ppb)	1.3 ml
Standard 6	White	Ready to use	12 µg/L (ppb)	1.3 ml
Wash buffer salt Tween		Salt for dissolving		
Conjugate	Red	Ready to use		3 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 ml
Stop solution	Yellow	Ready to use		14 ml

## 5. Materials required but not provided

### 5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- laboratory mincer/grinder, mortar, ultra-turrax or homogenizer
- centrifuge + centrifugal vials
- water bath (100 °C / 212 °F)
- ice bath (0 °C / 32 °F)
- shaker
- optional: paper filter and funnel
- graduated pipettes
- variable 20 - 200 µl and 200 - 1000 µl micropipettes
- optional: multichannel pipette or 8 channel pipette

### 5.2. Reagents:

- distilled or demineralized water

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant bag provided and further store at 2 - 8 °C (36 - 46 °F).

The standards, the conjugate and the substrate/chromogen are light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers. (exception: sample buffer, wash buffer (salt), stop solution)

## 8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen prior to test implementation
- an absorption smaller than 0.6 ( $A_{450\text{ nm}} < 0.6$ ) for the zero standard

## 9. Preparation of samples

The samples should be stored in a cool place, protected from light. The samples should be brought to room temperature before analysis (20 - 25 °C / 68 - 77 °F). For quality control it is recommended to analyze control samples.

### Remark:

**The samples should be measured at the day of extraction.**

**The sample buffer can also be used for the sample preparation of vitamin B12 containing samples in combination with the RIDASCREEN® FAST Vitamin B12 (R2103).**

## 9.1. Milk, milk powder and food for special medical purpose

### Milk:

- homogenized milk can be directly applied in the test, dilute if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- add 50 µl per well in the assay

### Milk powder, food for special medical purpose:

- suspend 1 g sample in approx. 5 ml distilled water, fill up to 10 ml with distilled water (dilution factor = 10)
- shake / mix for 10 min
- heat the diluted sample for 3 min at 100 °C (212 °F) in a water bath, cool down quickly (ice bath), dilute the supernatant or filtrate if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- add 50 µl per well in the assay

## 9.2. Grain and cereals

A representative sample should be ground and thoroughly mixed before proceeding with the extraction procedure (particle size < 250 µm).

- weigh 1 g of ground sample into a suitable vial, suspend with approx. 5 ml sample buffer and fill up to 10 ml with sample buffer (dilution factor = 10)
- shake / mix for 10 min
- centrifuge (or alternatively filter) in order to remove undissolved components
- take supernatant, if supernatant is still turbid, filtrate
- dilute the supernatant or filtrate if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- add 50 µl per well in the assay

## 9.3. Fortified flour

A representative sample should be ground and thoroughly mixed before proceeding with the extraction procedure (particle size < 250 µm)

- weigh 1 g of the homogenized flour sample into a suitable vial, suspend with approx. 5 ml sample buffer and fill up to 10 ml with sample buffer (dilution factor = 10)
- shake / mix for 10 min

- centrifuge (or alternatively filter) in order to remove undissolved components
- take supernatant, if supernatant is still turbid, filtrate
- dilute the supernatant or filtrate if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- add 50 µl per well in the assay

#### 9.4. Vitamin powder, -mixture, and -tablets

##### Vitamin powder / -mixture:

- dissolve 1 g of vitamin powder / -mixture in 5 ml sample buffer and fill up to 10 ml with sample buffer (dilution factor = 10)

##### Tablets:

- crush the tablets and dissolve 1 g in 5 ml sample buffer and fill up to 10 ml with sample buffer (dilution factor = 10)

Continue with these samples as is described below:

- shake / mix for 10 min
- centrifuge (or alternatively filter) in order to remove undissolved components
- dilute the supernatant or filtrate if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- add 50 µl per well in the assay

#### 9.5. Vitamin juice

- the homogenized juice can be used directly in the test
- centrifuge (or alternatively filter) in order to remove undissolved components
- dilute if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- add 50 µl per well in the assay .

## 10. Test implementation

### 10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **sample buffer** is provided as a 2-fold concentrate. The concentrate is diluted 1:2 with distilled water (e.g. 50 ml sample buffer concentrate + 50 ml distilled

water). The diluted sample buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

A PBS-Tween buffer is required as **wash buffer**, please use the wash buffer salt (pouch) contained in the kit (see 4.). Dissolve the total content of the pouch in one liter of distilled water. The ready to use wash buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternatively: Dissolve the content of the pouch in 100 ml of distilled water to obtain a 10-fold concentrated wash buffer. This solution expires after approx. 8 - 12 weeks, stored at room temperature (20 - 25 °C / 68 - 77 °F). Use one part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use washing buffer.

## 10.2. Test procedure

Standards and samples must be tested on one microtiter plate at one time. Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µl of each standard or prepared sample to the wells in duplicate; use a new pipette tip for each standard or sample.
3. Add 50 µl of conjugate to each well (multichannel pipette or 8 channel pipette) and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Using a 8 channel pipette, fill the wells with ready to use wash buffer (250 µl/well) (10.1). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.
5. Add 100 µl of the substrate/chromogen to each well (multichannel pipette or 8 channel pipette); and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
6. Add 100 µl of stop solution to each well (multichannel pipette or 8 channel pipette); and measure the absorbance at 450 nm within 5 min.

## 11. Results

A special software, the RIDA<sup>®</sup>SOFT Win / RIDA<sup>®</sup>SOFT Win.net (Art. No. Z9996) is available for evaluation of the RIDASCREEN<sup>®</sup> enzyme immunoassays. The results are given in µg/kg respectively µg/L and µg/100 g respectively µg/100 ml.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semi-logarithmic graph paper against the folic acid concentration (µg/kg or µg/L).

In order to obtain the folic acid concentration actually contained in a sample, the concentration read from the standard curve must be further multiplied by the corresponding dilution factor.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG