RIDA[®] QUICK Zearalenon RQS

Art. No.: R5504

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RIDA[®]QUICK SCAN Methode: **Zea Quant** or **Zea75-500ppb**

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RIDA[®] QUICK Zearalenon RQS

Brief information

RIDA[®]QUICK Zearalenon RQS (Art. No.: R5504) is an immunochromatographic test for the determination of zearalenone in corn.

All reagents required for the assay are contained in the test kit.

Results are evaluated with RIDA[®]QUICK SCAN (RQS, Art. No.: ZG5005) or the RIDA[®]SMART APP software (Art. No. ZRSAM1000).

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples) approx. 10 min
test implementation (incubation time) 5 min

Limit of detection approx. 75 µg/kg (ppb; evaluation with RIDA[®]QUICK SCAN)
approx. 50 µg/kg (ppb; evaluation with RIDA[®]SMART APP)

Specificity: The RIDA[®]QUICK Zearalenon RQS test determines zearalenone in corn samples. The test has been validated for this matrix.

The specificity of the RIDA[®]QUICK Zearalenon RQS test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Related products

RIDASCREEN[®] Zearalenon (Art. No. R1401)

RIDASCREEN[®]FAST Zearalenon (Art. No. R5502)

RIDASCREEN[®]FAST Zearalenon SC (Art. No. R5505)

1. Intended use

RIDA[®]QUICK Zearalenon RQS is a quantitative immunochromatographic test in strip format for the determination of zearalenone in corn.

2. General

The mycotoxin zearalenone is formed by fungi of the genus *Fusarium*.

Zearalenone is a phytohormone which displays, apart from its anabolic properties, mainly estrogenic effects. Because of its estrogenic properties, zearalenone may induce fertility disorders in animals with clinical signs of hyperestrogenism - an aspect of a disease which although reported mainly in hogs, is described in other species such as cow, horse and sheep, too.

The potential health risk for humans induced by this mycotoxin, which is taken up with foods of vegetable or animal origin, is extensively discussed.

3. Test principle

The basis of the immunochromatographic assay in test strip format is an antigen-antibody reaction. A specific antibody against zearalenone binds the zearalenone molecules in the samples.

The control band (control line) should be present in all cases in order to prove that the test strip is valid. The control band intensity decreases with increasing zearalenone concentration in the sample. When a sample contains very high levels of zearalenone, there is only a very weak control band visible.

Results are evaluated with RIDA[®]QUICK SCAN (RQS) reader or the RIDA[®]SMART APP software.

4. Reagents provided

Each kit contains sufficient materials for 20 determinations. Each test kit contains:

Component	Cap color	Format	Content
20 x Test strip	-	Ready to use	One determination each (separately packed)
1 x Mobile solvent	White	Ready to use	13 ml
1 x RIDA [®] SMART APP cover	-	Ready to use	Lot-specific evaluation cover for the use with RIDA [®] SMART APP

5. Materials required but not provided

5.1. Equipment:

- laboratory mincer / grinder
- balance
- optional: shaker
- optional: centrifuge
- graduated cylinder
- 100 µl, 600 µl micropipettes (e.g. R-Biopharm FP 100, Art. No.: Z0007)
- RIDA[®]QUICK SCAN (Art. No.: ZG5005) or activated RIDA[®]SMART APP software (Art. No. ZRSAM1000) installed on compatible smartphone

5.2. Reagents:

- 70 % methanol (e.g. 70 ml 100 % methanol / 30 ml distilled water), exact adjustment necessary

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The test strips are sensitive to humidity. Humid test strips may influence the test results negatively; therefore keep the strips away from humidity. This has to be noted especially for already opened reaction strip packing.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.de.

7. Storage instructions

Store the kit at 2 - 8 °C / 36 - 47 °F. Do not freeze test kit components.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Sample preparation

Allow all components and liquids to reach room temperature before use (20 - 25 °C / 68 - 77 °F)!

The samples should be stored in a cool place, protected from light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

8.1. Sample preparation for evaluation with RIDA[®]QUICK SCAN reader or RIDA[®]SMART APP software

- weigh 10 g of ground sample into a screw cap tube and add 20 ml of 70 % methanol *)
- close the tube and shake the sample vigorously for approx. 5 min (manually or with shaker)
- keep tube closed to avoid methanol evaporation and let the solution settle for approx. 3 - 5 min or centrifuge
- mix 100 µl of the clear supernatant with 600 µl mobile solvent
- apply 100 µl of diluted sample onto the application area of the test strip

*) Sample size may be increased if required, but the volume of methanol (70%) must be adapted accordingly, e.g. 50 g in 100 ml of 70 % methanol.

9. Test procedure

- apply 100 µl of the sample extract (see 8. Sample preparation) onto the application area of the test strip (see Fig. 1)
- read the result after 5 min incubation (+ max. 15 sec) with RIDA[®]QUICK SCAN or RIDA[®]SMART APP


10. Evaluation

The right band in the test strip reaction field is a control band (control line) and must appear after each test procedure (see Fig. 1). If the band is missing, the test result is not valid because of improper test procedure or deterioration of the reagents. Repeat the test with a new strip. If the right band is missing again, please inform your local distributor.

Please note: If the sample is highly contaminated with zearalenone the visibility of the control band is very weak and the test may give an invalid result.

To obtain comparable results, the test strips must always be evaluated after 5 min.

Fig.1: Test strip RIDA[®]QUICK Zearalenon RQS

10.1 Evaluation with RIDA[®]QUICK SCAN (Art. No. ZG5005)

First, please read the instructions for use of the RIDA[®]QUICK SCAN reader attentively.

Please find a short description regarding the use of RIDA[®]QUICK SCAN in the test kit.

RIDA[®]QUICK SCAN method: **Zea Quant** or **Zea75-500ppb**

10.2. Evaluation with RIDA[®]SMART APP (Art. No. ZRSAM1000)

First, please read the instructions for use of the RIDA[®]SMART APP software attentively. The test strip evaluation has to be performed according to these guidelines.

Please find a short description of how to use the RIDA[®]SMART APP on our website www.r-biopharm.de.

For the test strip evaluation the RIDA[®]SMART APP cover is needed which is included in the test kit. The RIDA[®]SMART APP cover contains lot-specific

information and must be used only with the appropriate lot of the test kit. For further information please check the certificate of analysis enclosed in the kit. Before starting the test implementation and test evaluation please check if the correct RIDA[®]SMART APP cover is available.

RIDA[®]SMART APP software application: **50-1000ppb**

11. Sensitivity

The RIDA[®]QUICK Zearalenon RQS test is capable of detecting zearalenone contaminations of $\geq 75 - 500 \mu\text{g/kg}$ if test strip evaluation is performed with RIDA[®]QUICK SCAN. Zearalenon contaminations of $\geq 50 - 1000 \mu\text{g/kg}$ can be detected when evaluating with the RIDA[®]SMART APP.

For further information please contact your local distributor or R-Biopharm AG (sales@r-biopharm).

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG