**Vitamin C**  
**(L-Ascorbinsäure / L-Ascorbic Acid)**

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**IFUMA17A:** Determination of ascorbic acid by HPLC (1987; rev. 2005)

**Beutler, H. O. (1984)** in Methods of Enzymatic Analysis (Bergmeyer, H. U., ed.) 3rd ed., Vol. VI, pp. 376-385, Verlag Chemie, Weinheim, Deerfield Beach / Florida, Basel

## 11. Abkürzungen

AO	Ascorbat-Oxidase
E	Extinktion
$\Delta E$	Extinktionsdifferenz
IFU	International Federation of Fruit Juice Producers
IFUMA	International Federation of Fruit Juice Producers - Method of Analysis
KOH	Kaliumhydroxid
MTT	[3-(4,5-Dimethylthiazol-2-yl)-2,5-diphenyltetrazoliumbromid]
PMS	(5-Methylphenaziniummethylsulfat)

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# VitaFast<sup>®</sup> Vitamin C (L-Ascorbic Acid)

## Brief information

Easy to conduct enzymatic test in microtiter plate format for the quantitative determination of vitamin C (L-ascorbic acid) in food, pharmaceutical products and other sample material. Additionally L-dehydroascorbic acid can be detected (see application note). The kit contains all the required reagents including an L-ascorbic acid standard for the standard curve. 50 reactions (plus four calibration curves) can be carried out. For evaluation you will need a microtiter plate photometer (550 - 570 nm).

## Sample preparation

Homogenise, extract, filter, dilute

## Time required for five samples

Sample preparation:	approx. 10 - 30 min
Conducting the test:	approx. 60 min
Evaluation:	approx. 5 min

## Test characteristics

Test volume:	320 µl
Measurement in:	microtiter plate, 96 wells
Wavelength:	550 - 570 nm
Temperature:	20 - 25 °C
Weighed sample:	1 - 20 g (ml)
Sample volume:	200 µl
Blank measurement:	meta-phosphoric acid
Standard range:	7.8 - 250 mg / 100 g (ml)
Detection limit	= quantification limit
Quantification limit:	for 1 g (ml) weighed sample: 7.8 mg / 100 g (ml)
Recovery rates:	95 - 105 % (with freshly prepared ascorbic acid)

## 1. Intended use

**VitaFast®** Vitamin C (L-Ascorbic Acid) is a test in microtiter plate format for the quantitative determination of vitamin C (L-ascorbic acid) in foods, pharmaceutical products and other sample material. Furthermore the total amount of vitamin C (L-ascorbic acid + L-dehydroascorbic acid) can be detected (see application note).

Sample preparation is carried out according to the HPLC method presented by the International Federation of Fruit Juice Producers - Methods of Analysis (IFUMA) No. 17a (Rev.2005).

### Information on L-ascorbic acid:

**Ascorbic acid** is a clear, odourless, crystalline solid substance which has a sour flavour. There are four stereoisomeric variants but only L-(+)-ascorbic acid (also referred to as vitamin C) is biologically active. Its vitamin activity is the most relevant property of L-ascorbic acid: humans cannot produce vitamin C and require it as part of their nutrition. A deficiency can result in the manifestation of scurvy.

L-ascorbic acid and its salt, the anion L-ascorbate, will hereinafter be referred to as ascorbic acid and ascorbate, respectively.

## 2. Principle of the test

Two assays are performed to determine ascorbic acid (which is present in the sample extract as ascorbate) based on its property of being a reductant that can be oxidised by donating an electron to an oxidising agent (oxidant):

### 1. Sample blank assay - with Ascorbate Oxidase (+AO)

In the first assay the ascorbate fraction of the sample is oxidised to dehydroascorbate by ascorbate oxidase (AO) in the presence of oxygen (Equation (1)).

At pH 3.5, the remaining reductants in the assay reduce the tetrazolium salt MTT\* to MTT-Formazan in the presence of the electron carrier PMS\*\* (Equation (2)). Dehydroascorbate does not react with MTT / PMS; hence **only the non-ascorbate reductants** are measured in this assay.

### 2. Sample assay - without Ascorbate Oxidase (-AO)

In the second assay the AO step is omitted. Thus, the **total sum of reductants including ascorbate** is measured according to Equation (2).

\* MTT = [3-(4,5-dimethylthiazolyl-2)-2,5-diphenyltetrazolium bromide]


\*\* PMS = (5-methylphenazinium methosulphate)

The **formazan** formed in the assays is proportional to the sum of reductants and is measured at 550 - 570 nm in a microtiter plate photometer. The ascorbic acid content can be calculated based on the absorbance difference


of the sample assay minus the absorbance difference of the sample blank assay. Results are evaluated based on an ascorbic acid standard curve.

The following reactions take place:

### Sample blank assay - with Ascorbate Oxidase (+AO)

### Sample assay - without Ascorbate Oxidase (-AO)

## 3. Contents of the test kit

Bring all solutions to 20 - 25 °C prior to use and blend by carefully swirling around prior to pipetting. Each test kit contains:

1 x Microtiter plate with 96 wells

1 x Solution 1

15 ml phosphate / citrate buffer (pH approx. 3.5) with MTT

2 x Bottle 2 - red screw cap

Ascorbate oxidase (AO, lyophilisate), approx. 250 U each

**Fully dissolve** the content of **one Bottle 2** with **0.2 ml sterile water** from the test kit (see below). Leave at room temperature for 10 min, swirl several times. The solution (= **Solution 2**) will remain stable for eight weeks if stored at 2 - 8 °C. Alternatively, Solution 2 can be aliquoted and stored at -20 °C for three months.

2 x Bottle 3 - green screw cap

PMS (lyophilisate - Keep dark!)

**Fully dissolve** the content of **one Bottle 3** with **1.5 ml sterile water** from the test kit (see below). Leave at room temperature for 10 min, swirl several times. The solution (= **Solution 3**) will remain stable for eight weeks if stored at 2 - 8 °C. Alternatively, Solution 3 can be aliquoted and stored at -20 °C for three months.

5 x Bottle 4 - white screw cap

Ascorbic acid standard (lyophilisate)

**Fully dissolve** the content of **one Bottle 4** with **2 ml 1.5 % meta-phosphoric acid, pH 3.5 (= MPA)**, directly before use. The solution (= **Solution 4**) corresponds to the highest standard (250 mg / 100 g [ml]) and must be used **fresh only**. Discard the remaining solution after conducting the test.

1 x Sterile water (30 ml)

For dissolving the lyophilisate in Bottles 2 and 3.

Note: Push the coloured lid up, pull off right up to the glass rim and turn entire lid to remove it.

Remark: after the expiry date no guarantee of quality

## 4. Required devices and reagents, not provided

### 4.1 Devices and accessories

- Analytical scale
- Microtiter plate photometer (550 - 570 nm wavelength)
- pH meter
- Magnetic stirrer
- 100 ml beakers
- 100 ml volumetric flask (brown glass)
- 100 ml Schott glass bottles (brown glass)
- Funnel and pleated filter
- 1.5 and 2.0 ml reaction vials
- Micropipettes and pipette tips 10 - 200 µl; 100 - 1000 µl

### 4.2 Reagents

- Distilled or deionised water (**referred to as water hereinafter**)
- Meta-phosphoric acid (1.5 % w/v, pH adjusted to 3.5 with KOH, approx. 10 mol / l) (**referred to as MPA hereinafter**)
- HCl, 1 mol / l (dilute stock solution according to supplier's instructions)
- KOH, 10 mol / l (dissolve 56.11 g potassium hydroxide in approx. 70 ml water (cool while dissolving in water) and top up to 100 ml with water)
- KOH, 2 mol / l (dilute the 10 mol / l solution 1 : 5 with water)

## 5. Storage instructions

Store the test kit / the reagents at 2 - 8 °C. See Section 3 for storage instructions applicable to dissolved reagents. Store PMS in the dark.

## 6. Sample preparation

Samples are extracted in MPA (1.5 %, pH 3.5).

**Use the samples prepared in this way directly for the assay.**

1. Weigh 1 - 20 g (ml) of homogenized sample into a 100 ml beaker.
- 2.\* Dissolve in approx. 50 ml MPA using a magnetic stirrer (approx. 15 min or until completely dissolved).
3. Adjust pH to 3.5 - 4.0 with HCl (1 mol / l) or KOH (2 mol / l).
4. Transfer quantitatively to a 100 ml volumetric flask (brown glass).
5. Top up to 100 ml with MPA, mix and filter into a 100 ml Schott glass bottle (brown glass) through a pleated filter; discard the first millilitres.
6. Use clear solution for the test.
7. If necessary, dilute the sample extract with MPA in 1.5 ml reaction vials, depending on the ascorbic acid concentration (= **dilution factor F**).

\* Dilute liquid samples with approx. 50 ml MPA.

### Notes:

- Alternatively, samples can be weighed directly into the 100 ml volumetric flask, provided that the sample dissolved in MPA (1.5 %, pH 3.5) stays at pH 3.5 - 4. The pH of the dissolved sample depends on the matrix and on the sample weight. If the standard sample preparation procedure shows that the pH of a specific sample is within the range of 3.5 - 4, or can be adjusted to this range by adding a defined volume of HCl or KOH, the pH electrode is not necessary for step 3 and the sample can be weighed directly into the volumetric flask. Step 4 (quantitative transfer) is omitted.
- If the solution is not clear at step 6, a fine mesh filtration using a 0.22 - 0.45  $\mu\text{m}$  filter is recommended (syringe filter CME or polyester syringe filters).
- Matrices that are not easy to filter, such as milk powder, can be centrifuged for 5 min at  $\geq 4,000 \times g$ . Filter the entire supernatant through a pleated filter afterwards to eliminate remaining particles.

### Tips for sample dilution:

- To calculate the required dilution, the **expected content** (in mg / 100 g [ml]) is **multiplied** by the **sample weight** (in g) and **divided** by the concentration of the **second standard** (125 mg / 100 g (ml)).
- Generally, the following applies to the analysis of diluted sample extracts: the absorbance difference  $\Delta A$  measured in the sample assay (i. e. the total reductants assay without AO) should be  $\geq 0.1$ . If  $\Delta A$  is smaller, the dilution is too strong. In that case it is recommendable to repeat the analysis with a less strongly diluted sample extract or with an increased amount of weighed sample.

## 6.1 Examples of sample extract dilution

- **Fruit gum**, expected vitamin C content: 18 mg / 100 g  
**Weighed sample:** 7 g  
**Dilution factor F** of the sample extract: 1 (undiluted)
- **Multivitamin effervescent tablet**,  
expected vitamin C content: 1,300 mg / 100 g  
**Weighed sample:** 1 g  
**Dilution factor F** of the sample extract: dilute to a ratio of 1 : 10 (F = 10)
- **Iron and vitamin C tablet**, expected vitamin C content: 30 g / 100 g  
**Weighed sample:** 1 g  
**Dilution factor F** of the sample extract: dilute to a ratio of 1:240 (F = 240)

## 7. Test procedure / assay setup

### 7.1 Dilution of the standard

**Dissolve 1 x Bottle 4** (ascorbic acid lyophilisate) with 2 ml MPA (1.5 %, pH 3.5). The resulting solution (= Solution 4) corresponds to an ascorbic acid concentration of 250 mg / 100 g (ml) and can be used directly as the highest standard (**standard 4**).

For the lower standards, **dispense MPA**( 1.5 %, pH 3.5) into 1.5 or 2.0 ml reaction vials, as specified in Table 1. Continue as shown in Table 1.

**Table 1:** Dilution of the ascorbic acid standard

<b>Standard series mg / 100 g (ml)</b>	<b>MPA (1.5 %, pH 3.5)</b>		<b>Standard</b>
standard 4: 250	-		standard 4 = Solution 4
standard 3: 125	500 µl	+	500 µl from 4
standard 2: 31.3	600 µl	+	200 µl from 3
standard 1: 7.8	600 µl	+	200 µl from 2
blank: 0	blank = MPA		-

#### **Note:**

The standard series corresponds to the specified concentrations under standard test conditions, i. e. with a weighed sample of 1 g in 100 ml and without further dilution of the sample extract, the sample's ascorbic acid content can be read directly from the calibration curve. Only if the amount of weighed sample is increased, or if the extract is diluted must this be taken into account in the calculation (see Section 8).

## 7.2 Preparation of solutions

In order to determine the ascorbate fraction of the total reducing potential of the sample, ascorbic acid must be removed oxidatively by ascorbate oxidase in the sample blank assay. For each required well of the microtiter plate, **mix 95 µl of Solution 1** and **5 µl of Solution 2** (red screw cap) in a 1.5 ml reaction vial (see Table 2 for examples). The mixture must be prepared freshly prior to each assay setup and is discarded afterwards.

**Table 2:** Premix Solutions 1 + 2 for the required number of samples (a pipetting reserve is already included in the calculation below)

Number of samples	Solution 1 µl		Solution 2 µl	Total µl
1	190	+	10	200
2	285	+	15	300
3	380	+	20	400
4	475	+	25	500
5	570	+	30	600
6	665	+	35	700
7	760	+	40	800
8	855	+	45	900
9	950	+	50	1,000
10	1,045	+	55	1,100
15	1,615	+	85	1,700
20	2,090	+	110	2,200

### Note:

Solution 2 is not needed for the determination of the total sum of reductants. Hence, it is not necessary to mix Solutions 1 and 2 for the sample assay.

## 7.3 Assay setup

The individual steps for preparing the assays are shown in **Table 3**. You need to prepare two assays:

### 1. Sample blank assay (+AO)

The sample blank assay (determination of non-ascorbate reductants after oxidative removal of ascorbic acid) is prepared **with** ascorbate oxidase (step 1.1, Solution 1 + 2).

### 2. Sample assay (-AO)

The sample assay (determination of the total sum of reductants) is prepared **without** ascorbate oxidase (step 1.2, only Solution 1).

### Note:

- If the VitaFast<sup>®</sup> Vitamin C excel evaluation sheet is used for evaluation after the measurements (see 8.1.1.), it is recommended to begin pipetting by starting with the blank in well A1, followed in a downward direction by the standards in increasing concentration (i.e. standard 1 in well B1 all the way through standard 4 in well E1).

- All samples need to be performed twice: once in the sample blank assay and once in the sample assay.
- Standards need to be performed only once - without AO.
- Always rinse the pipette tip with the solution to be pipetted first and make sure it contains no air bubbles.
- Prior to each absorption measurement, mix the assay in the microtiter plate by gently swirling the plate on the tabletop.
- After the measurements, rinse the microtiter plate thoroughly with tap water followed by deionised or distilled water, and then leave to air dry before using again.

**Table 3:**

Pipetting scheme for the determination of non-ascorbate reductants (sample blank, +AO) and total reductants (sample, -AO).

Each column represents one microwell in which the single working steps starting from 1.1 or 1.2, respectively, are carried out in a downward direction all the way through step 8.

Step	Sample +AO	Sample -AO	Standard*	Notes
1.1 Solution 1 + 2	100 µl	-	-	Dispense onto the <b>bottom</b> of the wells.
1.2 Only Solution 1	-	100 µl	100 µl	
2. Standard*	-	-	200 µl	Add to the bottom of the wells, mix by pipetting up and down once.
3. Sample	200 µl	200 µl	-	
<b>4. Incubate in the dark for 20 min at 20 - 25 °C.</b>				
<b>5. Measurement A1</b>	A1 <sub>sample blank</sub>	A1 <sub>sample</sub>	A1 <sub>standard</sub>	Mix assay by swirling plate on the tabletop, measure absorbance at 550 - 570 nm immediately afterwards.
6. Solution 3 (green screw cap) → <b>Start of reaction</b>	20 µl	20 µl	20 µl	Add and mix by pipetting up and down three times; additionally use pipette tip to stir assay in well.
<b>7. Incubate in the dark for 30 min at 20 - 25 °C.</b>				
<b>8. Measurement A2**</b>	A2 <sub>sample blank</sub>	A2 <sub>sample</sub>	A2 <sub>standard</sub>	Mix assay by swirling plate on the tabletop, measure absorbance at 550 - 570 nm immediately afterwards.

\* Five separate wells: one for the blank and for each standard (1 - 4)

\*\* The reaction is complete when the measured absorbance remains constant; this is usually the case after 30 minutes.

## 8. Evaluation

The absorbance difference  $\Delta A = A2 - A1$  is calculated for the standards and samples in each assay. For each standard, the  $\Delta A$  from the **blank** is then **subtracted** from the  $\Delta A$  **standard** measured in the **sample assay**.

For each sample, the  $\Delta A$  from the **sample blank assay** is **subtracted** from the  $\Delta A$  measured in the **sample assay**. Table 4 shows a detailed evaluation of the test assay.

**Table 4:**

Evaluation scheme of the pipetting assays for non-ascorbate reductants (sample blank assay with AO) and total reductants (sample assay without AO)

	A1	A2	$\Delta A$
<b>Standard</b>			
blank	$A1_{\text{blank}}$	$A2_{\text{blank}}$	$(A2_{\text{blank}} - A1_{\text{blank}})$
7.8 mg / 100 g (ml)	$A1_{7.8}$	$A2_{7.8}$	$(A2_{7.8} - A1_{7.8}) - \Delta A_{\text{blank}}$
31.3 mg / 100 g (ml)	$A1_{31.3}$	$A2_{31.3}$	$(A2_{31.3} - A1_{31.3}) - \Delta A_{\text{blank}}$
125 mg / 100 g (ml)	$A1_{125}$	$A2_{125}$	$(A2_{125} - A1_{125}) - \Delta A_{\text{blank}}$
250 mg / 100 g (ml)	$A1_{250}$	$A2_{250}$	$(A2_{250} - A1_{250}) - \Delta A_{\text{blank}}$
<b>Sample</b>			
sample <sub>blank</sub>	$A1_{\text{sampleblank}}$	$A2_{\text{sampleblank}}$	$(A2_{\text{sampleblank}} - A1_{\text{sampleblank}})$
sample	$A1_{\text{sample}}$	$A2_{\text{sample}}$	$(A2_{\text{sample}} - A1_{\text{sample}}) - \Delta A_{\text{sampleblank}}$

The  $\Delta A$  values of the standards are plotted against the concentration of the standards in a **diagram**. Linear regression is used to determine the **linear equation** (3), which is then solved for ascorbic acid [mg / 100 g (ml)] (4). The values of  $\Delta A$  sample (=  $\Delta A$  without AO -  $\Delta A$  with AO) are then substituted into the equation in order to calculate the ascorbic acid content.

$$(3) \Delta A = \text{slope} \times \text{ascorbic acid [mg / 100 g (ml)]} + \text{y-intercept}$$

$$(4) \text{Ascorbic acid [mg / 100 g (ml)]} = \frac{(\Delta A \text{ sample} - \text{y-intercept}) \times F}{\text{slope} \times \text{weighed sample [g (ml)]}}$$

F = Dilution factor of the sample extract

### 8.1 Additional evaluation options

#### 8.1.1 Excel evaluation sheet


Along with the manual evaluation as per Section 8, the raw data can be entered in the VitaFast<sup>®</sup> Excel evaluation sheet for ascorbic acid. The results are then calculated automatically. The evaluation sheet is available on request.

## 8.2 Evaluation example

The following example illustrates the evaluation procedure based on multivitamin juice. The expected ascorbic acid content is 35 mg / 100 ml. The sample volume is 10 ml, and the sample extract was used in a 1 : 3 dilution. The raw data from the sample blank assay and the sample assay are shown in Table 5. The calibration curve calculated from this is shown in Fig. 1. The ascorbic acid content is calculated by substituting  $\Delta A_{\text{sample}}$  in Equation (4).

**Table 5:**  
Evaluation example of the pipetting assays

	A1	A2	$\Delta A$
<b>Standard</b>			
blank	0.047	0.050	$(A2_{\text{blank}} - A1_{\text{blank}}) = 0.003$
7.8 mg / 100 g (ml)	0.050	0.061	$(A2_{7.8} - A1_{7.8}) - \Delta A_{\text{blank}} = 0.008$
31.3 mg / 100 g (ml)	0.050	0.130	$(A2_{31.3} - A1_{31.3}) - \Delta A_{\text{blank}} = 0.077$
125 mg / 100 g (ml)	0.053	0.502	$(A2_{125} - A1_{125}) - \Delta A_{\text{blank}} = 0.446$
250 mg / 100 g (ml)	0.054	1.036	$(A2_{250} - A1_{250}) - \Delta A_{\text{blank}} = 0.979$
<b>Sample</b>			
Sample <sub>blank</sub>	0.085	0.109	$(A2_{\text{sampleblank}} - A1_{\text{sampleblank}}) = 0.024$
Sample	0.081	0.559	$(A2_{\text{sample}} - A1_{\text{sample}}) - \Delta A_{\text{sampleblank}} = 0.454$

**Fig. 1:**  
Standard curve for the determination of ascorbic acid  
(slope = 0.00404; y-intercept = -0.04078)

## ▪ Result for the ascorbic acid sample

$$(4) \text{ Ascorbic acid [mg / 100 ml]} = \frac{(0.454 - [-0.04078]) \times 3}{0.00404 \times 10 \text{ ml}} = 36.7$$

## 9. Test properties

### 9.1 Specificity

Under test conditions, the assay is specific for L-ascorbic acid. Additionally, isoascorbic acid (erythorbic acid, D-ascorbic acid), which does not occur naturally, is detected. In the European Union, isoascorbic acid is approved exclusively as a synthetic food additive (E 315) in certain meat and fish products and in heat-treated cream.

### 9.2 Quantification and detection limit

#### 9.2.1 Limit of detection

The detection limit equals the quantification limit (lowest standard).

#### 9.2.2 Limit of quantification

The quantification limit is defined as the concentration of the lowest standard.

Under standard test conditions (weighed sample of 1 g [ml], sample extract undiluted):

- 7.8 mg / 100 g (ml)

### 9.3 Recovery rates

With freshly prepared ascorbic acid (kept airtight and dark), recovery rates are 95 - 105 % on average. However, since the recovery of ascorbic acid from real food matrices depends on storage conditions, the quality of the solvent used for extraction / dilution, the age of the preparation / sample, its exposure to light / oxygen, and other factors, recovery rates below 95 % are not uncommon and do not indicate an insufficient enzymatic reaction efficiency.

## 10. Bibliography

The described sample preparation method complies with the requirements of the International Federation of Fruit Juice Producers - Methods of Analysis (IFUMA).

**IFUMA17A:** Determination of ascorbic acid by HPLC (1987; rev. 2005)

**Beutler, H. O. (1984)** in Methods of Enzymatic Analysis (Bergmeyer, H. U., ed.) 3rd ed., Vol. VI, pp. 376-385, Verlag Chemie, Weinheim, Deerfield Beach / Florida, Basel

## 11. Abbreviations

AO	Ascorbate oxidase
A	Absorbance
$\Delta A$	Absorbance difference
IFU	International Federation of Fruit Juice Producers
IFUMA	International Federation of Fruit Juice Producers - Method of Analysis
KOH	Potassium hydroxide
MTT	[3-(4,5 – dimethylthiazolyl-2)-2,5-diphenyltetrazolium bromide]
PMS	5-methylphenazinium methosulphate

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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is a registered trademark of the ifp Institut für Produktqualität GmbH.