**VitaFast®**  
**Folsäure / Folic Acid**  
**Dotierstandard / Spiking standard**

Art. No.: P3001

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## **VitaFast® Folic Acid Spiking standard**

### **Brief information**

Folic acid standard for the use of spiking food, feed and pharmaceutical products.

### **1. Reagents provided**

3 x folic acid spiking standard (solid)

Content: **see Quality Assurance Certificate**

Remark: after the expiry date no guarantee of quality.

### **2. Required reagents, not provided**

Redist. or deionized water for the preparation of the spiking standard solution.

### 3. Warning and precautions for the user

No precautions necessary.

### 4. Storage instructions

Store the reagents at 2 - 8 °C (35.6 - 46.4 °F).

### 5. Preparation of the spiking standard solution

- open the folic acid spiking solution bottle, place the lid down with the opening facing upwards
- add **x ml (x = see Quality Assurance Certificate „Spiking standard“ or label „Spiking standard“)** **redist. or deionized water**  
**content corresponds to 10 µg Folic Acid per 1 ml (or 1 µg per 100 µl)**
- close the bottle with the lid and dissolve the standard by shaking = spiking solution **(use directly)**
- add the needed volume to the sample (see example below) and prepare the sample according to the product information VitaFast® Folsäure / Folic Acid

The sample should be spiked with 50 % or 100 % of the vitamin content in the sample.

#### Example for a spiking procedure:

Estimated folic acid content in the sample: 400 µg / 100 g; take 1 g sample, this corresponds with a folic acid content of 4 µg / 1 g sample (the calculation refers to 1 g weighed sample in 40 ml).

#### 100 % spiking:

Spiking with 100 % (of the content of the sample):

Weigh 1 g sample + add 400 µl spiking solution (corresponds 4 µg absolute, respectively 400 µg / 100 g)

Calculation of the estimated total folic acid content:

Estimated content of the sample + content of the 100 % spiking  
400 µg / 100 g (estimated content) + 400 µg / 100 g (spike) =  
800 µg / 100 g folic acid (total content)

Prepare the spiked sample according to **6. Sample preparation** of the product information insert VitaFast® Folsäure / Folic Acid (Art. No. P1001).