QuickGEN Filtration

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QuickGEN Sample Preparation Filtration

1. Intended use

Fast and optimal DNA-extraction from bacteria and yeast out of filterable drinks without preenrichment.

2. Content

The kit contains reagents for 100 extractions

1 x QuickGEN buffer

2 x Washing solution

3. Materials required, but not provided

3.1 Instruments

- Membrane filtration equipment
- Pipettes
- "Vortex"

3.2 Reagents and plastic ware

- Polycarbonate-filter, 50 mm Ø, 0.4 µm
- for yeast DNA-extraction: Lyticase (Sigma L2524), in PCR-Kits with freeze-dried stripes the PCR-tubes contain lyticase (s. PCR Kit manual)
- sterile safe-lock reaction tubes 1.5 – 2.0 mL
- suitable filter tips
- single use gloves
- sterile water

4. Warnings

Careful use with personal protection according to good laboratory practice is recommended. Do not incorporate. Avoid skin and eye contact with all solutions.

5. Storage

Store all solutions at room temperature 15 – 30 °C (59 – 86 °F)

6. Indications of deterioration of reagents

In case of accurate handling deterioration unknown.

7. Preliminary preparations

Prepare Lyticase solution:

5U / μ L in lyticase buffer (50 % 1 x TE + 50 % Glycerol, pH-Wert 7.5 – 8.0).

Store solution at -20 °C (- 4 °F). **In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual).**

8. DNA-Extraction

Protocol:

1. filter the sample (0.33 L – 0.5 L) via a 0.4 μ m polycarbonate filter
2. wash the filter with 10 mL Washing solution
3. wash the filter with 30 mL sterile water
4. grasp the filter with sterile forceps and place it into the tube
5. add 250 μ L QuickGEN buffer
Optional: for yeast-DNA isolation add 15 U lyticase and mix well. In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual).
6. vortex the filter for 30 sec.
7. use 5 μ L of the DNA for PCR*

*** For PCR-analysis special PCR-Kits (QuickGEN) have to be used**