# **RIDA<sup>®</sup> Tetracyclin**

**Art. No. R3599**

Dotierlösung / Spiking Solution

In vitro Test

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Storage at 2 - 8 °C

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# RIDA<sup>®</sup> Tetracyclin Spiking Solution

## 1. Intended use

RIDA<sup>®</sup> Tetracyclin Spiking Solution can be used for the production of positive control samples suitable for the validation (e.g. determination of the recovery rate) of the RIDASCREEN<sup>®</sup> Tetracyclin test (R3505).

## 2. Content

Component	Cap Color	Format		Volume
Spiking Solution	White	Lyophilisate	tetracycline hydrochloride	1 µg
Reconstitution Buffer	Blue	Ready to use		11 ml

## 3. Warnings and precautions for the user

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 4. Storage instructions

Store the reagent at 2 - 8 °C (39.2 - 46.4 °F).

No quality guarantee is accepted after expiry of the component (see label).

## 5. Procedure

- equilibrate lyophilisate and reconstitution buffer to room temperature (20 - 25 °C / 68 - 77 °C) before use
- add 10 ml of reconstitution buffer to the lyophilisate, relock vial and vortex (concentration: 100 ng/ml)
- after reconstitution, the spiking solution can be stored at 2 - 8 °C (39.2 - 46.4 °F) for up to 6 months, but no longer than the shelf life (see label)
- avoid repeated warming and cooling; it is recommended to divide into aliquots directly after reconstitution
- if necessary, dilute the spiking solution with sample buffer (see instructions for use of RIDASCREEN® Tetracyclin (R3505) under 4. Reagents provided)
- spiking controls have to be prepared in the same media as used for the respective sample in the test
- spiking of samples should be conducted in such a way, that the concentration measured in the test is in the center of the standard curve; for calculation of the spike level, the dilution factor has to be taken into account
- the volume of spiking solution needed for the spiking of a certain sample amount can be calculated as follows:

$$\text{Volume of spiking solution} = \frac{\text{spike level} \times \text{sample quantity}}{\text{concentration of spiking solution}}$$

General example:

$$\text{Volume of spiking solution} = \frac{50 \frac{\text{ng}}{\text{g}} \times 1 \text{ g}}{1000 \frac{\text{ng}}{\text{ml}}} = 0.05 \text{ ml to 1 g sample}$$

- ideally, the ratio of spiking solution to sample amount should be between 1:100 and 1:10; if the ratio between spiking solution and sample amount is lower, the spiking solution has to be pre-diluted accordingly

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG
**R-Biopharm AG**