RIDA[®] QUICK Soya

Art. No. R7103

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RIDA[®] QUICK Soya

Brief information

RIDA[®] QUICK Soya (Art. No. R7103) is an immunochromatographic test for the qualitative detection of soya contamination (native and processed) on surfaces (swab test for the hygiene control in production and in laboratories) and in food (biscuit, margarine, chocolate, sausage).

All reagents required for the swab test are contained in the test kit.

The test kit contains 25 test strips (in a tube) for 1 determination each. Results are evaluated visually.

Time requirement: sampling for swab testapprox. 1 min
 preparation swab test.....5 min
 preparation sample 10 min
 test implementation (incubation time) 10 min

Detection limit: - on **surfaces** approx. 0.5 µg soya protein / 100 cm²
 - **soya flour in wheat flour** 0.5 mg/kg soya protein
 - in **processed food** approx. 10 mg/kg soya protein
 (depending on the matrix)

Specificity: The antibodies specifically detect heated soya proteins.
 (The sample extraction at **100 °C/212°F** ensures that
 all soya proteins denature and can be detected in the
 Lateral Flow Device).

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

Related products:

RIDA[®]QUICK Soya Extraction buffer (Art. No. R7113)

RIDA[®]QUICK Soya accessory pack (Art. No. Z7103)

RIDASCREEN[®]FAST Soya (Art. No. R7102)

Set of 3 processed Soya Assay Controls (Art. No. R7132)

SureFood[®] ALLERGEN Soya

1. Intended use

RIDA[®]QUICK Soya can be used as a swab test for the soya determination on surfaces in the hygiene control and for the qualitative detection of soya in raw material and processed food. The test has been developed for the detection of low amounts of soya (contamination).

2. General

The allergen can be present as an ingredient or as a contamination in raw and cooked products. According to the **regulation (EU) No. 1169/2011**, soya and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

Soybeans (mature seeds, raw) contain approx. 40 % proteins. Thus, soya is widely used as substitute for proteins from animal sources (e.g. in tofu, soya milk or soya yoghurt). Besides the high protein content, soya beans are also rich in fat (approx. 20 %) and are used for the productions of oils and fats.

The usage of soya as food and feed has strongly increased in the last decades, leading to an increase in allergic reactions as well.

3. Test principle

The basis of the immunochromatographic test is based on antibodies which are specific for the detection of soya proteins. Results are read visually. Generally, the higher the analyte level in the sample, the stronger the red color of the test band will be.

4. Reagents provided

Each kit contains sufficient reagent for 25 determinations.

Component	Cap color	Format	Volume
Dip stick	white	Ready-to-use	25 pieces
Extraction buffer	transparent	Ready-to-use	20 ml
Conjugate 1	yellow	Ready-to-use	4 ml
Conjugate 2	blue	Ready-to-use	4 ml
Reagent vial 2ml		Ready-to-use	25 pieces
Test tubes 5 ml		Ready-to-use	30 pieces
Swabs		Ready-to-use	26 pieces
Tip for Micropipete		Ready-to-use	50 pieces
Evaluation card		Ready-to-use	1 pieces

5. Materials required but not provided

- water bath, pot, electric kettle or similar
- vortexer
- pipettes
- RIDA[®]QUICK Soya accessory pack (Art. No. Z7103), contains 150 µl Pipette, Test tube holder and a Floating rack

For the analysis of food

- RIDA[®]QUICK Soya Extraction buffer (Art. No. R7113)

6. Warnings and precautions for users

Dust and dirty laboratory equipment lead to soya contamination of the assay. In order to avoid cross-contamination by cereal dust, please note the following points:

- wear gloves before starting and during the assay
- clean surfaces, glass vials, mincers and other equipment

The dip sticks are very sensitive to humidity that could turn the test useless. For this reason keep the strips away from humidity!

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.de.

7. Storage instructions

Store the unopened kit at 2 - 8 °C (36 - 46 °F). Do not freeze the kit.

No quality guarantee is accepted after the expiry date on the kit label.

8. Test implementation

8.1 Swab test

General: Use fresh pipette tips for each pipetting step

1. Take as many 5 ml Test tubes and 2 ml Reagent vials as samples to be analyzed.
2. Add 750 µl of Extraction buffer to 2 ml Tube.
3. Immerse the lower end of the swab into the Extraction buffer and swab a sampling area of 10 x 10 cm (wear gloves) thoroughly.
4. Insert the Swab into the 2 ml Reagent vial and cut/break the stick off.
5. Close the 2 ml Reagent vial and place it into a Floating rack (Art. No. Z7103).
6. Cook the 2 ml Reagent vial and the Floating rack for 5 min at 100°C/ 212°F. If no Floating rack is used, Tube may open and falsify the result.
7. Let the sample cool down for approx. 5 minutes.
8. Add 5 drops of Conjugate 1 (yellow cap) and Conjugate 2 (blue cap) into a new 5 ml Test tube.
9. Add 150 µl cooled sample, mix well and incubate for 5 min at room temperature.
10. Use gloves to place the Dip stick vertically into the 5 ml Test tube with the arrow end pointing down. Do not immerse the Dip stick beyond the maximum line.
11. Take out the Dip stick after exactly 5 min and read the result using the evaluation card.

8.2 Extraction of food using RIDA[®]QUICK Soya Extraction buffer

- Grind 5 g sample to powder
- Weigh in **1 g** of the sample into a resealable 15-50 ml reaction tube
- Add 7.5 ml** RIDA[®]QUICK Extraction buffer (Art. No. R7113)
- Mix the sample well by vortexing until the sample is homogenous.**
- Cook the solution for **10 min at 100 °C/212 °F** in a boiling water bath
- Cool down the sample (e.g. in an ice bath) and filter the extract

Then the test procedure with the dip sticks is carried out:

1. Add 5 drops of Conjugate 1 (yellow cap) and Conjugate 2 (blue cap) into a 5 ml Test tube.
2. Add 150 µl cooled sample, mix well and incubate for 5 min at room temperature.
3. Use gloves to place the Dip stick vertically into the 5 ml Test tube with the arrow end pointing down. Do not immerse the Dip stick beyond the maximum line.
4. Take out the Dip stick after exactly 5 min and read the result using the evaluation card.

Remark

The filtrate can be stored in a well closed vessel 1 day in the dark at 4 °C (39 °F).

9. Results and Sensitivity

Positive result: two colored bands

The sample is positive if two colored bands (the control band at the top and the lower test band) are visible within the result window.

Swab test: > 0.5 µg soya protein / 100 cm²

Food test: > 10 mg/kg soya protein

Negative result: only the upper control band

The sample is negative if only the upper control band is visible within the result window.

Swab test: < 0.5 µg soya protein / 100 cm²

Food test: < 10 mg/kg soya protein

Invalid result: no colored band

If no band is visible within the result window after performing the test, the test is considered invalid.

Limitations of the method

- The test strip has been developed for the detection of soya contamination.
- The limit of detection is dependent on sample type and extraction efficiency or the properties of the swabbed surface and the kind of contamination respectively.
- A negative result does not necessarily indicate the absence of soya as soya may be not homogenously distributed or the level of soya in the product is below the limit of detection.
- For high concentrations (approx. 12000 mg/kg soya protein) a high dose effect can be observed (Hook-effect). Therefore, a dilution series in these concentration is recommended.

In general

- Due to the multitude of food types, matrix effects cannot be excluded.
- In processed food, proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivities.

Recommendations

- For documentation, the upper part of the dip stick together with the test bands must be cut off.
- If the negative assay control sample is evaluated as positive then a contamination of the laboratory or laboratory equipment is likely.
- The ELISA RIDASCREEN[®]FAST Soya (Art. No. R7102) should be used for quantification.
- For confirmation of the result SureFood[®] PCR should be performed
- For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result, all cross reactivities and exemplarily analyzed matrices are described in the validation report
- Use also allergen-free and allergen containing (spiked) samples as test controls
- As a test control we recommend to use the Set of 3 processed Soya Assay Controls (R7132) (for the application note please contact sales@r-biopharm.de)

For further product information and application notes, please contact sales@r-biopharm.de.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG
R-Biopharm AG