RIDASCREEN[®] β -Agonists

Art. No. R1704

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RIDASCREEN® β -Agonists

Brief information

RIDASCREEN® β -Agonists (Art. No.: R1704) is a competitive enzyme immunoassay for the quantitative analysis of β -agonists in urine, serum, meat, liver, milk and feed.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: urine, serum (direct): centrifugation, dilution
 urine (SPE): hydrolysis, solid phase extraction (SPE),
 evaporation, reconstitution
 meat: homogenization, centrifugation, solid phase
 extraction (SPE), evaporation, reconstitution
 liver: homogenization, extraction, centrifugation,
 evaporation, reconstitution
 milk: extraction, centrifugation, evaporation,
 reconstitution
 feed: homogenization, centrifugation, extraction

Time requirement: sample preparation (for 10 samples)
 urine, serum (direct)approx. 30 min
 urine (SPE)approx. 4 h
 meatapprox. 2 h
 liverapprox. 1 h
 milkapprox. 45 min
 feedapprox. 1.5 h
 test implementation (incubation time) 1 h

Detection limit: urine (direct) approx. 200 ng/L (ppt)
(corresponding to the urine (SPE) approx. 150 ng/L (ppt)
standard substance) serum approx. 900 ng/L (ppt)
 meat approx. 100 ng/kg (ppt)
 liver approx. 130 ng/kg (ppt)
 milk approx. 45 ng/L (ppt)
 feed approx. 1000 ng/kg (ppt)

Recovery rate:	urine (direct)	approx. 94 %
(corresponding to the	urine (SPE)	approx. 113 %
standard substance)	serum.....	approx. 102 %
	meat.....	approx. 76 %
	liver	approx. 89 %
	milk	approx. 80 %
	feed.....	approx. 93 %

The specificity of the RIDASCREEN® β -Agonists test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the limit of detection and the recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Specificity:	Clenbuterol (standard substance)	approx. 100 %
	Salbutamol	approx. 90 %
	Cimbuterol	approx. 90 %
	Brombuterol	approx. 80 %
	Mabuterol	approx. 60 %
	Terbutaline.....	approx. 45 %
	Carbuterol	approx. 40 %
	Mapenterol	approx. 35 %
	Cimaterol	approx. 20 %

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

RIDA® β -Agonists Clenbuterol Spiking Solution	(R1799)
RIDASCREEN® Clenbuterol.....	(R1711)
Clenbuterol test control (positive)	(R1707)
Clenbuterol test control (negative).....	(R1708)
RIDASCREEN® Ractopamin	(R9901)
RIDA® Ractopamin Spiking Solution	(R9999)

1. Intended use

RIDASCREEN® β -Agonists is a competitive enzyme immunoassay for the quantitative analysis of β -Agonists in urine, serum, meat, liver, milk and feed.

2. General

β -Agonists, such as clenbuterol or salbutamol, are synthetic derivatives of the naturally occurring catecholamines.

It has been known that β -agonists are suitable for use as performance improvers within the field of livestock production. In particular, the meat/fat ratio in fattened animals can be improved or the growth may be accelerated. However, such compounds have not been approved in the EU for use as fattening adjuvants. In addition to lipolytic and anabolic effects, β -agonists have relaxing effects on non-striated musculature. Thus they can be used as antiasthmatic and tocolytic agents. It is possible that β -agonists residues, after use in illegal practice, may lead to a risk for consumers. Therefore the use of β -agonists use in food production was prohibited.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with antibodies directed against β -Agonists. β -Agonists standards or sample and β -Agonists conjugate are added. Free β -Agonists and β -Agonists conjugate compete for the β -Agonists antibody binding sites (competitive enzyme immunoassay). Any unbound conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the β -Agonists concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses):

Component	Cap Colour	Format		Volume
Microtiter plate	-	Ready to use		96 wells
Standard 1	White	Ready to use	0 ng/L	2 ml
Standard 2	White	Ready to use	62,5 ng/L	1 ml
Standard 3	White	Ready to use	125 ng/L	1 ml
Standard 4	White	Ready to use	250 ng/L	1 ml
Standard 5	White	Ready to use	500 ng/L	1 ml
Standard 6	White	Ready to use	1000 ng/L	1 ml
Standard 7	White	Ready to use	2000 ng/L	1 ml
Wash buffer	White	Concentrate	20x	30 ml
Conjugate	Red	Concentrate	100x	0.15 ml
Conjugate/sample buffer	White	Concentrate	4x	20 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		12 ml
Stop solution	Yellow	Ready to use		15 ml

5. Materials required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- mixer
- upside down shaker
- vortex
- evaporator
- column rack stand
- centrifuge
- graduated pipettes
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes

5.2. Reagents:

Urine (SPE) and meat

- Urine only: β -Glucuronidase/Arylsulfatase from *Helix pomatia* (Merck-Chemicals, Art. No.: 104114); dilute 1:10 (1+9) in demineralized water
- 0.1 M HCl
- 1 M Tris Base (Sigma-Aldrich, Trizma base, Art. No.: T1503)
- 100 % Methanol p.a.
- 0.1 M acetic acid
- 2 % ammonium hydroxide: add 6.25 ml of 32 % ammonia solution (Merck-Chemicals, Art. No.: 105426) to 93.75 ml of 100 % methanol p.a.
- RIDA[®] C18 column (R-Biopharm, Art. No.: R2002)

Liver and milk

- Acetonitrile

Feed

- 1 M NaOH
- 1 M HCl

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the substrate/chromogen prior to test implementation
- a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for the standard 1

9. Preparation of Samples

The samples should be stored in a cool place.

9.1. Urine (direct)

- centrifuge: 5 min / 3000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- dilute supernatant 1:5 (1+4) in sample buffer in a glass vial
e.g. 50 µl supernatant + 200 µl sample buffer and mix thoroughly
- use 25 µl per well in the assay

9.2. Urine (SPE)

If the urine is turbid or contains a precipitate, a centrifugation (5 min / 3000 g / room temperature (20 - 25 °C / 68 - 77 °F)) or filtration step (Filter paper (e.g. Macherey-Nagel, Art. No.: MN 615 A)) is mandatory.

- transfer 2 ml of urine into a clean glass vial
- add 20 µl of 1:10 diluted *Helix pomatia* β-Glucuronidase/ Arylsulfatase
- add 7 ml 0.1 M HCl and vortex
- incubate for 2 hours at 50 °C (122 °F)
- add 1 ml 1 M Tris Base and vortex
- continue with solid phase extraction under 9.4.

9.3. Meat

- homogenize completely
- transfer 2 g of homogenized meat into a clean glass vial
- add 7 ml 0.1 M HCl and vortex
- add 1 ml 1 M Tris Base and vortex
- centrifuge: 10 min / 2000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- transfer supernatant without disturbance of the pellet into a new glass vial
- continue with solid phase extraction under 9.4.

9.4. Solid phase extraction (SPE)

The supernatant is purified by means of RIDA[®] C18 column (Art. No. R2002):

- rinse the column with 1.2 ml 100 % methanol p.a.
- apply 1.2 ml 0.1 M acetic acid to the column
- apply 3 ml of the filtrate to the column (flow 1 ml/min)
- elute sample with 1 ml 2 % ammonium hydroxide solution in methanol
- apply positive pressure or vacuum for 30 s and dry the column for 2 min by floating it with air or nitrogen to collect all liquid
- evaporate the eluted sample to dryness under a mild nitrogen or air flow at 50 °C (122 °F)
- dissolve the dried residue in 0.6 ml sample buffer and vortex
- use 25 µl per well in the assay

9.5. Liver

- homogenize liver completely
- transfer 2 g of homogenized liver in a clean glass vial
- add 8 ml of acetonitrile and mix vigorously (vortex)
- mix for 10 min (upside down shaker)
- centrifuge: 10 min / 2000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- transfer 1 ml of the supernatant into a new glass vial
- evaporate to dryness under a mild nitrogen or air flow at 50 °C (122 °F)
- dissolve the dried residue in 400 µl sample buffer and vortex
- use 25 µl per well in the assay

9.6. Serum

- centrifuge: 5 min / 3000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- dilute supernatant 1:5 (1+4) in sample buffer and mix thoroughly
- use 25 µl per well in the assay

9.7. Milk

- transfer 2 ml of milk into a glass vial
- add 8 ml of acetonitrile and mix vigorously (vortex)
- mix for 10 min (upside down shaker)
- centrifuge: 10 min / 2000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- transfer 2 ml of the supernatant into a glass vial
- evaporate to dryness under a mild nitrogen or air flow at 50 °C (122 °F)
- dissolve the dried residue in 250 µl sample buffer and vortex
- use 25 µl per well in the assay

9.8. Feed

- grind sample (approx. 50 - 100 g) into a fine powder
- add 2 ml of 1 M HCl and 18 ml of demineralized water to 2 g of the ground feed
- vortex for 3 min
- shake for 15 min at room temperature (20 - 25 °C / 68 - 77 °F)
- centrifuge: 20 min / 2000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- transfer supernatant into a new vial and add 1 ml of 1 M NaOH
- adjust pH to 7.0 - 7.8 with 1 M HCl and vortex for 1 min
- centrifuge: 20 min / 2000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- use 25 µl per well in the assay

10. Test implementation

10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

Eventually occurring precipitates in the concentrates have to be dissolved by shaking at room temperature before dilution.

Please use the contained buffer concentrate as **conjugate and sample buffer**. Dilute buffer concentrate 1:4 (1+3) with demineralized water (e.g. 20 ml concentrate + 60 ml demineralized water). Diluted buffer can be stored at 2 - 8 °C (35 - 46 °F) until the shelflife of the concentrate expires.

Dilute **wash buffer** concentrate 1:20 (1+19) with demineralized water (e.g. 2 ml concentrate + 38 ml demineralized water, sufficient for 4 microtiter strips). Diluted buffer can be stored at 2 - 8 °C (35 - 46 °F) until the shelf life of the concentrate expires.

The **β-agonists conjugate** (bottle with red cap) is provided as a concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is needed should be reconstituted. Centrifuge the conjugate in the vial shortly (1 min / 1000 g / room temperature (20 - 25 °C / 68 - 77 °F)). For reconstitution, the concentrate is diluted 1:100 (1+99) in ready-to-use conjugate buffer (e.g. 30 μl conjugate concentrate + 2970 μl ready-to-use conjugate buffer, sufficient for 4 microtiter strips).

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 25 µl of standard or prepared sample to separate duplicate wells and add 75 µl of diluted conjugate of each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all wells with 300 µl of wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
4. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
5. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 min after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win / RIDA[®]SOFT Win.net (Art. No. Z9996), is available to evaluate the RIDASCREEN[®] enzyme immunoassays.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = B/B_0 (\%)$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the β-agonists concentration [ng/L].

In order to obtain the β -agonists concentration in ng/kg (ppt) actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factor is as follows:

urine (direct)	5
urine (SPE)	1
serum.....	5
meat.....	1
liver	2
milk	0.625
feed.....	11

For further information or applications please contact your local distributor or R-Biopharm AG.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG