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Advanced Settings

Integration time (s) 0.5
Acquisitions per Cycle 10

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
QuickGEN PCR Kit

Alicyclobacillus Differentiation

1. Intended use

Qualitative detection and differentiation of *Alicyclobacillus spp.*, *Alicyclobacillus acidoterrestris* and *Alicyclobacillus acidocaldarius*.

2. Test principle

The (TaqMan[®]) real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ; ATTO/DQ; HEX/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

To avoid false negative PCR-results an Inhibition control (ROX/DQ) is amplified together in one reaction vessel with the specific sequence. The system contains dUTP. Optional: the use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit).

3. Kit contents

The kit contains sufficient reagents for 48 PCR reactions:

1 x Premix	white cap
48 x Dye Strips (freeze-dried, incl. IC-DNA)	reaction tubes
1 x ddH ₂ O	colourless cap
1 x Control-DNA (freeze-dried)	yellow cap

4. Storage

The Control-DNA is freeze-dried and has to be solved in ddH₂O prior to use (see 6.1).

Do **not** freeze the lyophilized Dye Strips and lyophilized Control-DNA.

The PCR-reagents should be stored at 2 – 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) after opening. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

5. Materials required, but not provided

5.1. Instruments

Real-time PCR instrument for low profile tubes clear

Centrifuge for 1.5 – 2.0 mL reaction tubes

Centrifuge for strips

Pipettes

“Vortex”

5.2. Reagents and plastic ware

sterile ddH₂O

sterile filter tips

optional: Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction)

6. PCR

6.1. PCR-Setup

When using the kit for the first time, the freeze-dried Control-DNA has to be shortly centrifuged and carefully resolved:

- add 55 μL sterile ddH₂O to the freeze-dried Control-DNA
- after 15 minutes mix well

Before every use thoroughly mix all PCR components and centrifuge briefly.

PCR-reaction per sample:

PCR-Components	volume (μL)
Premix	15.0
Sample-DNA	5.0*
Total volume	20.0

* if using the Simplex easy DNA kit: add 2.5 μL DNA and 2.5 μL sterile ddH₂O

1. Pipette 15 μL of the Premix into each PCR-tube, making sure that, prior to the first filling, the tip has been moistened.
2. Add 5 μL sample DNA, add 5 μL of the Control-DNA for the positive control reaction, add 5 μL of the extraction control and 5 μL of ddH₂O for the negative control* reaction. Use a fresh tip with each DNA filling.
3. Close the tubes immediately and centrifuge them shortly.
4. Place the tubes in the PCR-machine and start run.

Very important: * Please fill up the negative control with 5 μL ddH₂O to avoid unspecific amplification.

Work swiftly to avoid warming up and keep away from light.

6.2 PCR-Program

Note:

For the use of UNG the thermal cycler program has to be changed according to manufacturers` instructions.

Software 3.3.2

Step	Time	Temp.	
Initial denaturation	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	25 sec	60 °C	

Advanced Settings

Integration time (s)	0.5
Acquisitions per Cycle	10

Software 3.4.8

Step	Time	Temp.	
Initial denaturation	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	35 sec	60 °C	

Advanced Settings

Integration time (s)	0.5
----------------------	-----

7. Evaluation

The evaluation has to be made according to the data analysis programme recommended by the real-time instrument manufacturer.

***Alicyclobacillus spp.*-DNA:** HEX-channel

***Alicyclobacillus acidoterrestris*-DNA:** FAM-channel

***Alicyclobacillus acidocaldarius*-DNA:** ATTO-channel

Inhibition Control-DNA: ROX-channel


A sample is ***Alicyclobacillus spp.* positive**, if there is a detectable fluorescence increase in the **HEX-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the ROX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is ***Alicyclobacillus acidoterrestris* positive**, if there is a detectable fluorescence increase in the **FAM-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the ROX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is ***Alicyclobacillus acidocaldarius* positive**, if there is a detectable fluorescence increase in the **ATTO-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the ROX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is **negative**, if there is no detectable fluorescence increase in the different channels and the positive controls have a positive fluorescence signal. The negative controls show no amplification. The Inhibition Control in the ROX-channel has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

analysis flowchart

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method. GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.