RIDASCREEN[®] FAST Citrinin

Art. No. R6302

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RIDASCREEN®FAST Citrinin

Brief information

RIDASCREEN®FAST Citrinin (Art. No. R6302) is a competitive enzyme immunoassay for the quantitative determination of citrinin in cereals and feed.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: extraction, filtration and dilution

Time requirement: sample preparation (for 10 samples)..... approx. 10 min
test implementation (incubation time).....25 min

Limit of detection: 15 µg/kg (ppb)
(corresponding to the standard substance)

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

TRILOGY® Dried Standard Citrinin (TS-904)

EASI-EXTRACT® CITRININ (R6302)

1. Intended use

RIDASCREEN®FAST Citrinin assay is a competitive enzyme immunoassay for the quantitative determination of citrinin in cereals and feed.

2. General

The mycotoxin Citrinin is formed by fungi of the species *Aspergillus*, *Penicillium*, and *Monascus* species. These fungi are able to produce Citrinin and/or Ochratoxin A, therefore both mycotoxins appear together very often.

Citrinin was described for the first time in 1953 by Japanese authors. They detected this mycotoxin in yellow coloured rice which was contaminated with *Penicillium citrinum*.

Citrinin occurs in cereals such as wheat, barley, rye, oats, corn and rice and shows the same nephrotoxic effects as ochratoxin A.

The moisture level of cereals is very important for the development and growth of the citrinin producing fungi, requiring minimum levels of 16.5 to 19.5 %.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with citrinin. Standards, respectively sample solution and anti-citrinin antibodies, are added. Free and immobilized citrinin compete for the citrinin antibody binding sites (competitive enzyme immunoassay). After a washing step, secondary antibodies labeled with peroxidase are added. These bind to the bound anti-citrinin antibodies. Any unbound enzyme conjugated secondary antibodies are then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate (secondary antibodies labeled with peroxidase) converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is performed photometrically at 450 nm. The absorbance is inversely proportional to the citrinin concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for as many as 43 analyses (plus 5 standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
Standard 1*	white	Ready to use	0 µg/L	1.3 mL
Standard 2*	white	Ready to use	15 µg/L	1.3 mL
Standard 3*	white	Ready to use	45 µg/L	1.3 mL
Standard 4*	white	Ready to use	135 µg/L	1.3 mL
Standard 5*	white	Ready to use	405 µg/L	1.3 mL
Conjugate	red	Ready to use		6 mL
Antibody	black	Ready to use		3 mL
Substrate/Chromogen Red Chromogen	brown	Ready to use		10 mL
Stop Solution	yellow	Ready to use		14 mL

*) The dilution factor 5 for the sample has been already considered. Therefore, the citrinin concentrations of samples can be read directly from the standard curve.

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- graduated cylinder (plastic or glass) 100 ml
- for preparing sample extract: filter funnel and 50 ml flask
- grinder (mill)
- optional: shaker
- filter paper: Whatman No. 1 or equivalent
- 50 µl, 100 µl and 1000 µl micropipettes

5.2. Reagents:

- 70 % methanol solution: mix 70 ml methanol (100 %) with 30 ml distilled or deionized water
- distilled or deionized water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The standards contain citrinin, particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

Decontamination of the glassware and citrinin solutions is best carried out using a sodium hypochlorite solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (36 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

The test kit can be regularly used at least up to the expiry date (indicated on the kit package), if stored correctly.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the reddish substrate/chromogen prior to test implementation
- a value of less than 0.6 absorbance units ($A_{450\text{ nm}} < 0.6$) for the zero standard

9. Preparation of Samples

The samples should be stored in a cool place, protected from light.

Bring all reagents and samples to room temperature (20 - 25 °C / 68 - 77 °F) before use.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- weigh 5 g of ground sample and add it to a suitable container with 12.5 ml of methanol (70 %) *)
- shake vigorously for 3 min (manually or with shaker)
- filter the extract through Whatman No. 1 filter
- dilute 1 ml of the obtained filtrate with 1 ml of distilled or deionized water
- use 50 µl of the filtrate per well in the test

*) sample size may be increased if required, but the volume of methanol/water must be adapted accordingly

10. Test implementation

10.1. Preliminary comments

1. Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.
2. The specific reaction starts with the addition of the specific antibody. Nevertheless, not more than three strips should be applied to the test, when a single step pipette is utilized. More strips (up to 6) can be applied, when a multistep pipette is used.
3. Return all reagents to 2 - 8 °C (36 - 46 °F) immediately after use.

The **citrinin standards** are provided ready to use. The dilution factor 5 for the sample has been considered when labeling. Therefore, the citrinin concentration of samples can be read directly from the standard curve.

10.2. Test procedure

Accurate washing is very important. Do not allow microwells to dry up totally and avoid prolonged intervals between the working steps. Reproducibility in any EIA is largely dependent upon the consistency with which the microwells are washed. Carefully follow the recommended washing sequence as outlined in the EIA test procedure.

Avoid direct sunlight during all incubations. Covering the microtiter plates is recommended.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Pipet 50 μ l of standard or prepared sample into separate wells; use a new pipette tip for each standard or sample.
3. Add 50 μ l of the antibody to each well. Mix gently by shaking the plate manually and incubate for 10 min (+/- 1) at room temperature (20 - 25 °C / 68 - 77 °F).
4. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Using a wash bottle or multichannel pipette, fill the wells with distilled or deionized water (250 μ l per well). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.
5. Add 100 μ l of the conjugate (secondary antibody) to the bottom of each well and incubate for 10 min (+/- 1) at room temperature (20 - 25 °C / 68 - 77 °F).
6. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Using a wash bottle or multichannel pipette, fill the wells with distilled or deionized water (250 μ l per well). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.
7. Add 100 μ l of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 5 min (+/- 0.5) at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
8. Add 100 μ l of stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win.net (Art. No. Z9996), is available to evaluate the RIDASCREEN[®] enzyme immunoassays.

For single determinations we recommend logit/log evaluation and for double or multiple determinations cubic spline should be used.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the citrinin concentration [$\mu\text{g}/\text{kg}$]. The citrinin concentration in $\mu\text{g}/\text{kg}$ corresponding to the extinction of each sample can be read from the calibration curve.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG