RIDASCREEN[®] Nitrofurantoin (AOZ)

Art. No. R3703

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RIDASCREEN® Nitrofurantoin (AOZ)

Brief information

RIDASCREEN® Nitrofurantoin (AOZ) (Art. No.: R3703) is a competitive enzyme immunoassay for the quantitative analysis of AOZ in shrimp, meat (chicken, pork, beef), liver (beef and pork), fish, whole egg and milk.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit. The test kit is sufficient for 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization (milk samples: precipitation), derivatization, extraction, centrifugation, evaporation and degreasing

Time requirement: sample preparation for 10 samples and test procedure

option I (short):

sample preparation part 1 approx. 1.5 h
incubation 3 h
sample preparation part 2 approx. 1.5 h
test implementation (incubation time) 1 h 15 min

option II (long):

sample preparation part 1 approx. 1.5 h
incubation (overnight) approx. 16 h
sample preparation part 2 approx. 1.5 h
test implementation (incubation time) 1 h 15 min

Detection limit: shrimp, fish, milk approx. 50 ng/kg (ppt)
(corresponding to the standard substance) meat, liver, whole egg approx. 100 ng/kg (ppt)

Recovery rate: in spiked samples
(corresponding to the standard substance) whole egg approx. 90 - 100 %
shrimp, fish approx. 85 - 100 %
meat (chicken, beef), liver, milk approx. 80 - 100 %
meat (pork) approx. 75 - 100 %

The specificity of the RIDASCREEN® Nitrofurantoin (AOZ) test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Specificity:	NP-AOZ (standard substance).....	100 %
	NP-AHD, NP-AMOZ.....	< 0,01 %
	Nitrofurantoin	< 0,01 %
	Furazolidone	1,2 %
	AMOZ, AHD, SEM	< 0,01 %

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. It lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

RIDASCREEN® Nitrofurantoin (AMOZ)	(R3711)
RIDASCREEN® Nitrofurantoin (AHD).....	(R3713)
RIDASCREEN® Nitrofurantoin (SEM)	(R3715)
RIDA® Nitrofurantoin (AOZ) Spiking solution.....	(R3798)
RIDA® Nitrofurantoin (AMOZ) Spiking solution.....	(R3799)
RIDA® Nitrofurantoin (AHD) Spiking solution.....	(R3796)
RIDA® Nitrofurantoin (SEM) Spiking solution	(R3797)

1. Intended use

RIDASCREEN[®] Nitrofurantoin (AOZ) is a competitive enzyme immunoassay for the quantitative analysis of AOZ in shrimp, meat (chicken, pork and beef), liver (beef and pork), fish, whole egg and milk.

2. General

Nitrofurans are synthetic broad-spectrum antibiotics, which are frequently used in animal production due to their excellent antibacterial and pharmacokinetic properties. They have also been used as growth promoters during the production of shrimp, poultry and pigs. Long term animal experiments have shown that the parent compounds and their metabolites have carcinogenic and mutagenic characteristics. This led to the prohibition of nitrofurans for the treatment of animals used for food production. In 1993, the EU banned the nitrofurans furaltadone, nitrofurantoin and nitrofurazone for use in animals used as sources of food, and in 1995 the use of furazolidone was also prohibited. The analysis of nitrofurans is based on the detection of the tissue bound metabolites of nitrofurans. The parent compounds are difficult to detect accurately since they are metabolized very rapidly after treatment. The tissue bound nitrofurantoin metabolites however are present for a long time after administration and they are used to detect nitrofurantoin abuse.

Parent compound	Metabolite	Abbreviation
Nitrofurantoin	1-Aminohydantoin	AHD
Furaltadone	3-Amino-5-morpholinomethyl-2-oxazolidinone	AMTZ
Furazolidone	3-Amino-2-oxazolidinone	AOZ
Nitrofurazone	Semicarbazide	SEM

Prior to analysis, the metabolites have to be derivatized by incubation with 2-Nitrobenzaldehyde into NP-AHD, NP-AMTZ, NP-AOZ and NP-SEM.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-AOZ antibodies. AOZ standards or sample solution, AOZ conjugate and anti-AOZ antibodies are added. Free AOZ and AOZ conjugate compete for the AOZ antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-AOZ antibodies are also bound by the immobilized capture antibodies. Any unbound conjugate is then removed in a washing step.

Substrate/chromogen is added to the wells and incubated. Bound conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm; the absorption is inversely proportional to the AOZ concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses).

Component	Cap color	Format	Volume
Microtiter plate	-	Ready to use	96 wells
2-Nitrobenzaldehyde		solid	100 mg
Standard 1	White	Ready to use	0 ng/L 1.3 ml
Standard 2	White	Ready to use	25 ng/L 1.3 ml
Standard 3	White	Ready to use	50 ng/L 1.3 ml
Standard 4	White	Ready to use	100 ng/L 1.3 ml
Standard 5	White	Ready to use	200 ng/L 1.3 ml
Standard 6	White	Ready to use	400 ng/L 1.3 ml
Wash buffer salt Tween		For dissolving	
Conjugate	Red	Ready to use	6 ml
Antibody	Black	Ready to use	6 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use	10 ml
Stop solution	Yellow	Ready to use	14 ml

5. Materials required but not provided

5.1. Equipment:

Equipment	Shrimp	Meat, Fish, Liver, whole egg	Milk
Microtiter plate spectrophotometer (450 nm)	•	•	•
Centrifuge	•	•	•
Shaker (Vortex)	•	•	•
Mixer (ultraturrax)	•	•	•
Incubator	•	•	•
Evaporator	•	•	•
Magnetic stirrer	•	•	•
Pasteur pipettes	•	•	•
Graduated pipettes: variable 20 - 200 µl and 200 - 1000 µl micropipettes	•	•	•

5.2. Reagents:

Reagent	Shrimp	Meat, Fish, Liver, whole egg	Milk
1 M HCl	•	•	•
1 M NaOH	•	•	•
0.1 M K ₂ HPO ₄	•	•	•
Ethyl acetate	•	•	•
n-hexane (oder n-Heptane)	•	•	•
Dimethylsulfoxide (DMSO)	•	•	•
Carrez I: 0,36 M Potassium ferrocyanide (II) x 3 H ₂ O			•
Carrez II: 1,04 M Zinc sulfate x 7 H ₂ O			•

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen prior to test implementation
- a value of less than 0.6 absorbance units ($A_{450\text{ nm}} < 0.6$) for standard 1

9. Preparation of Samples

The samples should be stored in a cool place, protected against light.

Prepare 10 mM 2-nitrobenzaldehyde (see 4.) in dimethylsulfoxide (DMSO) (this solution has to be prepared directly before use, i.e. dissolve 7.6 mg 2-nitrobenzaldehyde in 5 ml DMSO).

9.1. Shrimp, meat (chicken, pork and beef), liver (bovine and porcine), fish, and whole egg (freshly whisked egg)

- homogenize the sample, use stomacher or mixer

Shrimp:

- mix 1 g of the homogenized sample with 4 ml distilled water, 0.5 ml 1 M HCl and 100 µl 10 mM 2-nitrobenzaldehyde (in DMSO) by shaking vigorously.

Note: pH value at this point should be between 1 and 2, if necessary adjust it by adding more 1 M HCl

- continue as described in 9.3.

meat (beef), liver (bovine, porcine), fish and whole egg (fresh whisked egg):

- mix 1 g of the homogenized sample with 3.9 ml distilled water, 0.5 ml 1 M HCl and 200 µl 10 mM 2-nitrobenzaldehyde (in DMSO) by shaking vigorously

- continue as described in 9.3.

meat (chicken, pork):

- mix 1 g of the homogenized sample with 3.8 ml distilled water, 0.5 ml 1 M HCl and 300 µl 10 mM 2-nitrobenzaldehyde (in DMSO) by shaking vigorously
- continue as described in 9.3.

9.2. Milk

- transfer 5 ml of milk into a centrifugal glass vial
- add 250 µl Carrez I and 250 µl Carrez II (see 5.2.)
- mix thoroughly (vortex) and centrifuge: 10 min / 3000 g / 4 - 12 °C (39 - 54 °F). If a refrigerated centrifuge is not available, chill sample to approx. 8 °C (46 °F) prior to centrifugation.
- mix 1.1 ml of the supernatant with 3.8 ml distilled water, 0.5 ml 1 M HCl and 300 µl 10 mM 2-nitrobenzaldehyde (in DMSO) by shaking vigorously
- continue as described in 9.3.

9.3. Derivatisation

- incubate at 50 °C (122 °F) for 3 hours or at 37 °C (98.6 °F) overnight (approx. 16 h)
- add 5 ml 0.1 M K₂HPO₄, 0.4 ml 1 M NaOH and 5 ml ethyl acetate, shake vigorously for 30 sec.
- optional: sample incubation up to 50 °C (122 °F) can be used for better phase separation
- centrifuge: 10 min / 3000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- transfer 2.5 ml of the ethyl acetate layer (upper layer) into a new vial and evaporate to dryness at 60 °C
- dissolve the residue in 1 ml n-hexane (or n-heptane) and mix properly with 1 ml sample buffer (see 10.1.)
- centrifuge: 10 min / 3000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- use 50 µl of the lower, aqueous phase per well in the assay

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

A PBS-Tween buffer is needed as **sample and wash buffer**, please use the wash buffer salt (see 4.) contained in the kit. Dissolve the entire buffer salt in one liter of distilled water. The ready to use buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternative: Dissolve the contents of the pouch in only 100 ml of distilled water to obtain a 10fold concentrated buffer. This solution expires after approx. 8 - 12 weeks, store at room temperature (20 - 25 °C / 68 - 77 °F). Use 1 part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use buffer.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µl of each standard or prepared sample to separate duplicate wells. Use a new pipette tip for each standard or sample.
3. Add 50 µl of the conjugate (red cap) to the bottom of each well.
4. Add 50 µl of the antibody (black cap) to each well, mix gently by shaking the plate manually and incubate for 1 h at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl of substrate/chromogen (brown cap) to each well. Mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution (yellow cap) to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 min after addition of stop solution.

11. Results

A special software, the RIDASOFT® Win.NET (Art. No. Z9996), is available to evaluate the RIDASCREEN® enzyme immunoassays.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = B/B_0 (\%)$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the AOZ concentration [ng/kg].

In order to obtain the AOZ concentration in ng/kg actually contained in a sample, the concentration read from the standard curve must be further multiplied by the corresponding dilution factor.

When working in accordance with the regulation stated, **the dilution factor is 2.**

For further information or applications please contact your local distributor or info@r-biopharm.de

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG