RIDASCREEN[®] FAST Sesame

Art. No. R7202

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RIDASCREEN[®] FAST Sesame

Brief information

RIDASCREEN[®] FAST Sesame (Art. No. R7202) is a sandwich enzyme immunoassay for the quantitative analysis of sesame or parts of sesame in instant soup, chocolate dessert, bakery products and crackers.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit. The test kit is sufficient for 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization, extraction and centrifugation

Time requirement: sample preparation (for 10 samples)approx. 20 min
test implementation (incubation time)30 min

Limit of detection: 0.14 mg/kg (ppm) Sesame; between 0.08 – 0.20 mg/kg (ppm) depending on the matrix

Limit of quantification: 2.5 mg/kg (ppm) Sesame

Specificity: The antibody specifically detects proteins from sesame.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments. Further information is described in the updated validation report.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded from the website <http://www.r-biopharm.com/products/food-feed-analysis>.

Related products

bioavid Lateral Flow Sesame (Art. No. BL609-10/-25)

SureFood® PCR ALLERGEN ID Sesame

1. Intended use

RIDASCREEN®FAST Sesame (Art. No. R7202) is a sandwich enzyme immunoassay for the quantitative analysis of sesame or parts of sesame (as ingredient components or contamination) in instant soup, chocolate dessert, bakery products and crackers.

2. General

Sesame seeds (*Sesamum indicum*) are used as raw or roasted ingredients in many different foods.

The allergen can be present as an ingredient or as a contamination in raw and cooked products. According to the regulation (EU) No. 1169/2011, Sesam and products thereof must be declared on food labels. Similar regulations exist e.g. in Canada and Australia.

Sesame seed food allergy is an increasingly recognized health burden, the symptoms vary from allergic eczema to life threatening anaphylactic reactions. Until now, seven different proteins of sesame were identified as allergens (Ses i 1-7). The only effective treatment to protect sensitized patients from allergic symptoms is a strict avoidance of sesame containing food.

3. Test principle

The wells of the microtiter strips are coated with specific antibodies to sesame-proteins. By adding standards and samples to the wells, sesame proteins present will bind to the specific antibodies. In a washing step components not bound are removed. Then antibody conjugated to peroxidase is added. This antibody conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. The detection of sesame protein takes place by adding substrate/chromogen. The enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the sesame concentration of the sample. The result is expressed in mg/kg sesame.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses).

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use	-	48 wells
Allergen extraction buffer	Green	Concentrate	10x	100 ml
Standard 1*	Transparent	Ready to use	0 mg/kg	1.3 ml
Standard 2*	Transparent	Ready to use	2.5 mg/kg	1.3 ml
Standard 3*	Transparent	Ready to use	5.0 mg/kg	1.3 ml
Standard 4*	Transparent	Ready to use	10.0 mg/kg	1.3 ml
Standard 5*	Transparent	Ready to use	20.0 mg/kg	1.3 ml
Wash buffer	Brown	Concentrate	10x	100 ml
Conjugate	Red	Ready to use	-	6 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use	-	10 ml
Stop solution	Yellow	Ready to use	-	14 ml

*) The dilution factor 20 for the sample has already been considered when labeling. Therefore, the sesame concentrations of samples can directly be read from the standard curve.

5. Materials required but not provided

5.1. Equipment

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal vials
- water bath
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenisator
- graduated pipettes
- paper filter
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes
- if possible multichannel pipette and/or multistepper pipette

5.2. Reagents

- distilled or deionized water
- skim milk powder (SMP, food quality)

6. Warnings and precautions

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiry date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen prior to test implementation
- a value of less than 1.2 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 5

9. Preparation of Samples

Working devices such as a mill, glass vials or spatulas must be cleaned before and after each sample preparation to remove any remains of sesame and to avoid contamination.

The **Allergen Extraction buffer** is provided as a **10fold concentrate**. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1 + 9) with distilled water before use (e.g. 100 ml buffer concentrate + 900 ml distilled water). The diluted Allergen Extraction buffer is stable at 20 – 25 °C (68 -77 °F) for approx. 4 weeks.

Add 5 g skim milk powder (SMP) per 100 ml diluted Allergen Extraction buffer. Prepare only the quantity of SMP-buffer that is actually needed. The SMP-buffer

should only be heated once to 60°C (140 °F) and it should be used within 24 hours.

9.1. Sample preparation

- grind 5 g of the sample carefully and mix thoroughly
- weigh 1 g of the sample and add 20 ml final diluted Allergen Extraction buffer containing SMP (the extraction buffer should already have a temperature of approx. 60 °C (140 °F))
- or add 19 ml of the final diluted Allergen Extraction buffer containing SMP to 1 ml liquid sample (the extraction buffer should already have a temperature of approx. 60 °C (140 °F))
- mix intensively and incubate for 10 min at 60 °C (140 °F), afterwards cool down (e.g. in an ice bath)
- centrifuge: 10 min, at least 2500 g, if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl of the supernatant or filtrate per well in the assay

Remark

Further dilutions should be prepared with the finally diluted Allergen Extraction buffer with SMP.

The sample extracts can be stored at 4 °C (39 °F) for 1 day.

10. Test implementation

10.1. Test Preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **wash buffer** is provided as a 10fold concentrate. Before dilution of the buffer concentrate, dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1 + 9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted wash buffer is stable at 20 - 25 °C (68 - 77 °F) for approx. four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than 3 strips (24 wells) at a time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101 or Mikrotiter Assembly Breakable Strip 1x8, Thermo Scientific) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standard and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard or prepared sample to separate duplicate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells.
4. Fill all the wells with 250 µl wash buffer (see 10.1). Repeat two more times.
5. Add 100 µl of conjugate to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
6. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
7. Add 100 µl substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
8. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDASOFT® Win.NET (Art. Nr. Z9996) is available for evaluation of the RIDASCREEN® enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit. The document "Compliance Criteria" provides criteria for evaluating standard curves. For quality assurance assay controls should be used.

In comparison with the certificate, higher values of the absorbance ($A_{450\text{ nm}}$) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or sesame contamination.

Please note

When working according to the described sample preparation, the dilution factor is 20. The allergen concentration can be read directly from the standard curve (the sample dilution factor of 20 is already taken into account for the standard concentrations (see 4. *)).

In general

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

Allergen containing samples that have been heat treated show a reduced recovery because the proteins denature and are no longer recognized by the antibody. The reduction in recovery depends strongly on the temperature and the duration of the heat treatment. If samples are heat treated at high temperature the recovery can be significantly reduced.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross-reactivities and exemplary analyzed matrices are described in the updated validation report.

The protein content and the protein composition may vary considerably between different sesame species. Therefore, different varieties may produce different results, since exemplary varieties were used for calibration. The protein content of the standard material is described in the validation report and can be used for the calculation to protein.

Recommendation:

In order to ensure a high analytical performance it is recommended to:

- Adjust the pH to a neutral value for extremely acidic or alkaline samples
- Use allergen-free and allergen containing (spiked) samples as test controls
- Carry out spiking experiments for an accurate and correct procedure
- Analyze each sample material in duplicates
- Perform SureFood[®] PCR to confirm results

– Contact sales@r-biopharm.de if automates (e.g. ThunderBolt[®] / Bolt) are used

Further product information and application notes, please contact sales@r-biopharm.de.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG