**CONGEN**

**SureFast<sup>®</sup> Animal+Plant  
Control 3plex**

Art. No. F4053

100 rxn

**User Manual**

**April 2018**

** Inhalt /  Content**

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## 1 General Information

### 1.1 Description

The SureFast® Animal+Plant Control 3plex real time PCR controls the DNA-preparation. In addition a differentiation between the DNA of vertebrates and plants is possible. The test contains an Internal Control DNA (ICD) as an internal control of sample preparation procedure and to determine possible PCR-inhibition.

The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of three fluorescence emissions at 522 nm, 553 nm and 670 nm (FAM, VIC/HEX and Cy5) at the same time. The technical validation of instruments was performed on Roche LightCycler® 480 II/Roche cobas z 480 Analyzer, Qiagen Rotor-Gene Q, Applied Biosystems 7500, Bio-Rad CFX96, BMS MIC, Agilent Mx3005P and Agilent AriaDx.

### 1.2 Limit of Detection

The SureFast® Animal+Plant Control 3plex real-time PCR has a limit of detection of  $\leq 500$  DNA copies. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

### 1.3 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 $\mu$ l	Yellow
2	Taq Polymerase	1 x 11 $\mu$ l	Red
D	Internal Control DNA	2 x 1800 $\mu$ l	Orange
3	Positive Control	1 x 200 $\mu$ l	Light Blue

Store all reagents at  $-20^{\circ}\text{C}$  and protected from light.

### 1.4 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFood® PREP Basic or SureFood® PREP Advanced)
- real-time PCR instrument with three detection channels (522 nm, 553 nm and 670 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

**1.5 Setup**

	<b>Blockcycler/MIC/LightCycler® 480 II</b>	<b>Rotor-Gene</b>
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	35	35
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup (exemplary)	Detection: End of Extension Phase  Detection System Plant: Various devices <b>FAM</b> -channel, Quencher: BHQ LightCycler 480 II 465 nm - 510 nm  Detection System Internal Control DNA (ICD): Various devices <b>VIC/HEX</b> -channel, Quencher: BHQ LightCycler 480 II 533 nm - 580 nm  Detection System Vertebrates: Various devices <b>Cy5</b> -channel, Quencher: BHQ LightCycler 480 II 618 nm - 660 nm	
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: <a href="http://www.congen.de/en/company/downloads">http://www.congen.de/en/company/downloads</a>		

**2 Qualitative Analysis**

**2.1 Protocol**

**2.1.1 DNA-Preparation**

For DNA-preparation the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended.

The test assay contains an Internal Control DNA (ICD), which can either be used as PCR inhibition control or as extraction control for the sample preparation procedure and as a PCR inhibition control.

If the ICD is used only as a PCR inhibition control, 1 µl of the ICD should be added to the Master-Mix.

If the ICD is used as an extraction control for the sample preparation procedure and as PCR inhibition control, 20 µl of the ICD should be added during extraction procedure. The ICD should always be added to the specimen-lysis buffer mixture and must not be added directly to the specimen.

**2.1.2 Preparation of the master-mix**

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and blank of extraction. The test assay contains an Internal Control DNA (ICD), which can either be used as PCR inhibition control or as extraction control.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

April 2018

Example for the calculation and preparation of 10 reactions for **ICD as extraction and PCR inhibition control**:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix each master-mix well and centrifuge shortly before use.

Example for the calculation and preparation of 10 reactions for **ICD only as PCR inhibition control**:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Internal Control RNA	1.0 µl	11.0 µl
Taq Polymerase	0.1 µl	1.1 µl
<b>Total volume</b>	<b>21 µl</b>	<b>231 µl</b>

Mix each master-mix well and centrifuge shortly before use.

### 2.1.3 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries.
- Pipette 5 µl of the Positive Control into the designated tubes/wells or capillaries.
- Centrifuge all tubes/wells or capillaries shortly at low speed.
- Place tubes/wells or capillaries into the real-time PCR instrument and start the run according to the setup.

### 2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

With this sensitive real-time PCR kit it is possible to detect minimal contaminations by the target DNA (for example airborne pollen, human DNA from the operator). These can lead to positive signals of the negative controls. Experience has shown that in this case the negative control has a Cp-value of > 33. The positive control should have a Cp-value around Cp 17-20.

Plant DNA is detected in the FAM-channel and vertebrates DNA is detected in the Cy5-channel. The internal amplification/extraction control will be detected in the VIC/HEX-channel.

A sample is stated **positive** for the respective parameter (plant/vertebrates), if the sample DNA shows amplification in the respective channel. A sample is stated **negative** for the respective parameter (plant/vertebrates), if the sample DNA shows no amplification in the respective channel and the ICD signal (VIC-channel) of the sample is **positive**.

If the sample DNA and the internal control DNA are **negative** the sample contains PCR-inhibiting substances or the sample preparation was not successful. Under these circumstances an evaluation of the samples is not possible. DNA isolation and purification for the sample need to be improved (see also table below).

April 2018

result in the respective channel			
FAM-channel plant	VIC/HEX-channel ICD	Cy5 - channel vertebrates	result
<b>positive</b>	<b>positive</b>	negative	<b>plant DNA detected</b>
negative	<b>positive</b>	<b>positive</b>	<b>vertebrates DNA detected</b>
negative	negative	negative	invalid

**Note:** The SureFast® Animal+Plant Control 3plex real-time PCR detects DNA of vertebrates. That's the reason for the detection of human DNA (*Homo sapiens*).

### 3 Further Information

#### 3.1 Product Information

- Validation Report

#### 3.2 Technical Support

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).