# RIDA<sup>®</sup> 17 $\beta$ -Östradiol

**Art. No. R2399**

Dotierlösung / Spiking Solution

In vitro Test

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# RIDA<sup>®</sup> 17β-Östradiol Spiking Solution

## 1. Intended use

RIDA<sup>®</sup> 17β-Östradiol Spiking Solution can be used for the production of positive control samples suitable for the validation (e.g. determination of the recovery rate) of the RIDASCREEN<sup>®</sup> 17β-Östradiol test (R2301).

## 2. Content

Component	Cap Colour	Format	Volume
Spiking Solution	White	Ready to use	200 ng 17β-estradiol / ml methanol 1 ml

## 3. Warnings and precautions for the user

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 4. Storage instructions

Store the reagent at - 20 °C (- 4 °F).

No quality guarantee is accepted after expiry of the component (see label).

## 5. Procedure

- equilibrate spiking solution to room temperature (20 - 25 °C / 68 - 77 °C) before use
- take an aliquot of the spiking solution, relock the bottle at once and continue storage at - 20 °C (- 4 °F)
- if necessary, dilute the spiking solution with methanol 100 % (v/v)
- to avoid any impact on sample preparation, the concentration of methanol in the spiked sample must not exceed 10 % (v/v)
- spiking controls have to be prepared in the same media as used for the respective sample in the test
- spiking of samples should be conducted in such a way, that the concentration measured in the test is in the center of the standard curve; for calculation of the spike level, the dilution factor has to be taken into account
- the volume of spiking solution needed for the spiking of a certain sample amount can be calculated as follows:

$$\text{Volume of spiking solution} = \frac{\text{spike level} \times \text{sample quantity}}{\text{concentration of spiking solution}}$$

General example:

$$\text{Volume of spiking solution} = \frac{50 \frac{\text{ng}}{\text{g}} \times 1 \text{ g}}{1000 \frac{\text{ng}}{\text{ml}}} = 0.05 \text{ ml to } 1 \text{ g sample}$$

- ideally, the ratio of spiking solution to sample amount should be between 1:100 and 1:10. If the ratio between spiking solution and sample amount is lower, the spiking solution has to be pre-diluted accordingly
- due to high vapour pressure, it is strongly recommended to pipette organic solvents using pipettes with air displacement technique or hamilton syringes; if single channel pipettes are used to transfer organic solvents, the pipette tip and the air inside the pipette needs to be saturated with organic solvent vapour before pipetting the desired volume; for this, aspirate and dispense the organic solvent at least 3 times before the desired volume is transferred

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG
**R-Biopharm AG**