# **RIDASCREEN<sup>®</sup> $\beta$ -Lactoglobulin**

**Art. No.: R4901**

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- Analysis of gluten or  $\beta$ -Lactoglobulin in supplemental enzymes

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# RIDASCREEN<sup>®</sup> $\beta$ -Lactoglobulin

## Brief information

RIDASCREEN<sup>®</sup>  $\beta$ -Lactoglobulin (Art. No. R4901) is a competitive enzyme immunoassay for the quantitative analysis of  $\beta$ -lactoglobulin in hydrolyzed milk products including hypoallergenic baby food.

For the detection of native and processed  $\beta$ -lactoglobulin in beverages containing whey, milk or milk powder as well as in juice, wine, beer and in food like sausages, dressings, bakery products, ice, chocolate, yoghurt, etc. the use of RIDASCREEN<sup>®</sup> FAST  $\beta$ -Lactoglobulin (Art. No. R4902) is recommended.

To detect casein or caseinates in food, the use of RIDASCREEN<sup>®</sup> FAST Casein (Art. No. R4612) is recommended.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for analysis.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples) .....approx. 30 min  
test implementation (incubation time) .....2 h 45 min

Limit of detection: 2.1 mg/kg  $\beta$ -lactoglobulin

Limit of quantification: 5 mg/kg  $\beta$ -lactoglobulin

Specificity: was determined by analyzing the cross reactivities to corresponding substances

$\beta$ -Lactoglobulin .....	100 %
$\alpha$ -, $\beta$ -, $\kappa$ -Casein.....	< 1 %
$\alpha$ -Lactalbumin.....	< 1 %
bovine serum albumin .....	< 1 %

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. corn bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

### **Related products**

RIDASCREEN® FAST Milk (Art. No. R4652)

RIDASCREEN® FAST Casein (Art. No. R4612)

RIDASCREEN® FAST  $\beta$ -Lactoglobulin (Art. No. R4902)

RIDA® Extractor 2 (Art. No. R4613)

Bioavid Lateral Flow Milch/Milk (Art. No. BL613-10 / BL613-25)

### **1. Intended use**

RIDASCREEN®  $\beta$ -Lactoglobulin is a competitive enzyme immunoassay for the quantitative analysis of  $\beta$ -lactoglobulin in hydrolyzed milk products including hypoallergenic baby food. The assay is calibrated to hydrolyzed  $\beta$ -lactoglobulin check samples; the standards consist of native  $\beta$ -lactoglobulin. The assay detects native and processed protein and fragments thereof (during the heating of milk, the native  $\beta$ -lactoglobulin is partially denatured).

### **2. General**

Cow's milk contains approx. 3.2 % proteins (2.6 % caseins and 0.6 % whey proteins). Approx. half of the whey fraction consists of  $\beta$ -lactoglobulin.

The allergen can be present as an ingredient or as a contamination in raw and cooked products. According to the **regulation (EU) No. 1169/2011**, milk and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

According to the Commission Directive 2006/141/EC on infant formulae, the immunoreactive protein shall be less than 1 %.

### 3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with  $\beta$ -lactoglobulin (BLG). Standards, respectively sample solutions and anti-BLG antibodies are added. Free and immobilized BLG compete for the antibody binding sites. After washing, secondary antibodies labeled with peroxidase (enzyme conjugate) are added. These bind to the antibody-BLG-complexes. Any unbound conjugate is then removed by a washing step. Enzyme substrate and chromogen are added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue product. Addition of the stop solution causes a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the  $\beta$ -lactoglobulin concentration in the sample.

### 4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses).

Component	Cap color	Format		Volume
<b>Microtiter plate</b>	-	Ready to use		96 wells
<b>Standard 1</b>	Transparent	Ready to use	0 ng/ml	1.3 ml
<b>Standard 2</b>	Transparent	Ready to use	10 ng/ml	1.3 ml
<b>Standard 3</b>	Transparent	Ready to use	30 ng/ml	1.3 ml
<b>Standard 4</b>	Transparent	Ready to use	90 ng/ml	1.3 ml
<b>Standard 5</b>	Transparent	Ready to use	270 ng/ml	1.3 ml
<b>Standard 6</b>	Transparent	Ready to use	810 ng/ml	1.3 ml
<b>Wash buffer</b>	Brown	<b>Concentrate</b>	<b>10x</b>	100 ml
<b>Conjugate</b>	Red	Ready to use		12 ml
<b>Antibody</b>	White	<b>Concentrate</b>	<b>11x</b>	0.7 ml
<b>Substrate</b>	Green	Ready to use		7 ml
<b>Chromogen</b>	Blue	Ready to use		7 ml
<b>Stop solution</b>	Yellow	Ready to use		14 ml

### 5. Reagents required but not provided

#### 5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal vials

- shaker
- water bath
- laboratory mincer/grinder, mortar, ultra-turrax or homogenizer
- graduated pipettes
- variable 20 - 200 µl and 200 - 1000 µl -micropipettes

## 5.2. Reagents:

- distilled or deionized water

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The colorless chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- any bluish coloration of the colorless chromogen prior to test implementation
- a value of less than 0.8 absorbance units ( $A_{450\text{ nm}} < 0.8$ ) for standard 1

## 9. Preparation of Samples

Working devices such as a mill, glass vials or spatulas must be cleaned before and after each sample preparation to remove any remains of milk and to avoid contamination.

The **wash buffer** for sample dilution, antibody dilution and washing is provided as a 10fold concentrate. Before use the buffer concentrate has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted buffer expires after approx. four weeks at 20 – 25 °C (68 - 77 °F).

- weigh 2 g of a sample in a centrifugal vial
- add distilled water up to a final volume of 50 ml and rotate the vial for 10 min (head over end) at room temperature (20 - 25 °C / 68 - 77 °F)
- dilute the extract 1:20 (1+19) with diluted buffer (e.g. 50 µl sample + 950 µl buffer); the final dilution therefore is 1:500
- use 50 µl per well in the assay

**The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 1 day.**

## 10. Test implementation

### 10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **antibody** (bottle with white cap) is provided as a concentrate. Since the diluted antibody solution has a limited stability, only the amount which actually is needed should be reconstituted. Before pipetting, the antibody solution should be shaken carefully. For reconstitution, the antibody concentrate is diluted 1:11 (1+10) in diluted buffer (see above), e. g. 200 µl concentrate + 2 ml buffer, ready to use sufficient for 4 microtiter strips.

### 10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps. It is recommended to pipette the antibody, the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µl of each standard solution or prepared sample to separate duplicate wells.
3. Add 50 µl of the diluted antibody to each well. Mix gently by shaking the plate manually and incubate for 2 h at room temperature (20 - 25 °C / 68 - 77 °F).
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete

- removal of liquid from the wells. Fill all the wells with 250 µl diluted washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
5. Add 100 µl of conjugate to each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F).
  6. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
  7. Add 50 µl of substrate and 50 µl of chromogen to each well. Mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
  8. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

## 11. Results

A special software, the RIDA<sup>®</sup>SOFT Win / RIDA<sup>®</sup>SOFT Win.net (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN<sup>®</sup> enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

The β-lactoglobulin concentration in ng/ml (ppb) corresponding to the absorbance of each sample is read from the calibration curve. This concentration must be further multiplied by the dilution factor 500, the concentration of β-lactoglobulin actually contained in a sample is then converted into mg/kg (ppm).

### **Calculation example:**

150 ng/ml β-lactoglobulin is read from the calibration curve. After multiplication with the dilution factor 500 the resulting concentration is 75000 ng/g (75 mg/kg) β-Lactoglobulin.

### **In general:**

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like casein or lactose for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

Native, intact  $\beta$ -lactoglobulin is underestimated in the test kit compared to hydrolyzed  $\beta$ -lactoglobulin (in hydrolyzed products more accessible epitopes are available for antibody binding).

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross reactivities and exemplary analyzed matrices are described in the Validation report.

### **Recommendation:**

In order to ensure a high analytical performance we recommend:

- to analyze each sample material in duplicates
- to use also allergen-free and allergen- containing (spiked) samples as test controls
- in case of extremely acid or alkaline samples, the pH should be adjusted to a neutral pH
- due to the multitude of food types, matrix effects cannot be excluded.
- to carry out spiking experiments for an accurate and correct procedure (for more information refer to the validation report)
- further dilute samples with absorbances < standard 6 and analyze again (note the additional dilution factor for calculation!)

**The product information folder with further information is available on request from your local distributor or R-Biopharm AG.**

### **Further application notes:**

- Analysis of gluten or  $\beta$ -lactoglobulin in supplemental enzymes

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG