# **RIDASCREEN<sup>®</sup> FAST $\beta$ -Lactoglobulin**

**Art. No. R4912**

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# RIDASCREEN<sup>®</sup> FAST $\beta$ -Lactoglobulin

## Brief information

RIDASCREEN<sup>®</sup> FAST  $\beta$ -Lactoglobulin (Art. No. R4912) is a sandwich enzyme immunoassay for the quantitative analysis of  $\beta$ -lactoglobulin in rice crispies, chocolate, and sausage.

For the detection of  $\beta$ -lactoglobulin in hydrolysed milk products, including hypoallergenic baby food, the use of the competitive enzyme immunoassay RIDASCREEN<sup>®</sup>  $\beta$ -Lactoglobulin (Art. No. R4901) is recommended.

For the detection of caseins or caseinates in food, the use of RIDASCREEN<sup>®</sup> FAST Casein (Art. No. R4612) or RIDASCREEN<sup>®</sup> FAST Milk is recommended.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit. The test kit is sufficient for 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples) ..... approx. 45 min  
test implementation (incubation time) ..... 30 min

Limit of detection: 0.04 mg/kg (ppm)  $\beta$ -lactoglobulin (0.02 – 0.07)  
depending on the matrix

Limit of quantification: 0.167 mg/kg (ppm)  $\beta$ -lactoglobulin

Specificity: the antibodies specifically detect  $\beta$ -lactoglobulin of cow's milk. There is a cross reactivity to sheep's, goat's, and buffalo's milk

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. corn bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments. During the development of the test kit more than 70 potential cross reactants are tested. Further information is described in the updated validation report.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

### **Related products**

RIDASCREEN®FAST Milk (Art. No. R4652)

RIDASCREEN®FAST Casein (Art. No. R4612)

RIDASCREEN®  $\beta$ -Lactoglobulin (Art. No. R4901)

RIDA® Extractor 2 (Art. No. R4613)

Bioavid Lateral Flow Milch/Milk (Art. No. BL613-10 / BL613-25)

Set of 4 MoniQA Milk Reference controls (Art. No. MQA 122016)

### **1. Intended use**

RIDASCREEN®FAST  $\beta$ -Lactoglobulin (Art. No. R4912) is a sandwich enzyme immunoassay for the quantitative analysis of native and processed  $\beta$ -lactoglobulin in rice crispies, chocolate, and sausage.

### **2. General**

Cow's milk contains 3.2 % proteins which consist of 10 %  $\beta$ -lactoglobulin (leading protein of whey) and 80 % caseins. The most important allergen especially for children is  $\beta$ -lactoglobulin. For adults the caseins become more dominant. Milk can induce allergic reactions at infancy. Whey or milk powder is often added to food products, therefore it is recommended to determine  $\beta$ -lactoglobulin in food. The allergen can be present as an ingredient or as a contamination in raw and cooked products. According to the **regulation (EU) No. 1169/2011**, milk and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

### **3. Test principle**

The wells of the microtiter strips are coated with specific antibodies to  $\beta$ -lactoglobulin. By adding standards and samples to the wells,  $\beta$ -lactoglobulin

present will bind to the specific antibodies. In a washing step components not bound are removed. Then antibody conjugated to peroxidase is added. This antibody conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. The detection of  $\beta$ -lactoglobulin takes place by adding substrate/chromogen solution. The enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the  $\beta$ -lactoglobulin concentration of the sample. The result is expressed in mg/kg  $\beta$ -lactoglobulin.

#### 4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
3 x Extractor 2	Blue	<b>Concentrate</b>	<b>2x</b>	30 ml
Allergen extraction buffer	Green	<b>Concentrate</b>	<b>10x</b>	100 ml
Additive 1	Blue			2 g
Standard 1*	Transparent	Ready to use	0 mg/kg	1.3 ml
Standard 2*	Transparent	Ready to use	0.167 mg/kg	1.3 ml
Standard 3*	Transparent	Ready to use	0.5 mg/kg	1.3 ml
Standard 4*	Transparent	Ready to use	1.5 mg/kg	1.3 ml
Standard 5*	Transparent	Ready to use	4.5 mg/kg	1.3 ml
Wash buffer	Brown	<b>Concentrate</b>	<b>10x</b>	100 ml
Conjugate	Red	Ready to use		7 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 ml
Stop solution	Yellow	Ready to use		14 ml

\*) The dilution factor of 100 for the sample has already been considered when labeling. Therefore, the  $\beta$ -lactoglobulin concentration of samples can directly be read from the standard curve.

## 5. Reagents required but not provided

### 5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal vials
- water bath (60 °C / 140 °F and 100 °C / 212 °F)
- paper filter
- variable 20 - 200 µl and 200 - 1000 µl micropipettes

### 5.2. Reagents:

- distilled water
- 1 M hydrochloric acid (HCl)
- 1 M sodium hydroxide (NaOH)
- for the extraction of pine nuts: bovine serum albumin (BSA, Serva, fraction V, protease free, Art. No. 11926.01)

## 6. Warnings and precautions for the users

This test should be performed only by trained laboratory staff. The instructions for performing the test must be followed strictly.

The Extractor 2 is harmful to health. It contains mercaptoethanol. It should be worked under a chemical hood and skin contact should be avoided (use gloves).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen prior to test implementation
- a value of less than 1.2 absorbance units ( $A_{450\text{ nm}}$ ) for standard 5

## 9. Preparation of Samples

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use. Working devices such as a mill, glass vials or spatulas must be cleaned before and after each sample preparation to remove any remains of milk and to avoid contamination. Make sure to heat the A-AEP in time in the 60 °C (140 °F) water bath (the 100 °C (212 °F) water bath is used for the extraction of the samples).

The **Allergen extraction buffer** is provided as a 10fold concentrate. Before dilution of the buffer concentrate completely dissolve any crystals in a water bath at 37 °C (99°F) and mix well. Following this, dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate 900 ml distilled water). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for approx. **four** weeks.

For the preparation of the **Allergen extraction buffer containing Additive 1 (A-AEP)** weigh 1.35 g of Additive 1 in a glass beaker and add 15 ml 1 M NaOH. Stir until the Additive 1 is solved. Then fill 700 ml diluted Allergen extraction buffer (see above) in a measuring cylinder. Add the 15 ml Additive 1 solution by stirring constantly, transfer eventually residues of the Additive 1 solution into the measuring cylinder by rinsing with diluted Allergen extraction buffer. Adjust the Additive 1 containing Allergen extraction (A-AEP) buffer to pH 9 with 1 M HCl and fill up to 750 ml with diluted Allergen extraction buffer. 750 ml A-AEP buffer is sufficient for 45 samples. The buffer can be used for approx. 3 weeks at room temperature 20 - 25 °C (68 - 77 °F) (do **not** store in the refrigerator). Use clean bottles as particles of dust can initiate crystallization. Discard the buffer if crystals are present.

The **Extractor 2** is provided as 2fold concentrate and has to be diluted 1:2 (1+1) with distilled water (e.g. 30 ml Extractor 2 + 30 ml dist. water). The complete diluted Extractor 2 is sufficient for 45 samples and can be used for approx. 3 month at room temperature (20 - 25 °C).

## 9.1. Sample extraction with Extractor 2

- homogenize a representative amount of the sample (5 - 50 g)

### **solid samples:**

- weigh in 1 g sample (for pine nut add also 0.5 g BSA) and add 4 ml prepared Extractor 2 (see 9.), mix vigorously, close the vial and cook it for 10 min at 100 °C in a water bath
- let the sample cool down shortly
- pre-heat the A-AEP to 60 °C (140 °F)
- add 16 ml heated (60 °C) A-AEP (see 9.) to the cooked sample

### **liquid samples:**

- take 1 ml and add 4 ml prepared Extractor 2 (see 9.), mix vigorously, close the vial and cook it for 10 min at 100 °C in a water bath
- let the sample cool down shortly
- pre-heat the A-AEP to 60 °C (140 °F)
- add 15 ml heated (60 °C) A-AEP to the cooked sample

### **further prepare solid and liquid samples as follows:**

- mix vigorously (shaker)
- extract for 10 min at 60 °C in a water bath  
cool down (e.g. ice water), centrifuge for 10 min / 2500 g if possible at 4 °C (39 °F) and/or filter (alternative: transfer 2 ml of the extract into a reaction vial and centrifuge at high speed for 10 min in a microcentrifuge)
- dilute the particle free supernatant or the filtrate 1:5 (1+4) with diluted Allergen Extraction buffer, without Additive 1 (see 9.), e.g. 100 µl + 400 µl (1:100 final)  
**Remark:** use this diluted supernatant immediately (within 30 min) in the assay. A longer time period may influence the recovery.
- use 100 µl per well in the assay

**The sample extracts before the last dilution step can be stored at room temperature for 4 hours.**

**Sample extracts with Extractor 2 can also be used both in the RIDASCREEN®FAST Casein (R4612) and in the RIDASCREEN®FAST Milk (R4652) assays.**

## 10. Test implementation

### 10.1. Test preparation for the ELISA

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **wash buffer** is provided as a 10fold concentrate. Before dilution of the buffer concentrate completely dissolve any crystals in a water bath at 37 °C (99 °F) and mix well. Following this, dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted wash buffer is stable at 20 - 25 °C (68 - 77 °F) for approx. four weeks.

### 10.2. Test procedure

Do not use more than 3 strips (24 wells) at a time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101 or Mikrotiter Assembly Breakable Strip 1x8, Thermo Scientific) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution 1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Add 100 µl of each standard or prepared sample to separate duplicate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
4. Add 100 µl of the conjugate to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 -77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete

removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.

6. Add 100 µl of the reddish substrate/chromogen solution to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after the addition of stop solution.

## 11. Results

A special software, the RIDASOFT® Win.NET (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN® enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit. Please note the Compliance Criteria, if the curve shape differs from the Certificate.

In comparison with the certificate, higher values of the absorbance (A450 nm) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or allergen contamination.

A further dilution and new detection of the samples is recommended for absorbance values (A450 nm) > standard 5.

### **Please note:**

When working in accordance with the regulation stated, the dilution factor is 100. The allergen concentration can be read directly from the standard curve - the sample dilution factor of 100 is already taken into account for the standard concentrations (see 4. \*).

For sample dilutions of more than 1:100 the further dilution factor must be considered for the calculation of the β-lactoglobulin concentration.

### **Calculation example**

Cow's milk was added to the food sample (remarks to the protein content of cow's milk, see 2. General).

0.5 mg/kg β-lactoglobulin is read from the calibration curve (a multiplication with the dilution factor is not necessary. This corresponds to approx. 5 mg/kg total milk protein ( $0.5 \text{ mg/kg} \times 3.2 \% / 0.3 \% = 5.33 \text{ mg/kg}$ ) and a total milk content of approx. 156 mg/kg ( $100 \% \times 5 \text{ mg/kg} / 3.2 \%$ ).

## **In general:**

Samples tested negative could still contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components such as casein and / or lactose.

As  $\beta$ -lactoglobulin represents about 10 % of the total milk protein, a 0.2 ppm  $\beta$ -lactoglobulin containing sample corresponds to approx. 2 ppm milk protein.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food, proteins may be altered or fragmented, this may have an impact on the recovery.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross-reactivities and exemplary analysed matrices are described in the updated validation report.

## **Recommendation:**

In order to ensure high analytical performance:

- In case of extremely acid or alkaline samples, the pH should be adjusted to a neutral pH
- Use allergen-free and allergen- containing (spiked) samples as test controls  
Due to the multitude of food types, matrix effects cannot be excluded. To ensure an accurate result spike experiments are recommended
- Each sample material should be analyzed in duplicates.
- The recovery for fragmented proteins is reduced in sandwich ELISA, such samples should be analyzed with a competitive ELISA test systems like the RIDASCREEN<sup>®</sup>  $\beta$ -Lactoglobulin (R4901)
- For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result, all cross reactivities and exemplarily analyzed matrices are described in the validation report
- For details using the ThunderBolt<sup>®</sup> / Bolt automation please contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

**The product information folder with further information is available on request from your local distributor or R-Biopharm AG.**

## **Futher Applications**

- Extraction for heavily liquid absorbing matrices using RIDA<sup>®</sup>Extractor 2 (Art. No. R4613)

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG