# **RIDASCREEN<sup>®</sup> FAST Gliadin sensitive**

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Enzyme immunoassay for the quantitative of gliadins and corresponding prolamins

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Please contact for questions and further information:

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RIDA® and RIDASCREEN®  
are registered trademarks of R-Biopharm AG  
Manufacturer: R-Biopharm AG, Darmstadt, Germany

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## Produktangebot

RIDASCREEN® Gliadin (R7001)  
RIDASCREEN®FAST Gliadin (R7061)  
RIDASCREEN® Gliadin competitive (R7021)  
Cocktail (patented) (R7006/R7016)  
Cocktail ECO (R7080)  
RIDA® Extraction Solution (colorless) (R7098)  
Set of 3 processed Gliadin Assay Controls (R7012)  
RIDA®QUICK Gliadin (R7003/R7004/R7005)  
SureFood® PCR Allergen Gluten

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# RIDASCREEN®FAST Gliadin sensitive

## Brief information

RIDASCREEN®FAST Gliadin sensitive (R7051) is a sandwich enzyme immunoassay for the determination of prolamins from gluten containing cereals in gluten-free declared food (cookies, rice flour, oat flour, sausages, quinoa flour, chocolate, millet flour).

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit. The test kit is sufficient for 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization and extraction
Standard material:	The RIDASCREEN® standard material is calibrated to the standard of the Prolamin Working Group.
Time requirement:	sample preparation Cocktail (patented) (for 10 samples)..... approx. 2 h Cocktail ECO (for 10 samples) ..... approx. 35 min. test implementation (incubation time) ..... 30 min.
Limit of detection:	0.2 mg/kg (ppm) gliadin bzw. 0.4 mg/kg (ppm) gluten (depending on the matrix)
Limit of quantification:	1.25 mg/kg (ppm) gliadin or to 2.5 mg/kg (ppm) gluten
Specificity:	The monoclonal antibody R5 reacts with the gliadin-fractions from wheat and corresponding prolamins from rye and barley.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments. Further information is described in the updated validation report.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework

conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded from the website <http://www.r-biopharm.com/products/food-feed-analysis>.

## Related products

RIDASCREEN® Gliadin (R7001)

RIDASCREEN®FAST Gliadin (R7061)

RIDASCREEN® Gliadin competitive (R7021)

Cocktail (patented) (R7006/R7016)

Cocktail ECO (R7080)

RIDA® Extraction Solution (colorless) (Art. No. R7098)

Set of 3 processed Gliadin Assay Controls (Art. No. R7012)

RIDA®QUICK Gliadin (R7003/R7004/R7005)

SureFood® PCR Allergen Gluten

## 1. Intended use

RIDASCREEN®FAST Gliadin sensitive (R7051) is a sandwich enzyme immunoassay for the determination of prolamins from gluten containing cereals in as gluten-free declared food (cookies, rice flour, oat flour, sausages, quinoa flour, chocolate, millet flour).

The sample preparation using the **Cocktail (patented)** (R7006/R7016) is the official R5-Mendez method according to the Codex Alimentarius and the AOAC. The faster sample preparation using the environmental-friendly **Cocktail ECO** (R7080) is convenient for the screening of samples. The Cocktail ECO has an extraction efficiency of approx. 70 - 110% compared to Cocktail (patented).

## 2. General

The use of wheat flour and gluten in foodstuffs is extremely common because of their heat stability and useful effects on e.g. texture, moisture retention and flavour. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley.

Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

According to the Codex Alimentarius (CODEX STAN 118/1979) two categories for labeling of food according to the gluten content exist:

1.) Food products which contain less than 20 mg/kg can be labeled as "**gluten-free**".

2.) Food products labeled as "**very low gluten**" can have a gluten content above 20 and up to 100 mg/kg.

The threshold of 20 mg/kg has been adopted by many national legislations in many countries.

### 3. Test principle

The basis of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific R5 antibodies against gliadins. By adding the standard or sample solution to the wells, present gliadin will bind to the specific capture antibodies. The result is an antibody-antigen-complex. Components not bound by the antibodies are then removed in a washing step. Then R5 antibody conjugated to peroxidase is added. This conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. Enzyme substrate/chromogen are added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the gliadin concentration of the sample.

### 4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
<b>Microtiter plate</b>	-	Ready to use	-	96 Kavitäten
<b>Buffer</b>	Transparent	Ready to use	-	2 x 100 ml
<b>Standard 1</b>	Transparent	Ready to use	0 ng / ml gliadin	1.3 ml
<b>Standard 2</b>	Transparent	Ready to use	2.5 ng / ml gliadin	1.3 ml
<b>Standard 3</b>	Transparent	Ready to use	5 ng / ml gliadin	1.3 ml
<b>Standard 4</b>	Transparent	Ready to use	10 ng / ml gliadin	1.3 ml
<b>Standard 5</b>	Transparent	Ready to use	20 ng / ml gliadin	1.3 ml
<b>Wash buffer</b>	Brown	<b>Concentrate</b>	<b>10x</b>	100 ml
<b>Conjugate</b>	Red	Ready to use	-	10 ml
<b>Substrate/Chromogen</b> Red chromogen Pro	Brown	Ready to use	-	10 ml
<b>Stop solution</b>	Yellow	Ready to use	-	14 ml

## 5. Materials required but not provided

### 5.1. Equipment

- microtiter plate spectrophotometer (450 nm)
- centrifuge, centrifugal vials (e.g. Brand 10742512)
- shaker
- laboratory mincer / grinder, pestle and mortar, ultra-turrax or homogenizer
- water bath (50 °C / 122 °F)
- graduated pipettes
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes

### 5.2. Reagents

- distilled or deionized water
- gluten-free **skim milk powder** (food quality)
- **Cocktail (patented)** (R7006/R7016, 105 ml/1000 ml) and **ethanol solution (80 %)**: i.e. add 120 ml ethanol p.a. to 30 ml distilled water and shake well  
or
- **Cocktail ECO** (R7080), and **ethanol solution (80 %)**: i.e. add 120 ml ethanol p.a. to 30 ml distilled water and shake well

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instructions for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen prior to test implementation
- a value of less than 1.2 absorbance units ( $A_{450\text{ nm}}$ ) for standard 5

## 9. Preparation of Samples

### 9.1. Preliminary comments

Airborne cereal dust and dirty laboratory equipment lead to gliadin contamination of the assay. Therefore, before starting and during the assay wear gloves.

- clean surfaces, glass vials, mincers and other equipment with 40 % ethanol or 2-propanol (see 5.2.)
- carry out the sample preparation in a room isolated from the ELISA procedure
- check for gliadin contamination of reagents and equipment with the test strips RIDA<sup>®</sup>QUICK Gliadin (R7003/R7004/R7005)
- it is recommended to work **under a chemical hood**, because of  $\beta$ -mercaptoethanol in the Cocktail (patented)
- $\beta$ -mercaptoethanol can disturb the ELISA, therefore dilute the samples **at least 1:500** (recommendation 1:500 for samples with approx. 10 mg/kg gluten and 1:2500 for samples with approx. 50 mg/kg gluten)

### 9.2. Extraction with the Cocktail (patented) (R7006 / R7016, official AOAC method)

Homogenize well a sufficient amount (at least 5 g or 5 ml) of sample (grind it thoroughly to powder and mix well or mix well the solution respectively).

- **liquid food samples:** use 0.25 ml of the homogenized sample and add 2.5 ml of the Cocktail (patented), close the vial and mix well
- **other food samples (e.g. soy and quinoa containing samples):** weigh 0.25 g of the homogenized sample and add 2.5 ml of the Cocktail (patented), close the vial and mix well
- **tannin and polyphenol containing food samples (e.g. chocolate, coffee, cocoa, chestnut flour, buckwheat, millet and spices):** weigh 0.25 g of the homogenized sample, add 0.25 g of skimmed milk powder and add 2.5 ml of the Cocktail (patented), close the vial and mix well
- **meat and sausages:** in these matrices gliadin may be not distributed evenly, therefore, weigh 50 g sample and homogenize: weigh 0.25 g of the

homogenized sample and add 2.5 ml of the Cocktail (patented), close the vial and mix well

- **oat samples:** gliadin may not be distributed evenly, furthermore the samples are difficult to homogenize. Therefore, homogenize 200g, then carry out the extraction with at least the fourfold amount of reagents: weigh 1 g of the homogenized sample and add 10 ml of the Cocktail (patented), close the vial and mix well

**Please further extract all samples as described in the following:**

- incubate for 40 min at 50 °C (122 °F)
- let the sample cool down and then mix it with 7.5 ml 80 % ethanol (see 5.2.) (for oat samples: 30 ml 80 % ethanol)
- close the vial and shake for 1 h up side down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F)
- centrifuge: 10 min, at least 2500 g, at room temperature (20 - 25 °C / 68 - 77 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged with high speed for 10 min in reaction caps by using a microcentrifuge)
- transfer the supernatant in a screw top vial
- dilute the sample 1:12.5 (1+11.5 / 80 µl + 920 µl) with buffer: the final dilution factor is 500
- use **immediately** 100 µl per well in the assay

**Remark:**

The supernatant obtained after the centrifugation or the filtrate can be stored in a tightly closed vial in the dark at room temperature (20 - 25 °C / 68 - 77 °F) up to eight weeks.

**9.3. Extraction with Cocktail ECO (faster and environmental friendly)**

The sample preparation is described in the test kit insert of the Cocktail ECO (R7080). For this fast and environmental friendly sample preparation, it is **not** necessary to work under a chemical hood.

**10. Test implementation**

**10.1. Test preparation**

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The wash buffer is provided as a 10fold concentrate. Before dilution of the buffer concentrate, dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled

water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted wash buffer is stable at 20 - 25 °C (68 - 77 °F) for approx. four weeks.

## 10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than 3 strips (24 wells) at a time. For an assay using three to six strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101 or Mikrotiter Assembly Breakable Strip 1x8, Thermo Scientific) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard or sample to separate duplicate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
4. Add 100 µl of the conjugate to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 min after addition of stop solution.

## 11. Results

A special software, the RIDASOFT® Win.NET (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN® enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit. The document 'Compliance Criteria' provides criteria for evaluating standard curves. For quality assurance assay controls should be used.

In comparison with the certificate, higher values of the absorbance ( $A_{450 \text{ nm}}$ ) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or Gluten contamination.

A further dilution and new detection of the samples is recommended for absorbance values ( $A_{450 \text{ nm}}$ ) > standard 5.

The gliadin concentration in ng/ml (ppb) is read from the RIDASOFT® Win.NET calibration curve and must be further multiplied by the dilution factor of at least 500. This result is then multiplied by 2 in order to obtain the gluten concentration (gliadin represents 50 % of the proteins present in gluten, Codex Definition). The RIDASOFT® Win.NET (version 1.93 or newer) indicates the results in gliadin and gluten.

### Example:

The absorbance value of a sample corresponds to 10 ng/ml gliadin in the calibration curve. Multiplying by the recommended dilution factor 500 leads to 5000 ng/ml, corresponding to 5 mg/kg (ppm) gliadin, respectively 0.0005 % gliadin. To calculate the gluten content, it is necessary to multiply by factor 2 which results in 10 mg/kg gluten, respectively 0.001 % gluten. This sample is considered to be gluten-free, because the gluten concentration is below 20 mg/kg.

### In general:

Samples tested negative still could contain a gluten contamination below the limit of detection of the assay, or they might contain other components like starch for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross reactivities and exemplary analyzed matrices are described in the validation report.

The protein content and protein composition of wheat, rye and barley cultivars may differ. Varying results are thus to be expected for different cultivars.

### **Recommendation:**

In order to ensure a high analytical performance:

- Each sample material should be analyzed in duplicates
- Use also gluten free and gluten containing (spiked) samples as test controls
- Carry out spiking experiments for an accurate and correct procedure
- Strong acidic or alkaline samples have to be neutralized prior to analysis
- Confirm results with PCR (e.g. SureFood<sup>®</sup> PCR Allergen Gluten)
- In the production of foods such as beer or sourdough, proteins are fragmented. In sandwich ELISAs protein fragments lead to a reduced recovery, such samples should be analyzed with a competitive ELISA test systems like the RIDASCREEN<sup>®</sup> Gliadin competitive (R7021)
- It is not recommended to extrapolate below the limit of quantification
- For details using the automation (e.g. Thunder Bolt<sup>®</sup> / Bolt, GEMINI) please contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de)

### **Further application notes:**

- Sample preparation for processed food with the RIDA<sup>®</sup> Extraction Solution (colorless) (Art. No. R7098) - only after validation
- Sample preparation for raw materials with ethanol
- Sample preparation for polyphenol containing raw materials (e.g. chocolate, coffee, cacao, buckwheat) with fish gelatine and ethanol

Further product information and application notes, please contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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