
Flow Chart for SureFood® PREP Advanced Protocol 1

Art. No. S1053

October 2018



Flow Chart for SureFood® PREP Advanced Protocol 1

Art. No. S1053

October 2018

- (5) Purification of bound nucleic acids &
- (6) Drying of the Spin Filter



Place the **clear Spin Filter** into a **new clear Receiver Tube** (Code R)

Add **550 µL Pre-Wash Buffer** (Code P)

Centrifugation 1 min 12000 rpm

Discard filtrate

Place Spin Filter back into the Receiver Tube (Code R)

Add 550 µL Wash Buffer (Code W)




Centrifugation 1 min 12000 rpm

Discard filtrate

Place Spin Filter back into the Receiver Tube (Code R)



Add $550~\mu L$ Wash Buffer (Code W)




Discard filtrate

Place Spin Filter back into the Receiver Tube (Code R)

Centrifugation 2 min 12000 rpm





Place the **Spin Filter** into a **clear 1.5 ml Receiver Tube** (Code T)

Add **50 µL preheated Elution Buffer** (Code E)

Incubation 3 min 65°C

Centrifugation 1 min 10000 rpm

Discard Spin Filter

The eluted DNA is ready-to-use for the PCR $\,$